

Regular wire now customized

to the Setup Wire with Extension 3-3 to fit a spaced malocclusion


If anterior spacing is present, currently doctors might not be able to insert the first wire without using multiple appointments. To avoid this, the “mushroom” bends in the round initial wires will be distalized to allow for ligation of the wire, so spaces from canine to canine can be closed conveniently.

There is no need for the doctor to calculate. Our digital wire design will apply the new values for any cases with more than 1.5 mm of space in the anteriors.

Questions, please call our Customer Care,

Phone +(49) 5472 9491 182 or contact: info-incognito@mmm.com



 3M Oral Care
Orthodontic Products
2724 South Peck Road
Monrovia, CA 91016 USA
www.3MUnitek.com

Customer Care:
(800) 401-3001

TOP-Service für Lingualtechnik GmbH
A 3M Company
Schledehauser Straße 81
49152 Bad Essen
Germany
Tel.: +49 (0)5472 – 94 91 - 0
Fax: +49 (0)5472 – 94 91 - 19
E-Mail: info-incognito@mmm.com
www.incognito.net

3M, Incognito, and Unitek are
Trademarks of 3M or 3M Deutschland
GmbH. All other trademarks are
property of their respective holders.

© 3M 2015. All rights reserved.