

3M Science.
Applied to Life.™



Venous Leg Ulcers

Simplify your wound care journey
with 3M products

Venous leg ulcers are associated with a variety of conditions, including low to high amounts of wound exudate, microbial bioburden, the potential for peri-wound maceration and edema. Using the **Clean-Protect-Fill-Cover-Compress** approach to wound care, 3M can help you manage all of these challenges to help achieve better care, better outcomes, and better value.

1. Clean

Clean the wound thoroughly according to your facility's protocol.

2. Protect

Help protect peri-wound skin and wound edges from maceration.



3. Fill

Antimicrobial – Help reduce and control bioburden in a wound that shows signs of infection.



Low exudate

Absorbent – Help manage wound drainage and fill wound depth.



High exudate

4. Cover

Cover the wound to manage wound exudate, help protect from contamination and help provide a moist wound healing environment.



5. Compress

Help manage edema with effective compression.



Ordering Information

3M™ Cavilon™ No Sting Barrier Film

| PROTECT | Cat. No. | Size | Items/Box | Boxes/Case |
|---------|----------|-----------------------|-----------|------------|
| | 3343E | 1.0 ml peel-open wand | 25 | 4 |
| | 3344E | 1.0 ml wipe | 30 | 6 |
| | 3345E | 3.0 ml wand | 25 | 4 |
| | 3346E | 28.0 ml spray bottle | 12 | 1 |

3M™ Tegaderm™ Ag Mesh Dressing with Silver

| FILL | Cat. No. | Size | Dressings/Box | Boxes/Case |
|------|----------|-------------------------------|---------------|------------|
| | 90500 | 2 in x 2 in (5 cm x 5 cm) | 5 | 4 |
| | 90501 | 4 in x 5 in (10 cm x 12.7 cm) | 5 | 4 |
| | 90502 | 4 in x 8 in (10 cm x 20 cm) | 5 | 4 |
| | 90503 | 8 in x 8 in (20 cm x 20 cm) | 3 | 4 |

3M™ Tegaderm™ Alginate Ag Silver Dressing

| COVER | Cat. No. | Size | Dressings/Box | Boxes/Case |
|-------|----------|-----------------------------------|---------------|------------|
| | 90301 | 2 in x 2 in (5 cm x 5 cm) | 10 | 4 |
| | 90303 | 4 in x 5 in (10.16 cm x 12.7 cm) | 10 | 4 |
| | 90307 | 1 in x 12 in (2.54 cm x 30.48 cm) | 5 | 4 |

3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing

| COVER | Cat. No. | Size | Dressings/Box | Boxes/Case |
|-------|----------|---|---------------|------------|
| | 90600 | 2 in x 2 in (5 cm x 5 cm) | 10 | 4 |
| | 90601 | 4 in x 4 in (10 cm x 10 cm) | 10 | 4 |
| | 90602 | 4 in x 8 in (10 cm x 20 cm) | 5 | 6 |
| | 90603 | 8 in x 8 in (20 cm x 20 cm) | 5 | 6 |
| | 90604 | 3.5 in x 3.5 in (8.9 cm x 8.9 cm) (Fenestrated) | 10 | 4 |
| | 90605 | 4 in x 24 in (10 cm x 60 cm) (Roll) | 6 | 6 |

3M™ Tegaderm™ Silicone Foam Border Dressing

| COVER | Cat. No. | Overall Dressing Size | Dressings/Box | Boxes/Case |
|-------|----------|----------------------------------|---------------|------------|
| | 90640 | 3 in x 3 in (7.5 cm x 7.5 cm) | 10 | 6 |
| | 90641 | 4 in x 4 in (10 cm x 10 cm) | 10 | 6 |
| | 90642 | 6 in x 6 in (15 cm x 15 cm) | 10 | 4 |

3M™ Coban™ 2 Layer Compression System

| COMPRESS | Cat. No. | Size | Black Nylon Stocking | Rolls/Box | Boxes/Case |
|----------|--|--|----------------------------|----------------------------|------------|
| | 2094N | Comfort Layer - 4 in x 2.9 yd (10 cm x 2.7 m) Compression Layer - 4 in x 5.1 yd (10 cm x 4.7 m) | 1 | 2 rolls 1 of each layer | 8 |
| 20096 | Comfort Layer - 6 in x 3.8 yd (15 cm x 3.5 m) Compression Layer - 6 in x 4.9 yd (15 cm x 4.5 m) | - | 2 rolls 1 of each layer | 8 | |

3M™ Coban™ 2 Layer Lite Compression System

| COMPRESS | Cat. No. | Size | Black Nylon Stocking | Rolls/Box | Boxes/Case |
|----------|----------|--|----------------------|----------------------------|------------|
| | 2794N | Comfort Layer - 4 in x 2.9 yd (10 cm x 2.7 m) Compression Layer - 4 in x 5.1 yd (10 cm x 4.7 m) | 1 | 2 rolls 1 of each layer | 8 |

For more information, contact your 3M Critical & Chronic Care Solutions Representative or call the 3M Health Care Customer Helpline at 1-800-364-3577. These products can be ordered from your local distributor.



Available in Canada from:

Critical & Chronic Care Solutions Division

3M Canada Company

P.O. Box 5757
London, ON N6A 4T1
1-800-364-3577
www.3M.ca/skinwound

3M Health Care

2510 Conway Avenue
St. Paul, MN, USA
55144-1000

3M Deutschland GmbH

Health Care Business
Carl-Schurz-Str. 1
41453 Neuss
Germany

3M, 3M Science. Applied to Life., Coban, Cavilon and Tegaderm are trademarks of 3M. Used under license in Canada. Please recycle. Printed in Canada. © 2015, 3M. All rights reserved. 1511-04014 E

3M™ Tegaderm™ High Integrity and 3M™ Tegaderm™ High Gelling Alginate Dressings

| Cat. No. | Size | Dressings/Box | Boxes/Case |
|----------|-----------------------------|---------------|------------|
| 90112 | 4 in x 4 in (10 cm x 10 cm) | 10 | 5 |
| 90114 | 4 in x 8 in (10 cm x 20 cm) | 5 | 4 |
| 90120 | 12 in (30.4 cm) (Rope) | 5 | 4 |
| 90212 | 4 in x 4 in (10 cm x 10 cm) | 10 | 5 |
| 90220 | 12 in (30.4 cm) (Rope) | 5 | 4 |