

Inlays or Onlays

Implants

Clear Aligners



PRODUCT SELECTOR

See product Instructions for Use for complete information.


Choose our best options for your procedure



Edentulous Cases

Larger Reconstruction

1-2 Unit Restorations



Inlays or Onlays


Implants

Larger Reconstruction

1-2 Unit Restorations

Clear Aligners

Edentulous Cases



- Impregnum™ Putty Super Quick (LB)
- Impregnum™ Putty Super Quick (HB or MB)
- Impregnum™ Putty Super Quick (LB)
- Impregnum™ Putty Super Quick (HB or MB)
- Impregnum™ Putty Super Quick (LB)
- Impregnum™ Putty Super Quick (HB or MB)

- Imprint™ 4 Light
- Imprint™ 4 Light Putty
- Imprint™ 4 Light Heavy
- Imprint™ 4 (Penta™) Super Quick (LB)
- Imprint™ 4 (Penta™) Super Quick (HB)
- Impregnum™ Super Quick (LB)
- Impregnum™ Super Quick (HB or MB)
- Impregnum™ Super Quick (MB)

- Imprint™ 4 Super Quick Light
- Imprint™ 4 Super Quick Heavy
- Impregnum™ Super Quick (LB)
- Impregnum™ Super Quick (HB or MB)
- Impregnum™ Super Quick (MB)
- Impregnum™ Super Quick (LB)
- Impregnum™ Super Quick (HB or MB)
- Impregnum™ Super Quick (MB)

- Imprint™ 4 Super Quick Light
- Imprint™ 4 Super Quick Heavy
- Impregnum™ Super Quick (LB)
- Impregnum™ Super Quick (HB or MB)
- Impregnum™ Super Quick (MB)
- Impregnum™ Super Quick (LB)
- Impregnum™ Super Quick (HB or MB)
- Impregnum™ Super Quick (MB)

- Imprint™ 4 Light
- Imprint™ 4 Light Putty
- Imprint™ 4 Light Heavy
- Imprint™ 4 (Penta™) Super Quick (LB)
- Imprint™ 4 (Penta™) Super Quick (HB)
- Imprint™ 4 (Penta™) Super Quick (MB)

- Imprint™ 4 Super Quick Light
- Imprint™ 4 Super Quick Heavy
- Impregnum™ Super Quick (LB)
- Impregnum™ Super Quick (HB or MB)
- Impregnum™ Super Quick (MB)
- Impregnum™ Super Quick (LB)
- Impregnum™ Super Quick (HB or MB)
- Impregnum™ Super Quick (MB)

3M™ Imprint™ 4 Penta™ Putty
VPS Impression Material

3M™ Imprint™ 4
VPS Impression Material

3M™ Impregum™ Super Quick
Polyether Impression Material



3M's family of impression materials

3M, Impregum, Imprint and Penta are trademarks of 3M
or 3M Deutschland GmbH. Used under license in Canada.
© 3M 2018. All rights reserved. 70-2013-0736-3



3M Science.
Applied to Life.™

Impregum™ Penta™ (MB)
Working/Setting Time: 1:45/4:15*

*Setting time from start of mixing 6:00 min.


3M™ Impregum™ Penta™
Polyether Impression Material

3M™ Imprint™ 4 Super Quick
VPS Impression Material

3M™ Impregum™ Soft
Polyether Impression Material

3M™ Imprint™ 4 Penta™ Putty
VPS Impression Material

Imprint™ 4 Penta™ Putty
Working/Setting Time: 1:30/2:30
* Shown with Imprint™ 4 Light




Imprint™ 4 (Penta™) Heavy
Working/Setting Time: 2:00/2:00
Imprint™ 4 Regular
Working/Setting Time: 1:00/2:00
Imprint™ 4 Light
Working/Setting Time: 1:00/2:00



3M™ Impregum™ Super Quick
Polyether Impression Material

Impregum™ Super Quick (LB)
Working/Setting Time: 0:45/2:00
Impregum™ Super Quick (MB)
Working/Setting Time: 0:45/2:00
Impregum™ Penta™ Super Quick (HB)
Working/Setting Time: 0:45/2:00



3M™ Impregum™ Soft
Polyether Impression Material

Impregum™ Penta™ Soft (HB)
Working/Setting Time: 1:45/4:15*
Impregum™ Penta™ Soft (MB)
Working/Setting Time: 1:45/4:15*
Impregum™ Penta™ Soft (LB)
Working/Setting Time: 1:45/4:15*




Impregum™ Penta™ (MB)
Working/Setting Time: 1:45/4:15*

3M™ Impregum™ Penta™
Polyether Impression Material

*Setting time from start of mixing 6:00 min.

Imprint™ 4 (Penta™) Super Quick Heavy
Working/Setting Time: 1:15/1:15
Imprint™ 4 Super Quick Regular
Working/Setting Time: 0:35/1:15
Imprint™ 4 Super Quick Light
Working/Setting Time: 0:35/1:15



3M™ Imprint™ 4
VPS Impression Material

3M™ Imprint™ 4 Super Quick
VPS Impression Material