

3M 拋棄式包布系列

3M Steri-Drape Surgical Drape Packs

衛署醫器製字第 003310 號

3M 拋棄式包布設計不但使用符合國際規範之手術包布材質、並依臨床對手術感染防護要求，統合 ” 隔離、無菌區、無菌表面 ” 之概念設計出功能性包布協助降低手術感染機會。

3M 拋棄式包布材質融合 3M 不織布、黏膠、與薄膜等技術打造下列材質特性以符合病患手術安全與感染防護：

- 僅供單次使用
- 滅菌包裝
- 有效阻絕液體
- 具抗菌性
- 具防滲透性
- 韌性強
- 低棉絮
- 兼顧手術環境安全與環保

使用說明

警告：不可透過防護單進行電擊，會形成電弧或無效電擊。

注意：不可重複滅菌使用。

手術包布：

液體防護等級：Level IV (符合 ANSI/AAMI PB70:2003)

Pack hydrostatic $\geq 100\text{cmH}_2\text{O}$ (測試依據 EN 20811:1992:EN13795-2:2002)

Resistance of material penetration of synthetic blood (測試依據 ASTM F1670)

Linting $< \log 4.0$ (測試依據 ISO9073-10:2003)

手術衣：

液體防護等級：Level III (符合 ANSI/AAMI PB70:2003)

Pack hydrostatic $\geq 100\text{cmH}_2\text{O}$ (測試依據 EN 20811:1992:EN13795-2:2002)

Linting $< \log 4.0$ (測試依據 ISO9073-10:2003)

操作方式：

1. 依各醫院及手術常規進行病患皮膚消毒。
2. 確認皮膚乾燥完全
3. 使用前移除黏膠範圍上背膠膠帶。
4. 將 3M 無菌防水防護單鋪陳於病人開刀區。

儲存條件：最好儲存於室溫下，避免過度的溫度與濕度。

保存期間：2 年

3M 拋棄式包布系列型號如下，各科室可依據手術特性、方式、與病患需求選擇合適的拋棄式包布單張(Individual drape)或包布配件 (Drape Accessory)。

產品規格

Drape Packs

Catalog Number	Pack Description	Item Included (cm x cm)
O-001T	Standard Orthopedic Pack (標準骨科包)	9099T(3)- OP Tape (膠條) 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-001T1	Standard Orthopedic Pack with Surgical Gown (標準骨科包含手術衣)	9099T(3)- OP Tape (膠條) 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-001T2	Standard Orthopedic Pack with Surgical Gown/Udrape (標準骨科包含手術衣/U 型隔離巾)	9099T(3)- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)

		7693T- Surgical Gown(拋棄式手術衣) length: 148 1015T- U- Drape (U 型隔離巾) 120x130 (± 5)
O-001T3	Standard Orthopedic Pack with Surgical Gown (2 pieces 9087) (標準骨科包含手術衣/2 中單)	9099T(3)- OP Tape (膠條), 10x55 (± 5) 9087T(2) -Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-001T4	Standard Orthopedic Pack (2 pieces 9087, 2 pieces 9088) (標準骨科包含 2 中單/2 大單)	9099T(3)- OP Tape (膠條), 10x55 (± 5) 9087T(2) -Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T(2) - Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-001T5	Standard Orthopedic Pack with U Drape (標準骨科包含 U 型隔離巾)	9099T(3)- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5)
O-002T	Orthopedic Pack without Mayo Stand Cover (骨科包無梅歐式單)	9099T(3) - OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5)

		9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-002T1	Orthopedic Pack without Mayo Stand Cover w/ Surgical Gown (骨科包無梅歐式單含手術衣)	9099T(3) - OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-002T2	Orthopedic Pack without Mayo Stand Cover w/ Surgical Gown/U drape (骨科包無梅歐式單含手術衣/U 型隔離巾)	9099T(3) - OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148 1015T- U- Drape (U 型隔離巾) 120x130 (± 5)
O-002T3	Orthopedic Pack without Mayo Stand Cover w/ Surgical Gown/(2 pieces 9087) (骨科包無梅歐式單含手術衣/2 中單)	9099T(3) - OP Tape (膠條), 10x55 (± 5) 9087T(2) -Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-002T4	Orthopedic Pack without Mayo Stand Cover w/ (2 pieces	9099T(3) - OP Tape (膠條), 10x55 (± 5) 9087T(2) -Adhesive Drape Sheet (中單) 183x183 (± 5)

	9087,2 pieces 9088) (骨科包無梅歐式單含 2 中單 /2 大單)	9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T(2)- Adhesive Drape Sheet (大單) 250x152 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-002T5	Orthopedic Pack without Mayo Stand Cover w/ U Drape (骨科包無梅歐式單含 U 型隔 離巾)	9099T(3) - OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5)
O-003T	Orthopedic Pack without Mayo Stand and Extremity Cover (骨科包無梅歐式單及腳套)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-003T1	Orthopedic Pack without Mayo Stand and Extremity Cover w/ Surgical Gown (骨科包無梅歐式單及腳套, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-003T2	Orthopedic Pack without Mayo Stand and Extremity Cover w/ Surgical Gown (骨科包無梅歐式單及腳套,	9099T- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5)

	含手術衣/U 型隔離巾)	9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148 1015T- U- Drape (U 型隔離巾) 120x130 (± 5)
O-003T3	Orthopedic Pack without Mayo Stand and Extremity Cover w/ Surgical Gown/2 pieces 9087 (骨科包無梅歐式單及腳套, 含手術衣)	9099T, - OP Tape (膠條), 10x55 (± 5) 9087T(2) -Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-003T4	Orthopedic Pack without Mayo Stand and Extremity Cover w/(2 pieces 9087, 2pieces 9088) (骨科包無梅歐式單及腳套, 含 2 中單, 2 大單)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T(2) -Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T(2) - Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-003T5	Orthopedic Pack without Mayo Stand and Extremity Cover w/ U Drape (骨科包無梅歐式單及腳套, 含 U 型隔離巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5)
O-004T	Standard Aperture Orthopedic Pack (標準骨科洞巾包)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9041T- Leggings (腳套) 35x60 (± 5)

		9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-004T1	Standard Aperture Orthopedic Pack With Surgical Gown (標準骨科洞巾包含手術衣)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-004T2	Standard Aperture Orthopedic Pack With Surgical Gown/U Drape (標準骨科洞巾包含手術衣及 U 型隔離巾)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-004T3	Standard Aperture Orthopedic Pack With Surgical Gown/U Drape/ 2 pieces 9087 (標準骨科洞巾包含手術衣及 U 型隔離巾, 2 中單)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9087T(2) -Adhesive Drape Sheet (中單) 183x183 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-004T4	Standard Aperture Orthopedic	9099T(2) - OP Tape (膠條), 10x55 (± 5)

	Pack With U Drape (標準骨科洞巾包含 U 型隔離巾)	1015T- U- Drape (U 型隔離巾) 120x130 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-004T5	Standard Aperture Orthopedic Pack With U Drape/ 9087 (標準骨科洞巾包含 U 型隔離巾)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-005T	Aperture Orthopedic Pack without Mayo Stand Cover (骨科洞巾包無梅歐式單)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-005T1	Aperture Orthopedic Pack without Mayo Stand Cover w/ Surgical Gown (骨科洞巾包無梅歐式單含手術衣)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-005T2	Aperture Orthopedic Pack without Mayo Stand Cover w/	9099T(2) - OP Tape (膠條), 10x55 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5)

	Surgical Gown/U Drape (骨科洞巾包無梅歐式單, 含手術衣及 U 型隔離巾)	9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-005T3	Aperture Orthopedic Pack without Mayo Stand Cover w/ Surgical Gown/ (2 pieces 9086) (骨科洞巾包無梅歐式單, 含手術衣, 2 治療巾-9086)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 9086T(2)- Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-005T4	Aperture Orthopedic Pack without Mayo Stand Cover w/ (2 pieces 9086) (骨科洞巾包無梅歐式單, 含 2 治療巾- 9086)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-005T5	Aperture Orthopedic Pack without Mayo Stand Cover w/ U Drape (骨科洞巾包無梅歐式單, 含 U 型隔離巾)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9041T- Leggings (腳套) 35x60 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)

O-006T	Aperture Orthopedic Pack without Mayo Stand Cover and Extremity Cover (骨科洞巾包無梅歐式單及腳套)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-006T1	Aperture Orthopedic Pack without Mayo Stand Cover and Extremity Cover w/ Surgical Gown (骨科洞巾包無梅歐式單及腳套, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-006T2	Aperture Orthopedic Pack without Mayo Stand Cover and Extremity Cover w/ Surgical Gown/ U Drape (骨科洞巾包無梅歐式單及腳套, 含手術衣及 U 型隔離巾)	9099T- OP Tape (膠條), 10x55 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-006T3	Aperture Orthopedic Pack without Mayo Stand Cover and Extremity Cover w/ 2 pieces 9086 (骨科洞巾包無梅歐式單及腳套, 含 2 治療巾-9086)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-006T4	Aperture Orthopedic Pack without Mayo Stand Cover and Extremity Cover w/ U Drape	9099T- OP Tape (膠條), 10x55 (± 5) 1015T- U- Drape (U 型隔離巾) 120x130 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5)

	(骨科洞巾包無梅歐式單及腳套, 含 U 型隔離巾)	9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-007T	Orthopedic Bilateral Drape Pack (骨科雙洞巾包)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-007T1	Orthopedic Bilateral Drape Pack With Surgical Gown (骨科雙洞巾包含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T-Adhesive Drape Sheet (中單) 183x183 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-007T2	Orthopedic Bilateral Drape Pack With Surgical Gown/ 2 pieces 9086 (骨科雙洞巾包含手術衣, 2 治療巾-9086)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-007T3	Orthopedic Bilateral Drape Pack With Surgical Gown/ U Drape	9099T- OP Tape (膠條), 10x55 (± 5) 1015T(2) - U- Drape (U 型隔離巾) 120x130 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5)

	(骨科雙洞巾包含手術衣及 U 型隔離巾)	9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-007T4	Orthopedic Bilateral Drape Pack With U Drape (骨科雙洞巾包含 U 型隔離巾)	9099T- OP Tape (膠條), 10x55 (± 5) 1015T(2) - U- Drape (U 型隔離巾) 120x130 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-007T5	Orthopedic Bilateral Drape Pack With 2 pieces 9086 (骨科雙洞巾包含 2 治療巾 -9086)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-008T	Orthopedic Bilateral Drape Pack without Mayo Stand Cover (骨科雙洞巾包無梅歐式單)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-008T1	Orthopedic Bilateral Drape Pack without Mayo Stand Cover w/ Surgical Gown (骨科雙洞巾包無梅歐式單, 含	9099T- OP Tape (膠條), 10x55 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5)

	手術衣)	9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-008T2	Orthopedic Bilateral Drape Pack without Mayo Stand Cover w/ Surgical Gown (骨科雙洞巾包無梅歐式單, 含手術衣, 2 治療巾- 9086)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-008T3	Orthopedic Bilateral Drape Pack without Mayo Stand Cover w/ Surgical Gown/1015 (骨科雙洞巾包無梅歐式單, 含手術衣及 U 型隔離巾)	9099T- OP Tape (膠條), 10x55 (± 5) 1015T(2) - U- Drape (U 型隔離巾) 120x130 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
O-008T4	Orthopedic Bilateral Drape Pack without Mayo Stand Cover w/ 1015 (骨科雙洞巾包無梅歐式單, 含 U 型隔離巾)	9099T- OP Tape (膠條), 10x55 (± 5) 1015T(2) - U- Drape (U 型隔離巾) 120x130 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-008T5	Orthopedic Bilateral Drape Pack without Mayo Stand Cover w/ 1015 (骨科雙洞巾包無梅歐式單, 含 2 治療巾- 9086)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2)- Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9056T- Bilateral Foot Drape (雙腳洞巾單) 195x300 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228

		(± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
O-009T	Aperture Orthopedic Pack (骨科洞巾包)	9099T(2) - OP Tape (膠條), 10x55 (± 5) 9072T(2)- Drape Sheet/ Back Table Cover (器械桌單)114x150 (± 5) 9085T- Adhesive Towel Drape (治療巾單)100x85 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-001T	Standard Basic Surgical Pack (標準基本手術包)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) or 9089T- Adhesive Drape Sheet (9089T 大單) 270X183 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-001T1	Standard Basic Surgical Pack with Surgical Gown (標準基本手術包含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5)or 9089T- Adhesive Drape Sheet (9089T 大單) 270X183 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148

GS-001T2	Standard Basic Surgical Pack with 9090 Drape (標準基本手術包含特大單)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9090T- Adhesive Drape Sheet (特大單) 300x228 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-001T3	Standard Basic Surgical Pack with 9090 Drape/Surgical Gown (標準基本手術包含特大單及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9090T- Adhesive Drape Sheet (特大單) 300x228 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-002T	Standard Basic Surgical Pack without Mayo Stand Cover (標準基本手術包無梅歐式單)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5) or 9089T- Adhesive Drape Sheet (9089T 大單) 270X183 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-002T1	Standard Basic Surgical Pack without Mayo Stand Cover w/ Surgical Gown (標準基本手術包無梅歐式單, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9088T- Adhesive Drape Sheet (大單) 250x152 (± 5)or 9089T- Adhesive Drape Sheet (9089T 大單) 270X183

		(± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-002T2	Standard Basic Surgical Pack without Mayo Stand Cover w/ 9090 Drape (標準基本手術包無梅歐式單, 含特大單)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9090T- Adhesive Drape Sheet (特大單) 300x228 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-002T3	Standard Basic Surgical Pack without Mayo Stand Cover w/ 9090 Drape (標準基本手術包無梅歐式單, 含特大單及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(2) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 9090T- Adhesive Drape Sheet (特大單) 300x228 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-003T	Standard Laparotomy Pack (標準剖腹包)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T(4)- Adhesive Towel Drape (治療巾單) 75x75 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-003T1	Standard Laparotomy Pack with Surgical Gown (標準剖腹包含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T(4) - Adhesive Towel Drape (治療巾單) 75x75 (± 5)

		9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-003T2	Standard Laparotomy Pack 2 (標準剖腹包 2)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-003T3	Standard Laparotomy Pack 2 with Surgical Gown (標準剖腹包 2 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4) - Adhesive Towel Drape (治療巾單) 60x60 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-003T4	Standard Laparotomy Pack with 2 pieces 9087 (標準剖腹包含 2 中單)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T(2) - Adhesive Drape Sheet (中單) 183x183 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-003T5	Standard Laparotomy Pack with 2 pieces 9087/Surgical Gown (標準剖腹包含 2 中單及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T(2) - Adhesive Drape Sheet (中單) 183x183 (± 5)

	衣)	9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-003T6	Standard Laparotomy Pack without Towel Drape (標準剖腹包無治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-003T7	Standard Laparotomy Pack without Towel Drape w/ Surgical Gown (標準剖腹包無治療巾, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-004T	Laparotomy Pack without Mayo Stand Cover (剖腹包無梅歐式單)	9099T- OP Tape (膠條), 10x55 (± 5) 9084T(4)- Adhesive Towel Drape (治療巾單) 75x75 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-004T1	Laparotomy Pack without Mayo Stand Cover w/ Surgical Gown (剖腹包無梅歐式單, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9084T(4) - Adhesive Towel Drape (治療巾單) 75x75 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148

GS-004T2	Laparotomy Pack without Mayo Stand Cover 2 (剖腹包無梅歐式單 2)	9099T- OP Tape (膠條), 10x55 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4) Adhesive Towel Drape (治療巾單) 60x60 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-004T3	Laparotomy Pack without Mayo Stand Cover 2 w/ Surgical Gown (剖腹包無梅歐式單 2, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4) Adhesive Towel Drape (治療巾單) 60x60 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-004T4	Laparotomy Pack without Mayo Stand Cover w/ 2 pieces 9087 (剖腹包無梅歐式單, 含 2 中單)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T(2)- Adhesive Drape Sheet (中單) 183x183 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-004T5	Laparotomy Pack without Mayo Stand Cover w/ 2 pieces 9087/Surgical Gown (剖腹包無梅歐式單, 含 2 中單及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T(2) - Adhesive Drape Sheet (中單) 183x183 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-004T6	Laparotomy Pack without Mayo Stand Cover and Towel Drape (剖腹包無梅歐式單及治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or

		9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-004T7	Laparotomy Pack without Mayo Stand Cover and Towel Drape w/ Surgical Gown (剖腹包無梅歐式單及治療巾, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9063T- Laparotomy Drape (剖腹大單) 183x285 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-005T	Oval Aperture Surgical Pack (橢圓型洞巾包)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T(2) - Adhesive Drape Sheet (中單)183x183 (± 5) 9050T, 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-005T1	Oval Aperture Surgical Pack (橢圓型洞巾包含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T(2) - Adhesive Drape Sheet (中單)183x183 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-005T2	Oval Aperture Surgical Pack with 9086 Towel Drape (橢圓型洞巾包含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4)- Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-005T3	Oval Aperture Surgical Pack	9099T- OP Tape (膠條), 10x55 (± 5)

	with 9086 Towel Drape (橢圓型洞巾包含 9086 治療巾及手術衣)	9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4)- Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-005T4	Oval Aperture Surgical Pack without Towel Drape (橢圓型洞巾包無治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-005T5	Oval Aperture Surgical Pack without Towel Drape w/ Surgical Gown (橢圓型洞巾包無治療巾, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-006T	Oval Aperture Surgical Pack Without Mayo Stand Cover (橢圓型洞巾包無梅歐式單)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T(2) - Adhesive Drape Sheet (中單)183x183 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-006T1	Oval Aperture Surgical Pack without Mayo Stand Cover w/ Surgical Gown	9099T- OP Tape (膠條), 10x55 (± 5) 9087T(2) - Adhesive Drape Sheet (中單)183x183 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300

	(橢圓型洞巾包無梅歐式單, 含手術衣)	(± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-006T2	Oval Aperture Surgical Pack without Mayo Stand Cover w/ 9086 Towel Drape (橢圓型洞巾包無梅歐式單, 含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-006T3	Oval Aperture Surgical Pack without Mayo Stand Cover w/ 9086 Towel Drape/Surgical Gown (橢圓型洞巾包無梅歐式單, 含 9086 治療巾及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-006T4	Oval Aperture Surgical Pack without Mayo Stand Cover and Towel Drape (橢圓型洞巾包無梅歐式單及治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-006T5	Oval Aperture Surgical Pack without Mayo Stand Cover and Towel Drape (橢圓型洞巾包無梅歐式單及治療巾, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9050T- Adhesive Aperture Drape (橢圓型洞巾單) 228x300 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183

		(± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-007T	Standard Laparoscopy Pack (標準腹腔內視鏡包)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-007T1	Standard Laparoscopy Pack (標準腹腔內視鏡包含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-007T2	Standard Laparoscopy Pack with 9087 Drape (標準腹腔內視鏡包含中單)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-007T3	Standard Laparoscopy Pack with 9087 Drape/ Surgical Gown (標準腹腔內視鏡包含中單及	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡

	手術衣)	洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-007T4	Standard Laparoscopy Pack with 9086 Towel Drape (標準腹腔內視鏡包含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-007T5	Standard Laparoscopy Pack with 9086 Towel Drape/Surgical Gown (標準腹腔內視鏡包含 9086 治療巾及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-007T6	Standard Laparoscopy Pack without Towel Drape (標準腹腔內視鏡包無治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-007T7	Standard Laparoscopy Pack without Towel Drape w/ Surgical Gown	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡

	(標準腹腔內視鏡包無治療巾, 含手術衣)	洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-008T	Laparoscopy Pack without Mayo Stand Cover (腹腔內視鏡包無梅歐式單)	9099T- OP Tape (膠條), 10x55 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-008T1	Laparoscopy Pack without Mayo Stand Cover w/ Surgical Gown (腹腔內視鏡包無梅歐式單, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-008T2	Laparoscopy Pack without Mayo Stand Cover w/ Surgical Gown (腹腔內視鏡包無梅歐式單, 含中單)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T- - Adhesive Drape Sheet (中單) 183x183 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-008T3	Laparoscopy Pack without Mayo Stand Cover w/ Surgical Gown	9099T- OP Tape (膠條), 10x55 (± 5) 9087T- Adhesive Drape Sheet (中單) 183x183 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡

	(腹腔內視鏡包無梅歐式單, 含中單及手術衣)	洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-008T4	Laparoscopy Pack without Mayo Stand Cover w/ 9086 Towel Drape (腹腔內視鏡包無梅歐式單, 含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-008T5	Laparoscopy Pack without Mayo Stand Cover w/ 9086 Towel Drape/Surgical Gown (腹腔內視鏡包無梅歐式單, 含 9086 治療巾及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-008T6	Laparoscopy Pack without Mayo Stand Cover and Towel Drape (腹腔內視鏡包無梅歐式單及治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GS-008T7	Laparoscopy Pack without Mayo Stand Cover and Towel Drape w/ Surgical gown (腹腔內視鏡包無梅歐式單及治療巾, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 1031T- Laparoscopic Drape - Supine Position (腹腔內視鏡洞巾- 仰臥位) 285x182 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183

		(± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GS-009T	Thyroid Pack (甲狀腺手術包)	9051T- Thyroid Drape (甲狀腺切除大單)183x315 (± 5) 9077T- Back Table Cover (器械桌單) 200x150 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) Utility Drape (2) -治療巾 66x38 (± 5)
GS-009T1	Thyroid Pack with 1072 Aperture Pouch (甲狀腺手術包含 1072 集水袋)	9051T- Thyroid Drape (甲狀腺切除大單)183x315 (± 5) 9077T - Back Table Cover (器械桌單) 200x150 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) Utility Drape (2) -治療巾 66x38 (± 5) 1072T- Aperture Pouch Drape (集水袋) 23.6x46
GS-010T	Basic Set up Pack (基礎準備包)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T(4) - Adhesive Towel Drape (治療巾單) 75x75 (± 5) 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
CR-001T	Laparoscopy Pack (Lithotomy) (內視鏡截石術臥位包)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T- Adhesive Towel Drape (治療巾單)75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
CR-001T1	Laparoscopy Pack with Surgical Gown (Lithotomy) (內視鏡截石術臥位包含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9084T- Adhesive Towel Drape (治療巾單)75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內

		視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
CR-001T2	Laparoscopy Pack with 9086 Towel Drape (Lithotomy) (內視鏡截石術臥位包含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
CR-001T3	Laparoscopy Pack with 9086 Towel Drape/Surgical Gown (Lithotomy) (內視鏡截石術臥位包含 9086 治療巾及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
CR-001T4	Laparoscopy Pack without Towel Drape (Lithotomy) (內視鏡截石術臥位包無治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
CR-001T5	Laparoscopy Pack without Towel Drape w/ Surgical Gown (Lithotomy)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內

	(內視鏡截石術臥位包無治療巾, 含手術衣)	視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
CR-002T	Laparoscopy Pack without Mayo Stand Cover (Lithotomy) (內視鏡截石術臥位包無梅歐式單)	9099T- OP Tape (膠條), 10x55 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
CR-002T1	Laparoscopy Pack without Mayo Stand Cover w/ Surgical Gown (Lithotomy) (內視鏡截石術臥位包無梅歐式單, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9084T- Adhesive Towel Drape (治療巾單) 75x75 (± 5) or 9084T1(4)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
CR-002T2	Laparoscopy Pack without Mayo Stand Cover w/ 9086 Towel Drape (Lithotomy) (內視鏡截石術臥位包無梅歐式單, 含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
CR-002T3	Laparoscopy Pack without Mayo Stand Cover w/ 9086	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5)

	Towel Drape/Surgical Gown (Lithotomy) (內視鏡截石術臥位包無梅歐式單, 含 9086 治療巾及手術衣)	5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
CR-002T4	Laparoscopy Pack without Mayo Stand Cover and Towel Drape (Lithotomy) (內視鏡截石術臥位包無梅歐式單及治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
CR-002T5	Laparoscopy Pack without Mayo Stand Cover and Towel Drape w/ Surgical Gown (Lithotomy) (內視鏡截石術臥位包無梅歐式單及治療巾, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 1032T- Laparoscopic Drape - Lithotomy Position (腹腔內視鏡洞巾- 截石術臥位) 228x241 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GYN-001T	Standard Gynecology Pack (標準婦科包)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T(2) - Adhesive Drape Sheet (中單)183x183 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GYN-001T1	Standard Gynecology Pack with Surgical Gown (標準婦科包含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T(2) - Adhesive Drape Sheet (中單)183x183 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183

		(± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GYN-001T2	Standard Gynecology Pack with 9086 Towel Drape (標準婦科包含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GYN-001T3	Standard Gynecology Pack with 9086 Towel Drape/Surgical Gown (標準婦科包含 9086 治療巾及 手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 77693T- Surgical Gown(拋棄式手術衣) length: 148
GYN-001T4	Standard Gynecology Pack without Towel Drape (標準婦科包無治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GYN-001T5	Standard Gynecology Pack without Towel Drape w/ Surgical Gown (標準婦科包無治療巾, 含手術 衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GYN-002T	Standard Gynecology Pack without Mayo Stand Cover	9099T- OP Tape (膠條), 10x55 (± 5) 9087T- Adhesive Drape Sheet (中單)183x183 (± 5)

	(標準婦科包無梅歐式單)	1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GYN-002T1	Standard Gynecology Pack without Mayo Stand Cover (標準婦科包無梅歐式單, 含手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T(2) - Adhesive Drape Sheet (中單)183x183 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GYN-002T2	Standard Gynecology Pack without Mayo Stand Cover w/9086 Towel Drape (標準婦科包無梅歐式單, 含9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4),1080T, 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GYN-002T3	Standard Gynecology Pack without Mayo Stand Cover w/9086 Towel Drape/Surgical Gown (標準婦科包無梅歐式單, 含9086 治療巾及手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GYN-002T4	Standard Gynecology Pack without Mayo Stand Cover and Towel Drape (標準婦科包無梅歐式單及治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GYN-002T5	Standard Gynecology Pack without Mayo Stand Cover and Towel Drape w/ Surgical Gown (標準婦科包無梅歐式單及治	9099T- OP Tape (膠條), 10x55 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or

	療巾, 含手術衣)	9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7693T- Surgical Gown(拋棄式手術衣) length: 148
GU-001T	Standard Urology Pack (標準泌尿外科包)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T- Adhesive Drape Sheet (中單)183x183 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7697R- Surgical Gown(拋棄式加強型手術衣) length: 148
GU-001T1	Standard Urology Pack without Reinforced Surgical Gown (標準泌尿外科包無加強型手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9087T- Adhesive Drape Sheet (中單)183x183 (± 5) 1080T- Gynaecology Drape (婦科洞巾單) 195x179 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GU-001T2	Standard Urology Pack with 9086 Towel Drape (標準泌尿外科包含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7697R- Surgical Gown(拋棄式加強型手術衣) length: 148
GU-001T3	Standard Urology Pack with 9086 Towel Drape without Reinforced Surgical Gown (標準泌尿外科包含 9086 治療巾, 無加強型手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228

		(± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GU-001T4	Standard Urology Pack without Towel Drape (標準泌尿外科包無治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7697R- Surgical Gown(拋棄式加強型手術衣) length: 148
GU-001T5	Standard Urology Pack without Towel Drape and Reinforced Surgical Gown (標準泌尿外科包無治療巾及 手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GU-002T	Standard Urology Pack without Mayo Stand Cover (標準泌尿外科包無梅歐式單)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T- Adhesive Drape Sheet (中單)183x183 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7697R- Surgical Gown(拋棄式加強型手術衣) length: 148
GU-002T1	Standard Urology Pack without Mayo Stand Cover and Reinforced Surgical Gown (標準泌尿外科包無梅歐式單 及加強型手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 9087T- Adhesive Drape Sheet (中單)183x183 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GU-002T2	Standard Urology Pack without	9099T- OP Tape (膠條), 10x55 (± 5)

	Mayo Stand Cover w/ 9086 Towel Drape (標準泌尿外科包無梅歐式單, 含 9086 治療巾)	9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7697R- Surgical Gown(拋棄式加強型手術衣) length: 148
GU-002T3	Standard Urology Pack without Mayo Stand Cover and Reinforced Surgical Gown w/ 9086 Towel Drape (標準泌尿外科包無梅歐式單及加強型手術衣, 含 9086 治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 9086T(4) - Adhesive Towel Drape (治療巾單) 100x75 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GU-002T4	Standard Urology Pack without Mayo Stand Cover and Towel Drape (標準泌尿外科包無梅歐式單及治療巾)	9099T- OP Tape (膠條), 10x55 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5) 7697R- Surgical Gown(拋棄式加強型手術衣) length: 148
GU-002T5	Standard Urology Pack without Mayo Stand Cover /Towel Drape/Reinforced Surgical Gown (標準泌尿外科包無梅歐式單, 治療巾及加強型手術衣)	9099T- OP Tape (膠條), 10x55 (± 5) 1081T- TUR/Urological Drape (泌尿外科單) 164x152 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) or 9078T- Drape Sheet/Back Table Cover (器械桌單) 200x183 (± 5)
GU-003T	Cystoscopy Pack (膀胱鏡檢包)	9099T- OP Tape (膠條), 10x55 (± 5) 9005T- Cysto Drape (膀胱鏡檢單) 152x228 (± 5) 9040T- Leggings (腳套, 一對) 75x114 (± 5) 9069T- Drape Sheet/Back Table Cover (器械桌單) 90x75 (± 5)

		7697R- Surgical Gown(拋棄式加強型手術衣) length: 148
PED-001T	Pediatric Hernia Pack (小兒疝氣包)	9099T- OP Tape (膠條), 10x55 (± 5) 9059T- Pediatric Hernia Sheet with Aperture (小兒疝氣大單) 279x183 (± 5) 9084T1(2)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 9077T- Back Table Cover (器械桌單) 200x150 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5)
PED-001T1	Pediatric Hernia Pack without Mayo Stand Cover (小兒疝氣包無梅歐式單)	9099T- OP Tape (膠條), 10x55 (± 5) 9059T- Pediatric Hernia Sheet with Aperture (小兒疝氣大單) 279x183 (± 5) 9084T1(2)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) 9077T- Back Table Cover (器械桌單) 200x150 (± 5)
TD-200T	Towel Drape Pack (4 pieces) (治療巾包, 4 件裝)	9095T(4)- Absorbent Pad (高吸收治療巾) 50x70 (± 5), Blue Wrap
TD-200T1	Towel Drape Pack (4 pieces) (治療巾包, 4 件裝)	9086T(4)- Adhesive Towel Drape (治療巾單) 100x75 (± 5) Blue Wrap
TD-200T2	Towel Drape Pack (2 pieces) (治療巾包, 2 件裝)	9095T(2)- Absorbent Pad (高吸收治療巾) 50x70 (± 5), Blue Wrap
TD-200T3	Towel Drape Pack (2 pieces) (治療巾包, 2 件裝)	9086T(2)- Adhesive Towel Drape (治療巾單) 100x75 (± 5) Blue Wrap
TD-200T4	Towel Drape Pack (1 pieces) (治療巾包, 1 件裝)	9095T(1) - Absorbent Pad (高吸收治療巾)50x70 (± 5), Blue Wrap
TD-200T5	Towel Drape Pack (1 pieces) (治療巾包, 1 件裝)	9086T(1)- Adhesive Towel Drape (治療巾單) 100x75 (± 5) Blue Wrap
TD-200T6	Towel Drape Pack (4 pieces) (治療巾包, 4 件裝)	9084T(4)- Adhesive Towel Drape (治療巾單)75x75 (± 5) , Blue Wrap
TD-200T7	Towel Drape Pack (2 pieces) (治療巾包, 2 件裝)	9084T(2) - Adhesive Towel Drape (治療巾單)75x75 (± 5) Blue Wrap
TD-200T8	Towel Drape Pack (1 pieces) (治療巾包, 1 件裝)	9084T(1) - Adhesive Towel Drape (治療巾單)75x75 (± 5), Blue Wrap
TD-200T9	Towel Drape Pack (4 pieces) (治療巾包, 4 件裝)	9084T1(4) -Adhesive Towel Drape (治療巾單) 60x60 (± 5) Blue Wrap
TD-200T10	Towel Drape Pack (2 pieces) (治療巾包, 2 件裝)	9084T1(2)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) Blue Wrap

TD-200T11	Towel Drape Pack (1 pieces) (治療巾包, 1 件裝)	9084T1(1)- Adhesive Towel Drape (治療巾單) 60x60 (± 5) Blue Wrap
C-300T	C-Arm Drape Pack (C-Arm 套)	4081T- C-Arm Drape (C-Arm 套)104x150 Cotton String Rubber Band(10)
C-300T1	C-Arm Drape Pack (C-Arm 套)	4082T- C-Arm Drape (C-Arm 套) 90x150 Cotton String Rubber Band(10)
C-300T2	Microscope Drape Pack (顯微鏡套)	4083T-,Microscope Drape (顯微鏡套) 90x135 Cotton String Rubber Band(10)
SA-001	Basic Spinal Anesthesia Pack (基本腰椎麻醉包)	Blue Wrap 1092T-Small Drape with adhesive aperture (小洞巾) 57x66 (± 5) or 1093T- Small Drape with adhesive aperture (小洞巾) 75x75 (± 5) Hand towel Swabstick(4)
SA-002E	Basic Spinal Anesthesia Pack -Epidural (基本腰椎麻醉包-硬腦膜外)	Blue Wrap 1092T-Small Drape with adhesive aperture (小洞巾) 57x66 (± 5) or 1093T- Small Drape with adhesive aperture (小洞巾) 75x75 (± 5) Hand towel Swabstick(4) 1626W Tegaderm Dressing
SA-003	Spinal Anesthesia Pack (腰椎麻醉包)	Blue Wrap 1092T-Small Drape with adhesive aperture (小洞巾) 57x66 (± 5) or 1093T- Small Drape with adhesive aperture (小洞巾) 75x75 (± 5) Hand towel
SA-004E	Basic Spinal Anesthesia Pack - Epidural (腰椎麻醉包-硬腦膜外)	Blue Wrap 1092T-Small Drape with adhesive aperture (小洞巾) 57x66 (± 5) or 1093T- Small Drape with adhesive aperture (小洞巾) 75x75 (± 5)

		Hand towel 1626W Tegaderm Dressing
SA-005	Spinal Anesthesia Pack -V (基本腰椎麻醉包-V)	Blue Wrap 1092T-Small Drape with adhesive aperture (小洞巾) 57x66 (± 5) or 1093T- Small Drape with adhesive aperture (小洞巾) 75x75 (± 5) Hand towel Swabstick(4) 1626W Tegaderm Dressing
VT-001T	Surgical Drape Pack (外科手術包)	9099T(4) - OP Tape (膠條), 10x55 (± 5) 9045T- Adhesive Split Sheet (U 型大單), 228x260 (± 5) 9036T- Mayo Stand Cover (梅歐式單) 58x137 (± 5) 9053T- Extremity Drape (彈性洞巾單) 228x300 (± 5) 9042T- Leggings (腳套) 35x100 (± 5) 9079T -Drape sheet/Back Table Cover (器械桌單) 200x228 (± 5) 1019T- Long U-Drape (長 U 型隔離巾)120x180 (± 5) 9090T- Adhesive Drape Sheet (特大單) 300x228 (± 5) 9095T- Absorbent Pad (高吸收治療巾) 50x70 (± 5) 4801T- C-Arm Drape (C-Arm 套) 104x150 Cotton string Rubber band(10) Hand towel(2) 7698R- - Surgical Gown(拋棄式加強型手術衣) 156

製造廠名稱及地址：台灣明尼蘇達礦業製造股份有限公司(台北市敦化南路二段
95 號 6 至 11, 15 及 16 樓)委託宇仁醫療器材科技股份有限公司
(苗栗縣竹南鎮永貞路 1 段 181 巷 107 號)製造
藥商名稱：台灣明尼蘇達礦業製造股份有限公司
藥商地址：台北市敦化南路二段 95 號 6 至 11, 15 及 16 樓