

A full line of restorative solutions to fit your unique needs.

Versatility & Aesthetics



Filtek™ Supreme XTE Universal Restorative

- Exceptional aesthetics
- Excellent polish retention
- Exceptional handling
- Outstanding strength for anterior and posterior
- Wear resistance similar to enamel*
- True nanotechnology
- 36 shades, 4 opacities



Filtek™ Supreme XTE Flowable Restorative

- Excellent aesthetics
- Excellent cavity adaptation
- Excellent polish retention
- Low wear
- True nanotechnology
- 12 shades

Efficiency & Simplicity



Filtek™ Bulk Fill Posterior Restorative

- One-step placement
- Fast and easy to use
- Excellent handling
- Excellent wear resistance
- Stress relief
- True nanotechnology
- 5 shades

Value

Paradigm™ Nano Hybrid Universal Restorative








- Easy to polish
- Excellent for everyday anterior and posterior restorations
- Excellent handling - holds its shape – does not slump
- Excellent strength and wear resistance
- Easily blends with the surrounding teeth through outstanding chameleon effect.
- Available in 12 popular shades (including 2 opaque shades) that match Vita® Classical Shade Guide to cover all of your basic restorative needs.
- Nano Hybrid technology

*In vitro data on file

3M ESPE has you covered

for every class of direct restoration

	Versatility & Aesthetics		Efficiency & Simplicity	Value
	Filtek™ Supreme XTE Universal Restorative (Incremental up to 2mm)	Filtek™ Supreme XTE Flowable Restorative (Incremental up to 2mm)	Filtek™ Bulk Fill Posterior Restorative (Bulk up to 5mm)	Paradigm™ Nano Hybrid Universal Restorative (Incremental up to 2mm)
Class I 	✓	✓ Liner/Base	✓ Up to 4mm	✓
Class II 	✓	✓ Liner/Base	✓ Up to 5mm	✓
Class III 	✓	✓		✓
Class IV 	✓			✓
Class V 	✓	✓		✓
Minimally Invasive	✓	✓	✓	✓
Stress Bearing	✓		✓	✓
Single Shade Aesthetics	✓	✓	✓	✓
Dual/Multi-Shade Aesthetics	✓			

✓ = Recommended uses

Refer to Instructions for Use for full list of indications.

3M Australia Pty Limited ABN 90 000 100 096

Building A, 1 Rivett Rd, North Ryde NSW 2113
1300 363 878 www.3MESPE.com.au

3M New Zealand Limited

94 Apollo Drive, Rosedale Auckland 0632
Ph: 0800 80 81 82 www.3MESPE.co.nz

3M, ESPE, Filtek and Paradigm are trademarks of 3M or 3M ESPE. Please recycle.
Printed in Australia. © 3M 2014. All rights reserved. PB7100/1014

3M ESPE