

3M Science.
Applied to Life.™

See prophy powder in a whole new way.

Clinpro™ Glycine Prophy Powder



Subgingival plaque biofilm removal by spraying Clinpro™ Glycine Prophy Powder along the gingival margin with a 30° to 60° angle in a steady motion to facilitate subgingival cleaning up to 5 mm pocket depth.*

5 mm

Summary of advantages

- For supra- and subgingival plaque biofilm removal
- Low abrasiveness on dentin, root surfaces and restorations
- Proven to effectively remove biofilm in periodontal pockets
- Removes light to moderate extrinsic stains
- Gentle on restorations, orthodontic brackets, implants and sealants
- Pleasant tasting
- Can be used in commercially available air-polishing equipment
- Effective for areas difficult to reach with paste

Comfortable cleaning for your patients ... confidence in fast and effective biofilm removal for you.

You take special care with your preventive measures to ensure your patients' comfort and continued oral health. Now there's a way to enhance your prophy procedure and advance their care.


Unlike conventional air-polishing powders such as sodium bicarbonate, Clinpro™ Glycine Prophy Powder is made with soft and gentle water-soluble glycine technology for a thorough but gentle cleaning—even subgingivally.

- A fast and comfortable cleaning for your patients
- Safe removal of plaque biofilm from soft root surfaces, dental restorations, implants and around orthodontic brackets
- Gentler on dental tissues than conventional procedures using currettes and scalers

Clinpro™ Glycine Prophy Powder from 3M is a gentle prophy powder that has the longest clinical history of glycine air-polishing powders—used clinically since 2001.



Around orthodontic brackets



On implants and restorations



Around gingival margins

Benefits for you and your patients:

	Clinpro™ Glycine Prophy Powder	Conventional Sodium Bicarbonate Powders	Conventional Prophy Paste
Used to remove subgingival plaque biofilm up to 5 mm pocket depth	X		
Used to remove supragingival plaque biofilm	X	X	X
Made with less abrasive glycine technology	X		
Safe on restorations, implants, orthodontic brackets and sealants	X		X
Water soluble	X	X	X
Pleasant, non-salty taste	X		X
Can be used for patients with sodium restrictions	X		X
Removes light to moderate extrinsic stains	X	X	X
Removes heavy discoloration		X	
Used for supportive maintenance in periodontitis therapy	X		
Used for supportive maintenance for peri-implant therapy	X		



Ordering Information

Item #	Product Information
67008	Clinpro™ Glycine Prophy Powder: 2 bottles – 5.6 oz/160 grams each
12616	C&P Powder: 10 oz bottle

Clinpro™ prophylaxis solutions help your patients, and you, achieve success. Make 3M in-office and at-home preventive products part of your patient care protocol.

www.3M.com/ProphyPowder



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
Phone 1-888-363-3685

3M, ESPE and Clinpro are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Printed in USA. © 3M 2016. All rights reserved. 70-2013-0761-1