



Polish with diamonds. Skip the paste.

Sof-Lex™ Diamond Polishing System

How much time and effort do you spend creating beautiful smiles? Whether you currently use a rubberised finishing and polishing system or an intraoral diamond polish, the process can be time-consuming. And, even with your best effort, the gloss may not last. 3M has a simple solution for both problems, using two of our innovative technologies.

Restore with Filtek™ Supreme XTE Universal Restorative. Brilliant aesthetics is just one reason why doctors use this nanocomposite. Thanks to 3M's true nanotechnology, it is easy to polish and offers excellent polish retention.

Polish with the Sof-Lex™ Diamond Polishing System. Forget the messy paste. Our pre-polishing spiral prepares the restoration for final gloss, while our diamond-impregnated polishing spiral gives your restorations that gorgeous paste-like gloss. The system offers the convenience of a rubberised system while also adapting to all tooth surfaces.

You'll be happy to know that while the spirals are effective, they're also kinder to gingival tissues*—and maintain the integrity and anatomy of your restorations!

When patients leave your office smiling, you'll marvel at how simple it's become to create beautiful, natural-looking aesthetics.

*Compared to other finishing and polishing tools.



Filtek™ Supreme XTE Universal Restorative polished with Sof-Lex™ Diamond Polishing System providing a clear reflection due to the high gloss finish obtained.

You can create a diamond paste-like gloss with just two steps.



1
Sof-Lex™
Pre-Polishing Spiral
(beige)



2
Sof-Lex™
Diamond Polishing Spiral
(pink)

Summary of advantages

- Imparts paste-like gloss in the convenience of a rubberised system
- Unique, flexible shape adapts to all tooth surfaces
- Fast and easy to use
- Multi-use, can be sterilised and reused
- High, long-lasting gloss when used with Filtek™ Supreme XTE Universal Restorative

Don't you love it when you find a product that works great, but is also cost effective?

Then you'll love this! The multi-purpose Sof-Lex™ spirals can be sterilised and reused, making them an economical choice for your practice.

Because abrasive particles are embedded throughout the rubberised spirals, all surfaces of the spiral design can be used: side, top, bottom and edge.

And forget the need to stop and switch between traditional points, cups, discs and brushes, too. The flexible spirals adapt to all tooth surfaces—anterior and posterior, convex and concave.



Sof-Lex™ Pre-Polishing Spiral

The pre-polishing spiral (beige) smooths and removes scratches in restorations that develop during contouring and prepares the surface for high-gloss polishing.



Sof-Lex™ Diamond Polishing Spiral

The diamond polishing spiral (pink) contains diamond particles to create a smooth, high gloss polish.

You're going to love how easy it is to create a beautiful smile!

Ordering Information

| Item # | Product Information |
|---------|--|
| 5092-IM | Sof-Lex™ Diamond Polishing System—Introductory Kit Kit includes: 5 Sof-Lex™ Pre-Polishing Spirals (Beige); 5 Sof-Lex™ Diamond Polishing Spirals (Pink); 1 Mandrel; Instructions for Use; Technique Card |
| 5090 | Sof-Lex™ Pre-Polishing Spirals—Refill Includes: 15 Sof-Lex Pre-Polishing Spirals (Beige) |
| 5091 | Sof-Lex™ Diamond Polishing Spirals—Refill Includes: 15 Sof-Lex Diamond Polishing Spirals (Pink) |



3M Oral Care

3M Australia Pty Limited
Building A, 1 Rivett Rd
North Ryde NSW 2113
1300 363 454
www.3M.com.au

3M New Zealand Limited
94 Apollo Drive
Rosedale Auckland 0632
0800 808 182
www.3M.co.nz

3M, 3M Science. Applied to Life., Filtek and Sof-Lex are trademarks of 3M. © 3M 2016. All rights reserved. Always follow the instructions for use. 3M data on file supports paste-like gloss, multi-use and kinder to gingival tissue references.