

**3M** Science.  
Applied to Life.™

# The surgical clipper with a pivoting head.

Safe for patients. Easy for you.



**3M™ Surgical Clipper 9661L**  
with **Pivoting Head** for use with the  
**3M™ Surgical Clipper Blade 9660**



## The clipper that is safe, reliable and easy to use.

Start with a unique pivoting head. Add a long-lasting battery with a short recharge time, a battery indicator light, an ergonomically designed handle and submersibility. What's the result? A better clipping experience for both you and your patient.



\* When used with the 3M™ Surgical Clipper 9661L, the 3M™ Surgical Clipper Blade 9690 is non-pivoting.



### Experience the difference a pivoting head makes.

When used with the 9660 blade, the unique pivoting head conforms to body contours, making it easy to provide a comfortable clipping experience for patients. You can efficiently cut through all types of head and body hair, wet or dry, in a single pass.



### Low nick rate.

The 3M™ Surgical Clipper 9661L's trusted 9660 blade offers performance advantages where it matters most – safety. The away-from-skin cutting position of the clipper blade minimizes nicks and cuts, which helps reduce the risk of SSIs. A nick rate of less than 2% was observed in a clinical study.<sup>1</sup>



### Lithium-ion battery enhances efficiency.

Clip uninterrupted for a longer time. Charge for a shorter time. Powered by a Li-ion battery, the 3M surgical clipper 9661L provides 160 minutes of run time on a four-hour recharge.



### Submersible for easy cleaning.

Meets IPX7 rating of submersible in water no deeper than one meter for up to 30 minutes.



### Ergonomic shape is easy to hold.



# 3M™ Surgical Clipper 9961L with Pivoting Head

Features	9660 Blade	9690 Blade
Easy-to-attach blade.	✓	✓
Submersible for easy cleaning.	✓	✓
Cuts through all types of head and body hair, wet or dry, in a single pass.	✓	
Individually sealed single-use blade.	✓	✓
Quickly and effectively removes hair in trauma and neuro cases.		✓
Away-from-the-skin cutting blade minimizes potential for nicks or cuts.	✓	✓
Li-ion battery – designed to hold a charge longer.	✓	✓
Ergonomically designed handle that's easy to hold.	✓	✓
Battery indicator light shows when it's time to recharge.	✓	✓
Approximately 160 minutes of run time on a four-hour recharge.	✓	✓
Pivoting head.	✓	

Catalog Number	Product Description	Items per Box	Boxes per Case
9660	3M™ Surgical Clipper Blade 9660	50	1
9661L	3M™ Surgical Clipper 9661L with Pivoting Head	1	1
9662L	3M™ Surgical Clipper Charger for 9661L	1	1
9667L	3M™ Surgical Clipper Starter Kit 9667L: 3M™ Surgical Clipper 9661L and 3M™ Surgical Clipper Charger 9662L	1	1
9690	3M™ Surgical Clipper Blade 9690 (for head and heavy body hair)	20	1

## For more information about 3M™ Surgical Clipper 9661L or other 3M Health Care products:

- Contact your local 3M representative
- Call the 3M Health Care Customer Helpline at 800.228.3957
- Visit [3M.com/medical](http://3M.com/medical)



**3M Health Care**  
**3M Infection Prevention Division**  
 3M Center, Building 275-4E-01  
 2510 Conway Avenue  
 St. Paul, MN 55144 USA

1.800.228.3957  
[3m.com/medical](http://3m.com/medical)

'3M study 05-013632.

Please recycle. Printed in the U.S.A.  
 © 3M 2016. 3M is a trademark of 3M.  
 All rights reserved.  
 70-2011-5613-3