

3M™ Medipore™ Family of Products



Easy does it

The soft, stretchy, conformable tapes, dressings and pre-cut pads that are easy to use and easy on skin.

3M

For over 50 years, healthcare professionals have depended on 3M™ Medical Tapes and Dressings for high-quality performance.

3M™ Medipore™ Family of Products

The 3M Medipore family of products is:

- Latex free
- Hypoallergenic
- Highly breathable
- Easy to handle—no curling or self-adhesion
- Water resistant
- Easy to apply and remove
- Gentle
- Soft and comfortable
- Conformable with excellent cross and diagonal stretch





3M™ Medipore™ H Soft Cloth Surgical Tape and Medipore™ Soft Cloth Surgical Tape

High adhesion, breathable soft cloth tape for a range of applications.

Features and benefits

- Easy-tear perforated rolls
- Gentle and breathable, to maintain skin integrity
- Excellent cross and diagonal stretch
- Hypoallergenic and latex free
- Water resistant

Suggested applications

- When repeated taping is needed
- Securing over areas where skin may stretch (e.g., edema, distention, hematoma formation or movement)
- Securing dressings and ostomy pouches to fragile skin
- Anchoring post-surgical dressings



3M™ Medipore™ Dress-It Pre-Cut Dressing Cover

High adhesion, pre-cut pads used for dressing secured.

Features and benefits

- Pre-cut, easy to use
- Gentle to skin
- Available in 5 sizes

Suggested applications

- As retentive sheeting over gauze



3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing

A gentle, conformable, absorbent dressing.

Features and benefits

- All-in-one sterile dressing – reduces application time
- Conforms to body contours, providing comfort
- Gentle, porous and breathable
- Unique, non-adherent absorbent pad

Suggested applications

- Post-surgical dressings—especially on areas where swelling or movement is expected
- Superficial and partial-thickness wounds
- Lacerations, abrasions and burns
- IV catheter sites
- Light to moderately draining wounds

3M™ Medipore™ Products

To learn more about the Medipore family of products, visit www.3M.com/SkinHealth, contact your 3M Skin & Wound Care Representative or call the 3M Helpline at 1-800-228-3957. Outside of the United States, contact the local 3M subsidiary.

***Disclaimer:** The reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coding for 3M products. Insurers' reimbursement policies can vary, however, and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level.

| Catalog Number | NDC/NHRC Number | Size (in × yd / cm × m) | Each/Box | Boxes/Case | HCPCS Code* |
|---|-----------------|-------------------------|----------|------------|-------------|
| 3M™ Medipore™ H Soft Cloth Surgical Tape | | | | | |
| 2861 | 8333-2861-01 | 1 × 10/2,5 × 9,1 | 24 | 1 | A4452 |
| 2862 | 8333-2862-01 | 2 × 10/5 × 9,1 | 12 | 1 | A4452 |
| 2863 | 8333-2863-01 | 3 × 10/7,6 × 9,1 | 12 | 1 | A4452 |
| 2864 | 8333-2864-01 | 4 × 10/10,1 × 9,1 | 12 | 1 | A4452 |
| 2866 | 8333-2866-01 | 6 × 10/15,2 × 9,1 | 12 | 1 | A4452 |
| 2868 | 8333-2868-01 | 8 × 10/20,3 × 9,1 | 6 | 1 | A4452 |
| 2862S | 8333-2862-03 | 2 × 2/5 × 1,8 | 48 | 1 | A4452 |
| 2864S | 8333-2864-03 | 4 × 2/10,1 × 1,8 | 24 | 1 | A4452 |
| 2866S | 8333-2866-03 | 6 × 2/15,2 × 1,8 | 16 | 1 | A4452 |

| Catalog Number | NDC/NHRC Number | Size (in × in / cm × cm) | Each/Box | Boxes/Case | HCPCS Code* |
|---|-----------------|--------------------------|----------|------------|-------------|
| 3M™ Medipore™ Soft Cloth Surgical Tape | | | | | |
| 2961 | 8333-2961-01 | 1 × 10/2,5 × 9,1 | 24 | 1 | A4452 |
| 2962 | 8333-2962-01 | 2 × 10/5 × 9,1 | 12 | 1 | A4452 |
| 2963 | 8333-2963-01 | 3 × 10/7,6 × 9,1 | 12 | 1 | A4452 |
| 2964 | 8333-2964-01 | 4 × 10/10,1 × 9,1 | 12 | 1 | A4452 |
| 2966 | 8333-2966-01 | 6 × 10/15,2 × 9,1 | 12 | 1 | A4452 |
| 2968 | 8333-2968-01 | 8 × 10/20,3 × 9,1 | 6 | 1 | A4452 |
| 2962S | 8333-2962-03 | 2 × 2/5 × 1,8 | 48 | 1 | A4452 |
| 2964S | 8333-2964-03 | 4 × 2/10,1 × 1,8 | 24 | 1 | A4452 |
| 2966S | 8333-2966-03 | 6 × 2/15,2 × 1,8 | 16 | 1 | A4452 |

| Catalog Number | NDC/NHRC Number | Size (in × in / cm × cm) | Each/Box | Boxes/Case | HCPCS Code* |
|---|-----------------|--------------------------|----------|------------|-------------|
| 3M™ Medipore™ Dress-It Pre-Cut Dressing Covers | | | | | |
| 2954 | 8333-2954-01 | 3 7/8 × 4 5/8 / 10 × 12 | 25 | 4 | A4452 |
| 2955 | 8333-2955-01 | 3 7/8 × 7 7/8 / 10 × 20 | 25 | 4 | A4452 |
| 2956 | 8333-2956-01 | 5 7/8 × 5 7/8 / 15 × 15 | 25 | 4 | A4452 |
| 2957 | 8333-2957-01 | 5 7/8 × 11 / 15 × 28 | 25 | 4 | A4452 |
| 2958 | 8333-2958-01 | 7 7/8 × 11 / 20 × 28 | 25 | 4 | A4452 |

| Catalog Number | NDC/NHRC Number | Dressing Size (in × in / cm × cm) | Pad Size (in × in / cm × cm) | Each/Box | Boxes/Case | HCPCS Code* |
|--|-----------------|-----------------------------------|------------------------------|----------|------------|-------------|
| 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing | | | | | | |
| 3562 | 8333-3562-01 | 2 × 2 3/4 / 5 × 7 | 1 × 1 1/2 / 2,5 × 4 | 50 | 4 | A6219 |
| 3564 | 8333-3564-01 | 2 3/8 × 4 / 6 × 10 | 1 × 2 3/8 / 2,5 × 6 | 50 | 4 | A6219 |
| 3566 | 8333-3566-01 | 3 1/2 × 4 / 9 × 10 | 1 3/4 × 2 3/8 / 4,5 × 6 | 25 | 4 | A6219 |
| 3568 | 8333-3568-01 | 6 × 6 / 15 × 15 | 4 × 4 / 10 × 10 | 25 | 4 | A6219 |
| 3569 | 8333-3569-01 | 3 1/2 × 6 / 9 × 15 | 1 3/4 × 4 / 4,5 × 10 | 25 | 4 | A6219 |
| 3570 | 8333-3570-01 | 3 1/2 × 8 / 9 × 20 | 1 3/4 × 6 / 4,5 × 15 | 25 | 4 | A6219 |
| 3571 | 8333-3571-01 | 3 1/2 × 10 / 9 × 25 | 1 3/4 × 8 / 4,5 × 20 | 25 | 4 | A6219 |
| 3573 | 8333-3573-01 | 3 1/2 × 13 3/4 / 9 × 35 | 1 3/4 × 11 3/4 / 4,5 × 30 | 25 | 4 | A6220 |

www.3M.com/SkinHealth



Skin & Wound Care Division
3M Health Care
 2510 Conway Ave.
 St. Paul, MN 55144
 USA
 1-800-228-3957
www.3M.com/healthcare

3M Canada
 Post Office Box 5757
 London, Ontario N6A 4T1
 Canada
 1-800-364-3577
www.3M.com/ca/healthcare

Please recycle. Printed in U.S.A.
 © 3M 2005, 2012. All rights reserved.
 70-2010-9054-8

3M and Medipore are trademarks of 3M.