




Tears Easily Works Reliably

Transparent, perforated plastic tape designed for dependability and ease of use

- Easy, straight, bi-directional tear allows clinicians to tailor the tape for various dressings and devices
- Easy to handle with gloves
- Good adhesion to skin and tubing
- Transparent, porous
- Hypoallergenic and latex-free
- Water-resistant

Suggested Applications

- Securing tubing and devices, (e.g., catheters, I.V. tubing)
- Anchoring bulky dressings
- When tape width must be customized



3M

3M™ Transpore™ Surgical Tape

Catalog #	NDC/NHRC No.	Size	Rolls/ Box	Boxes/ Case	HCPCS Code
1527-0	8333-1527-00	½ in. × 10 yd./1,25 cm × 9,1 m	24	10	A4452
1527-1	8333-1527-01	1 in. × 10 yd./2,5 cm × 9,1 m	12	10	A4452
1527-2	8333-1527-02	2 in. × 10 yd./5 cm × 9,1 m	6	10	A4452
1527-3	8333-1527-03	3 in. × 10 yd./7,6 cm × 9,1 m	4	10	A4452
1527S-1	8333-1527-11	1 in. × 1 ½ yd./2,5 cm × 1,37 m	100	5	A4452
1527S-2	8333-1527-12	2 in. × 1 ½ yd./5 cm × 1,37 m	50	5	A4452

Disclaimer: The reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coding for 3M products. Insurers' reimbursement policies can vary, however, and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level.

3M.com/MedicalTape



For more information, contact your 3M Health Care Sales Representative, call the 3M Health Care Customer Helpline at 1-800-228-3957, or visit **3M.com/MedicalTape**. This product can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.



Skin & Wound Care Division
3M Health Care
2510 Conway Ave.
St. Paul, MN 55144
USA
1-800-228-3957
www.3M.com/healthcare

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-364-3577
www.3M.com/ca/healthcare

Please recycle. Printed in U.S.A.
© 3M 2001, 2012. All rights reserved.
70-2008-5943-0

3M and Transpore are trademarks of 3M.