



Bend. Transform. Warm.

Our new warming blanket bends and conforms while providing uniform temperatures.

The new 3M™ Bair Hugger™ Multi-position upper body blanket is a customer-inspired innovation designed to maximize patient coverage with one easy-to-use warming solution—because degrees matter, and so does versatility. Its unique bendability enables clinicians to adapt the blanket to multiple surgical positions and procedures, to help maintain normothermia.



Transform the way you warm.

The blanket's unique bendability feature was inspired by the clinician's need to optimize patient body surface coverage in order to help maintain normothermia and drive positive patient outcomes.

Effective and Versatile

Engineered to deliver improved heat transfer in a wide range of surgical procedures¹

Bends and conforms while providing uniform temperatures

Offers optimal patient body surface coverage²

Conformable and Transformable

Offers fast and easy application

Offers greater ability to conform to the patient's body¹

Engineered to provide minimum loftiness

Suitable for perioperative use

Blanket Features

When deployed, the attached clear head drape and two neck vents keep warm air around an intubated patient's head and allow observation

Two resealable hose ports provide flexibility in positioning

Integrated tie strips and continuous adhesive strip can secure the blanket once placed

Material is soft, comfortable, lightweight, and radiolucent

¹ 3M data on file.

² 89% (72 of 81) of clinical users indicated the product offers optimal patient body surface coverage. 3M data on file.

Ordering Information

Item

Product Information

Model 62200

3M™ Bair Hugger™ Multi-position upper body warming blanket
Size: 78" x 24" (198 x 61 cm) Weight: 3.7oz (104g)
Units/Case: 10

To order, please contact your local 3M Bair Hugger sales representative or call Customer Service at 1-800-228-3957. Visit www.bairhugger.com.

