




## A full line of restorative solutions to fit your unique needs.




With a rich legacy of oral care products and services, 3M is one of the world's most trusted companies. Filtek™ restoratives have been used in more than 1 billion restorations worldwide. Our family of Filtek™ dental restoratives lets you choose just the right option to fit any patient situation—from highly esthetic, to highly simple for an easier day.

Versatility & Esthetics		Esthetics & Simplicity	Efficiency & Simplicity	
<p><b>3M™ Filtek™ Supreme Ultra Universal Restorative</b></p> <ul style="list-style-type: none"> <li>• Exceptional esthetics</li> <li>• Natural-looking fluorescence</li> <li>• Excellent polish retention</li> <li>• Exceptional handling</li> <li>• Outstanding strength for anterior and posterior</li> <li>• Wear resistance similar to enamel*</li> <li>• TRUE nanotechnology</li> <li>• 36 shades, 4 opacities</li> </ul>	<p><b>3M™ Filtek™ Supreme Ultra Flowable Restorative</b></p> <ul style="list-style-type: none"> <li>• Excellent esthetics</li> <li>• Excellent polish retention</li> <li>• Excellent adaptation</li> <li>• Excellent wear resistance</li> <li>• TRUE nanotechnology</li> <li>• 12 shades</li> </ul>	<p><b>3M™ Filtek™ Universal Restorative</b></p> <ul style="list-style-type: none"> <li>• Predictable, repeatable esthetics</li> <li>• Long-lasting polish</li> <li>• Excellent, creamy handling</li> <li>• Low shrinkage stress/stress relief</li> <li>• Excellent wear resistance</li> <li>• TRUE nanotechnology</li> <li>• 8 designer shades plus an Extra White (XW) cover the 19 VITA classical and bleach shades**</li> <li>• Pink Opaquer with excellent masking ability</li> </ul>	<p><b>3M™ Filtek™ One Bulk Fill Restorative</b></p> <ul style="list-style-type: none"> <li>• Increased opacity for improved esthetics</li> <li>• One-step placement up to 5 mm</li> <li>• Excellent handling</li> <li>• Excellent adaptation</li> <li>• Stress relief</li> <li>• Excellent wear resistance</li> <li>• TRUE nanotechnology</li> <li>• 5 shades</li> </ul>	<p><b>3M™ Filtek™ Bulk Fill Flowable Restorative</b></p> <ul style="list-style-type: none"> <li>• One-step placement up to 4 mm</li> <li>• Easy adaptation</li> <li>• Low shrinkage stress</li> <li>• Excellent wear resistance</li> <li>• 4 shades</li> </ul>
				

\*In vitro data on file.

\*\*3M internal data

We have you covered for every class of direct restoration.

	Versatility & Esthetics		Esthetics & Simplicity	Efficiency & Simplicity	
	<b>3M™ Filtek™ Supreme Ultra Universal Restorative</b> (Incremental up to 2 mm)	<b>3M™ Filtek™ Supreme Ultra Flowable Restorative</b> (Incremental up to 2 mm)	<b>3M™ Filtek™ Universal Restorative</b> (Incremental up to 2 mm)	<b>3M™ Filtek™ One Bulk Fill Restorative</b> (Bulk up to 5 mm)	<b>3M™ Filtek™ Bulk Fill Flowable Restorative</b> (Incremental up to 4 mm)
<b>Class I</b> 	✓	✓ (Liner/Base)	✓	✓ (up to 4 mm)	✓ (Liner/Base)
<b>Class II</b> 	✓	✓ (Liner/Base)	✓	✓ (up to 5 mm)	✓ (Liner/Base)
<b>Class III</b> 	✓	✓	✓		✓
<b>Class IV</b> 	✓		✓		
<b>Class V</b> 	✓	✓	✓		✓
<b>Minimally Invasive</b>	✓	✓	✓	✓	✓
<b>Stress Bearing</b>	✓		✓	✓	

✓ = recommended use

Refer to Instructions for Use for full list of indications.



[www.3M.com/Filtek](http://www.3M.com/Filtek)

3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA  
Phone 1-800-634-2249  
Web 3M.com/dental

3M Canada  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
Phone 1-888-363-3685

3M and Filtek are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2019. All rights reserved. 70-2013-0807-2