


# One Bulk Fill Restorative


## Anterior Restorations




1




2




3



4



5



6

# One Bulk Fill Restorative

## Anterior Restorations



Restoration Type	Increment Depth	Halogen Lights (with output of 550-1000 mW/cm <sup>2</sup> )	LED Lights (with output 1000- 2000 mW/cm <sup>2</sup> )
Anterior	≤ 3 mm.	20 sec.	10 sec.

3M, ESPE and Filtek are trademarks of 3M or 3M Deutschland GmbH.  
Used under license in Canada. © 3M 2016. All rights reserved.



**3M ESPE**  
Dental Products  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA



**3M Deutschland GmbH**  
Dental Products  
Carl-Schurz-Strasse 1  
D-41453 Neuss, Germany

