

3M Science.
Applied to Life.™

Esthetic Restoration Solution

**A beautiful smile
she can't wait to share.**



Because you want every patient to leave with a beautiful smile.

Patients love the way highly polished restorations look and feel. But creating a natural-looking gloss can take more of your time. So what do you do? Sacrifice esthetics for productivity? Or spend the extra time creating a beautiful smile?

The days of time-consuming esthetic restorations are gone. With our Esthetic Restoration Solution, you'll love how simple it is to create natural-looking restorations which are polished to a high gloss in no time.

Don't be surprised if word of your great work gets around... because a beautiful smile is worth sharing.



The truth is, the composite you use matters. You may spend a lot of time polishing composite restorations, only to find the initial gloss doesn't last long.



Class I

Or, you can create smiles your patients will love by using our Esthetic Restoration Solution. Thanks to 3M's nanotechnology, Filtek™ Supreme XTE Restorative is easy to polish and offers excellent polish retention.




Class II



Class III

The nanotechnology also makes Filtek Supreme XTE restorations strong, giving you the option of using it in any place or for any Class of restoration.



Class IV



Class V

Nanotechnology. For beautiful and strong restorations ... anywhere.

Diamond particles. For a gorgeous, paste-like gloss.

The esthetic appearance of a restoration depends on the shade, shape, texture and the gloss achieved by final polishing.

With the 3M™ ESPE™ Sof-Lex™ Diamond Polishing System, you can get a gorgeous, diamond paste-like gloss—with the convenience of a rubberized system!

The Sof-Lex™ Pre-Polishing Spiral (beige) smooths and removes scratches in restorations that develop during contouring and prepares the surface for high-gloss polishing. The Sof-Lex™ Diamond Polishing Spiral (pink) delivers the final high-gloss polish.

You'll find it reassuring to know that you no longer need to stop and switch between traditional points, cups, discs and brushes. The flexible spirals adapt to all tooth surfaces—anterior and posterior, convex and concave. They can be sterilized and reused, making them an economical choice for your practice.

A difference you can see!



Filtek™ Supreme XTE Universal Restorative polished with Sof-Lex™ Diamond Polishing System (left) vs a competitor composite polished with their finishing and polishing paste system. Notice a clearer reflection with the Sof-Lex™ Diamond Polishing System.

Now simple esthetics can be as easy as 1, 2, 3 ...

A factor in the success of this solution is 3M's unprecedented history in abrasive technologies. We leveraged that expertise to create a polishing solution that can help you maximize the esthetic benefits of Filtek™ Supreme XTE Universal Restorative with ease and simplicity.

Filtek™ Supreme XTE
Universal Restorative

1



Sof-Lex™
Pre-Polishing Spiral

2



Sof-Lex™ Diamond
Polishing Spiral

3



Ordering Information

Item No.	Product Information
	Filtek™ Supreme XTE Universal Restorative
4911S	Filtek Supreme XTE Universal Restorative Single Shade Introductory Kit – Capsule Kit includes: 40 – 0.2g capsules: 20 each of Body shades A2 and A3; Technique Guide; Instructions for Use
4910S	Filtek Supreme XTE Universal Restorative Single Shade Introductory Kit – Syringe Kit includes: 4 – 3g syringes: 2 each of Body shades A2 and A3; Technique Guide; Instructions for Use
4911P	Filtek Supreme XTE Universal Restorative Professional Kit – Capsule Kit includes: 160 – 0.2g capsules in plastic reusable case: 10 each of Dentine shades: A2, A3, A4 20 each of Body shades: A2, A3, A3.5, A4 10 each of Enamel shades: A1, A2, A3, W 10 each of Translucent shades: AT 1 Shade Selector Wheel; Technique Guide; Instructions for Use
4910P	Filtek Supreme XTE Universal Restorative Professional Kit – Syringe Kit includes: 12 – 3g syringes each of shades: Dentine shades: A2, A3, A4 Body shades: A2, A3, A3.5, A4 Enamel shades: A1, A2, A3, W Translucent shades: AT 1 Shade Selector Wheel; Technique Guide; Instructions for Use



Simple Esthetics

While Filtek™ Supreme XTE Universal Restorative is known for its wide range of shades and opacities, most single-shade restorations are made with one of 8 Body shades.

A1B A2B A3B A3.5B
B1B B2B C1B C2B

Filtek Supreme XTE Refills

Dentine shades	Syringe	Capsule
A1D	4910A1D	4911A1D
A2D	4910A2D	4911A2D
A3D	4910A3D	4911A3D
A4D	4910A4D	4911A4D
B3D		4911B3D
C4D	4910C4D	4911C4D
White	4910WD	4911WD
Body shades	Syringe	Capsule
A1B	4910A1B	4911A1B
A2B	4910A2B	4911A2B
A3B	4910A3B	4911A3B
A3.5B	4910A3.5B	4911A3.5B
A4B	4910A4B	4911A4B
B1B	4910B1B	4911B1B
B2B		4911B2B
B3B	4910B3B	4911B3B
C1B	4910C1B	4911C1B
C2B	4910C2B	4911C2B
C3B	4910C3B	4911C3B
D2B	4910D2B	4911D2B
White	4910WB	4911WB
Extra White	4910XWB	4911XWB
Enamel shades	Syringe	Capsule
A1E	4910A1E	4911A1E
A2E	4910A2E	4911A2E
A3E	4910A3E	4911A3E
B1E	4910B1E	4911B1E
B2E	4910B2E	4911B2E
D2E	4910D2E	4911D2E
White	4910WE	4911WE
Extra White	4910XWE	4911XWE
Translucent shades	Syringe	Capsule
Clear	4910CT	4911CT
Blue	4910BT	4911BT
Amber	4910AT	

Syringe shade refills include 1 – 3g syringe; Capsule shade refills: Dentine, Enamel and Translucent Shade Refills include 10 – 0.2g capsules; Body Shade Refills include 20 – 0.2g capsules; Instructions for Use.



Sof-Lex™ Diamond Polishing System

Item #	Product Information
5090	Sof-Lex™ Pre-Polishing Spirals—Refill Includes: 15 Sof-Lex Pre-Polishing Spirals (Beige)
5091	Sof-Lex™ Diamond Polishing Spirals—Refill Includes: 15 Sof-Lex Diamond Polishing Spirals (Pink)



3M South Africa (Pty) Ltd.
146a Kelvin Drive
Woodmead, 2128
3m.espe.za@mmm.com

www.3M.com/Dental

www.3M.com/Dental

3M, ESPE, Filtek and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH.
© 3M 2016. All rights reserved.