

3M Science.
Applied to Life.™

Raising the bar with stability and control.



3M™ Incognito™ Brackets with Tip Bar
Invisible, intelligent, individualised.

Incognito™
Appliance System

Introduction of the **3M™ Incognito™ Brackets with Tip Bar** continues to enhance lingual treatment, delivering targeted results through a new feature that improves angulation control. Extensions on the mesial and distal sides of the lingual bracket elongate the bracket slot to increase the stability of the interplay between the brackets and archwires.

*"I evaluated the new **Incognito Brackets with Tip Bar**, with mesial and distal extensions to the slot on the lower canines and upper 3-to-3. I found that the slot extension helps to **control and correct the angulation** of those teeth in the early stage of development. This product improvement will be essential to my future cases, especially considering the **reduction in chair and treatment time.**"*

Dr. Esfandiar Modjahedpour

Latest improvement in a series of ongoing innovation.

In conjunction with the existing form and function, the newly developed tip bar feature helps the doctor by supporting more precise finishing with increased control through added stability.

Lab testing showed that the Incognito Bracket with Tip Bar design for anterior brackets, used in combination with a 3M™ AlastiK™ Lingual Ligature, resulted in 75% more angulation movement compared to a conventional Incognito anterior bracket in combination with a powerchain overtie.¹

Targeted results.

With the addition of the Incognito Brackets with Tip Bar to the Incognito Appliance System, you accentuate the ability to control angulation.

Simply put, the high-precision bracket slots mean more controlled tooth movement. This new bracket feature is optionally available for upper arch 3-to-3 and lower arch canines, where tip control is most important. Incognito Brackets with Tip Bar are also compatible with the 3M™ Incognito™ Clear Precision Tray.




Incognito Brackets with Tip Bar



Traditional Incognito Bracket

Keeping up with the Incognito System.



Features and benefits.

Setting a new standard for aesthetic, comfortable and individually customised orthodontic options, the Incognito Bracket with Tip Bar design ...

- Enables utilisation of the Incognito Bracket with Tip Bar for as many teeth as possible, even for crowded cases.¹
- Offers same or better patient comfort,¹ ease of use,¹ and handling¹ as the conventional Incognito Bracket.
- Available on lower canines and upper 3-to-3.
- Provides even more angulation movement when combined with an AlastiK Lingual Ligature.¹

REFERENCES

1. Data on File

How to order Incognito Brackets with Tip Bar

Orders can be placed within 3M™ Unitek™ Treatment Management Portal | TMP.



3M Oral Care
A Division of 3M Australia Pty Limited
Building A, 1 Rivett Road,
North Ryde NSW 2113 Australia
Phone: 1300 363 484
3M.com.au/Oralcare

3M New Zealand Limited
94 Apollo Drive, Rosedale
Auckland 0632 New Zealand
Phone: 0800 441 622
3M.co.nz/Oralcare

Raise the bar today!
Visit 3M.com.au/Incognito or
contact your 3M Representative.

3M, AlastiK, Forsus, Incognito, and Unitek are trademarks of 3M. © 3M 2017. All rights reserved.
70-2013-0837-9 1701