


3M™ Bair Hugger™ Warming Blankets and Gowns

Intraoperative Blankets


XL Upper Body
Model 523




Lower Body
Model 525




Torso
Model 540



Dual Port Torso
Model 542




Multi-Position
Upper Body
Model 622



Specialty and Underbody Blankets

Adult Underbody
Model 545



Cath Lab
Model 560




Surgical Access
Model 570




Spinal Underbody
Model 575




Lithotomy Underbody
Model 585




Full Body Surgical
Model 610



Sterile Cardiac
Model 630



Full Access
Underbody
Model 635



Sterile Full Access
Underbody
Model 637



Sterile Cardiac
Access
Model 645



Pediatric Blankets and Gown

Pediatric Long
Model 530



Small Lower Body
Model 537



Pediatric Underbody
Model 555



Large Pediatric Underbody
Model 550



Pediatric
Model 310




Pediatric Warming Gown
Model 81501



Pre- and Post-Operative Blankets

Outpatient w/booties
Model 110



Outpatient
Model 111




Full Body
Model 300



Chest Access
Model 305



Multi-Access
Model 315



Perioperative Gowns

OR Warming Gown
Model 8XX01



Plus Warming Gown
Model 8XX02



Flex Warming Gown
Model 8XX03



3M™ Bair Hugger™ Warming Unit, Model 77500



- Four temperature settings allow delivery of effective forced-air patient warming
- Three temperature sensors and adjustable airflow provide safe, quick temperature response
- Tracks and reports information on forced-air warming duration, over-temperature conditions and fault codes, all accessible via the front panel
- Accessories for the Model 775 warming unit allow for mounting on a hard surface, an I.V. pole, bed rail or rolling cart

3M™ Bair Hugger™ Patient Adjustable Warming Unit, Model 87500

- Offers the ability to prewarm patients, assisting in the prevention of unintended hypothermia
- Personalized comfort as patients control the unit's airflow and temperature (ambient to 43°C)
- Easy "snap-fit" flexibility in mounting the unit on the wall, a bedrail or an IV pole
- Convenient, built-in hose storage

