

# 3M™ ESPE™ Clinpro™ Sealant

Introducing colour-change technology.




**COLOUR YOU CAN'T MISS**

**3M ESPE**

# 3M™ ESPE™ Clinpro™ Sealant

is the most exciting advancement in sealant technology in decades  
with an improvement that's hard to overlook.

YOU CAN SEE THE DIFFERENCE IN OUR  
COLOUR-CHANGE TECHNOLOGY.



Clinpro sealant goes on pink for easy-to-see application, and cures to a natural white. Applied with our convenient ultra-fine syringe tip, its low viscosity formula flows easily into pits and fissures.


WE ASKED DENTAL PROFESSIONALS WHAT  
WOULD MAKE THE IDEAL SEALANT

We listened and developed our new Clinpro sealant to meet their needs:

- Convenient direct delivery syringe for easy application
- Colour-change technology to visualize placement
- Low viscosity to flow easily into pits and fissures
- Contains and releases fluoride
- Bonds to enamel

DENTAL PROFESSIONALS PREFER IT OVER  
OTHER SEALANTS.

*77% of respondents rated the overall performance of Clinpro sealant to be superior to their current sealant.*



Source: 3M Dental Sealant Evaluation Summary Report. The Market Solutions Group, Inc., February 2001.



### 3M™ ESPE™ Clinpro™ Sealant Introductory Kit, Syringes:

- 12626 (2) 1.2 ml Clinpro Sealant Syringes  
(1) 3 ml 3M™ ESPE™ Scotchbond™ Etchant Syringe  
(20) Sealant Syringe Tips (Black)  
(24) Etchant Syringe Tips (Blue)

### 3M™ ESPE™ Clinpro™ Sealant Refill, Syringe

- 12627 (1) 1.2 ml Clinpro Sealant Syringe  
(10) Sealant Syringe Tips (Black)

**3M ESPE**  
Dental Products  
3M Canada Company  
Post Office Box 5757  
London, Ontario N6A 4T1

© 2010, 3M. All rights reserved.  
3M, Clinpro and Scotchbond are trademarks of 3M.  
ESPE is a trademark of 3M ESPE AG.  
Used under license in Canada.

3M ESPE Customer Care: 1-888-363-3685  
3M ESPE Website: [www.3MESPE.ca](http://www.3MESPE.ca)