

Dressing change considerations

- ▶ Observe the dressing frequently
- ▶ As the dressing absorbs, exudate will wick to the top of the dressing and discolouration may be noticeable
- ▶ When the exudate spreads to the edges of the dressing or the dressing begins to leak, a dressing change is indicated
- ▶ Dressing change frequency will depend on the type of wound, volume of exudate and clinical situation
- ▶ Change at least every 7 days, or as indicated by local protocol

3M™ Tegaderm™ Foam Adhesive Dressing



Ok



Observe



Change

3M™ Tegaderm™ Foam Dressing



Ok



Observe



Change