



Science.
Applied to Life.™



3M™ Tegaderm™ Foam Adhesive Dressing Mini Sizes

Smart. Adapts to a change in exudate levels.
Simple. To apply and remove.
Small. To fit those difficult to dress wounds.

3M™ Tegaderm™ Foam Adhesive Dressing

Mini Sizes

The **mini wrap** is a custom fitted dressing for easy application on the fingers and toes.

► **Unique**

There are no other dressings of this type currently available in the market.

► **Safe**

An aseptic application method and a bacteria proof covering over the foam. Wraps quickly and easily around the digit without the need for cutting.

► **Economic**

Dressing can remain in place for up to 7 days.

► **Comfortable**

It is low profile with no seams. It conforms well to digits allowing normal function and movement.

► **Versatile**

These dressings can be used anywhere on the body and have proved useful for dressing genitalia, ear pressure ulcers as well as fingers and toes.



90615



**Smart.
Simple.
Small.**



The **mini oval** is a miniature version of the Tegaderm foam dressing.

► **Smart**

Can be used on low to high exuding wounds, reducing the need for multiple dressings.

► **Safe**

An aseptic application method.

► **Economic**

Dressings can remain in place for up to 7 days.

► **Comfortable**

It has a low profile and conforms well to the body. The extended Tegaderm border is low profile and reduces the risk of rucking.

► **Gentle**

It is an island dressing with soft foam that protects the peri-wound area.

► **Practical**


The oval shape reduces the possibility of early dressing failure.

► **Versatile**

The mini oval has a multitude of uses: skin tags, mole removal, minor operations, bone assisted hearing aid wounds, pressure ulcers on ears and chin, maxillofacial surgery, malleoli, toe surgery, rhinoplasty to name but a few...

► **Four layers**

- 4 The Tegaderm layer is bacteria proof and viral proof
- 3 The moisture control layer that enables the dressing to evaporate more fluid when the wound is exuding heavily
- 2 The absorbent layer to manage extra fluid
- 1 The foam layer for quick wicking of exudate, preventing maceration



3M™ Tegaderm™ Foam Adhesive Dressing Mini Sizes

Ordering information

Dressing	3M code	Size	Dressing/box	NHS code	PIP code
	90614	6.9 x 7.6 cm mini oval	10	ELA492	360-9906
	90615	6.9 x 6.9 cm mini wrap	10	ELA493	360-9948



3M United Kingdom PLC
3M House
Morley Street
Loughborough
Leicestershire LE11 1EP
+44 (0)1509 611 611

3M Ireland Limited
The Iveagh Building,
The Park, Carrickmines,
Dublin D18 X015
Ireland
+353 (0)1 280 3555

www.3M.co.uk/healthcare

GF325
J327676

3M, Tegaderm, Coban, Medipore
and Micropore are trademarks of the
3M company. © 3M 2016.
All rights reserved.

3M Science.
Applied to Life.™