



Attest™

Sterilization Solutions

Upgrade to Dual Performance

Introducing the New 3M™ Attest™ Dual-reader BI System
with Any-Well Technology for Steam and VH2O2

STEAM

VH2O2





The Power of 2-in-1

Improve the efficiency of how you process BIs with 3M's latest innovation, so you can do more with less. Introducing our 3M™ Attest™ Auto-reader technology that incubates and delivers BI results for both steam and vaporized hydrogen peroxide sterilization in one reader.

Any-Well Technology


The new 24-minute Dual System works with Any-Well technology, which means Attest™ 1491 or 1492V (for steam) or 1295 (for vaporized hydrogen peroxide sterilization — VH₂O₂) BIs can be incubated in any well of the reader at the same time.

An Upgrade that Gives You More

Get the peace of mind and confidence that comes from knowing every patient receives safe equipment with a dual Auto-reader system that's:

Simple

Current customers can have their existing Attest™ Auto-reader model 490 (steam) or 490H (VH2O2) upgraded to start incubating both steam and VH2O2 biological indicators in one unit.




Consistent

3M offers 24-minute biological indicator results for both steam and vaporized hydrogen peroxide (VH2O2)—providing the ease of consistency.

24 minute BI results*

Familiar & Cost-Sensitive

Use the same trusted Attest™ systems with no change to 3M BIs or test packs. The upgrade is free, with no new hardware to buy.



Connected

Use the same 3M™ Web Application you use now for real-time BI status and results and for digital record keeping.

Your reader will still connect to leading Instrument Tracking Software Systems.



To get started with an upgrade of your current 3M™ Attest™ 490 or 490H Auto-reader, contact your 3M Representative, or sign up at go.3M.com/AttestUpgrade

Current Attest™ Auto-reader 390 users or new 3M customers can contact a 3M sales rep to purchase a new model 490 Auto-reader pre-configured with dual functionality.

Steam Sterilization Monitoring, Now Optimized


3M™ Attest™ Super Rapid Readout Biological Indicators use the same technology you've trusted for years, now optimized to deliver an even faster result.

ISO and FDA Compliant

Attest™ Super Rapid Readout BIs meet the performance requirements of ISO 11138-1:2017 and ISO 11138-3:2017 and satisfy FDA requirements for biological indicators.

Trusted Steam Sterilization Monitoring Portfolio


- Get the performance of the same trusted technology you are familiar with by using 3M™ Attest™ Super Rapid Readout process challenge devices (PCDs).
- The 3M™ Attest™ system is the #1 utilized steam sterilization test pack technology in the U.S., and is now available with the speed of 24-minute results.
- 3M™ Comply™ Chemical indicators, Bowie-Dick test packs and indicator tape provide the core four steam sterilization monitoring products.



*When used with 3M™ Attest™ Auto-reader 490 or the Attest™ Auto-reader 490H with software version 4.2.7 or greater.

Know. Then Go.

Now Every Load Monitoring is not just possible, it's practical. Take control with 24-minute results and get the peace of mind that comes from knowing that every patient will receive a high standard of care.



Results You Can Count on for Vaporized Hydrogen Peroxide Sterilization

The 3M™ Attest™ system for VH2O2 has the broadest range of FDA clearances in the industry for VH2O2 sterilizers and cycle types.** This provides a truly independent monitor for hydrogen peroxide sterilization assurance needs. The 1295 Biological Indicator system is FDA cleared for the following:

STERRAD® Vaporized Hydrogen Peroxide Systems

For use with STERRAD® 100S, STERRAD® 100NX (Standard, Flex, Express and Duo cycles), STERRAD® NX (Standard and Advanced cycles), STERRAD® 100NX with ALLClear™ Technology (Standard, Flex, Express and Duo cycles), STERRAD® NX with ALLClear™ Technology (Standard and Advanced cycles).

Amsco® V-PRO® Low Temperature Sterilization System

For use with V-PRO® maX (Lumen, Non Lumen, and Flexible cycles), V-PRO® 1 plus (Lumen and Non Lumen cycles), V-PRO® 1 (Lumen cycles), V-PRO® 60 (Lumen, Non Lumen, and Flexible cycles).

3M Provides a Complete VH2O2 System for Monitoring

- Use the key three monitoring products from 3M— BIs, CIs, and indicator tape—to leverage the full system of sterilization assurance products.
- A consistent blue-towards-pink color change for all VH2O2 indicator products including the CI on the biological indicator, internal process CI and external indicator tape, makes it easy to read and interpret results for all of your VH2O2 monitoring needs.

**As of 5/1/2018



The 3M™ Attest™ Brand Advantage



The 3M™ Attest™ brand advantage gives you a standardized, simplified and supported system for sterilization across steam, VH2O2 and EO modalities.

Fight SSIs from Every Angle

Our 3M portfolio of sterilization assurance products and sterilization monitoring expertise is your first line of defense in your fight to help protect patients from surgical site infections.



Load Control
Pack Control
Exposure Control
Equipment Control

Product	Description	Packaging	Shipping
490/490H	3M™ Attest™ Auto-reader 490/490H	1 unit/box	1 box/case
41482V	3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack	24 each/box	1 box/case
1491	3M™ Attest™ Super Rapid Readout Biological Indicator	50 each/box	4 boxes/case
1243	3M™ Comply™ (SteriGage™) Chemical Integrators	500 each/bag	2 bags/case
		100 each/bag	10 bags/case
1355	3M™ Comply™ Indicator Tapes	18mm roll	28 rolls/case
		24mm roll	20 rolls/case
00135LF	3M™ Comply™ Bowie-Dick Plus Test Pack	20 packs/box	4 boxes/case



Load Control
Pack Control
Exposure Control

Product	Description	Packaging	Shipping
490/490H	3M™ Attest™ Auto-reader 490/490H	1 unit/box	1 box/case
1295	3M™ Attest™ Rapid Readout Biological Indicator	30 each/bag	4 bags/case
1248	3M™ Comply™ Chemical Indicator Strip	250 each/bag	4 bags/case
1228	3M™ Comply™ Indicator Tape	1,9cm roll	24 rolls/case

To get started with an upgrade of your current 3M™ Attest™ 490 or 490H Auto-reader, contact your 3M Representative, or sign up at go.3M.com/AttestUpgrade.

For more information or support visit go.3M.com/3MAdvantage or call the 3M Health Care Help Line at 1-800-228-3957.



3M Infection Prevention Division
3M Health Care
3M Center, 2510 Conway Ave
St. Paul, MN 55144-1000
U.S.A.
1-800-228-3957
www.3M.com/medical

3M, Attest, Comply and SteriGage are trademarks of 3M. All other trademarks are property of their respective owners.

Please recycle. Printed in the U.S.A.
© 3M 2018. All rights reserved.
70-2011-7107-4