



Clarity™

Esthetic Orthodontic Solutions

Parts List

3M™ Clarity™ Ultra Self-Ligating Brackets



3M™ MBT™ Appliance System Prescription

MAXILLARY			IN/OUT		Part Number – .022 in.	
Tooth	Torque	Angulation	in.	mm	L	R
Central	+17°	4°	0.026	0.66	003-301	003-302
Lateral	+10°	8°	0.035	0.89	003-303	003-304
Cuspid Hk	0°	8°	0.021	0.53	003-307	003-308
Bicuspid Hk	-7°	0°	0.029	0.74	003-310	003-311

MANDIBULAR			IN/OUT		Part Number – .022 in.	
Tooth	Torque	Angulation	in.	mm	L	R
Anterior	-6°	0°	0.045	1.14	003-350	
Cuspid Hk	0°	3°	0.020	0.51	003-353	003-354
1st Bicuspid Hk	-12°	2°	0.030	0.76	003-357	003-358
2nd Bicuspid Hk	-17°	2°	0.034	0.86	003-361	003-362

Add Prefix "4" for 3M™ APC™ Flash-Free Adhesive System
Example: 4003-301

Patient Kits

	Part Number – .022 in.
U/L Lt/Rt 5x5 Hk	003-111
U Lt/Rt 5x5 Hk	003-116

Add Prefix "4" for APC Flash-Free Adhesive System

Instruments for the Ultra SL System

	Part Number
3M™ Unitek™ Debonding Instrument	804-175
3M™ Unitek™ Open-Close Instrument	804-155

For more information ...
 Visit 3M.com/ClarityUltra
 or contact your
 3M Oral Care
 representative.




REF 804-175



REF 804-155

3M Oral Care
 2510 Conway Avenue
 St. Paul, MN 55144-1000 USA
 Phone 1-800-423-4588
 Web 3M.com/ortho

3M Canada
 Health Care Division
 300 Tartan Dr.
 London, ON N5V 4M9 Canada
 Phone 1-800-443-1661

3M, APC, Clarity, MBT, and Unitek are trademarks of 3M.
 Used under license in Canada.
 © 3M 2018. All rights reserved.
 70-2013-0963-3 1807