



3M post and core solution. Complexity simplified.

For predictable results – in less time: The endodontic procedure is very complex. So the simpler and more reliable the endodontic restorative procedure is, the better – for both the dentist and the patient. With fewer clinical steps and fewer products, the 3M post and core solution lets you work faster, easier and more efficiently.

It helps you deliver predictable, robust outcomes. All you need is the new RelyX 3D glass fiber posts plus three other easy-to-use 3M products: RelyX Unicem 2 Automix cement, Scotchbond universal adhesive, and Filtek One bulk fill restorative. These four products have one common objective: Simplify the complex!



3M™ Filtek™ One
Bulk Fill Restorative




3M™ Scotchbond™
Universal
Adhesive



3M™ RelyX™
Fiber Post 3D
Glass Fiber Post



3M™ RelyX™ Unicem 2 Automix
Self-Adhesive Universal
Resin Cement



3M™ RelyX™ Unicem 2 Automix Self-Adhesive Universal Resin Cement offers excellent bond strength plus saves you precious time, as there is no need to etch, prime or bond. The automix syringe, together with the endo tip, lets you apply the ready-mixed cement directly in the root canal from bottom to top. In one application – without voids.

3M™ RelyX™ Fiber Post 3D Glass Fiber Post's coronal 3D macro retention offer higher mechanical retention to the core build-up material. Due to its microporous surface it also provides a strong bond in the root canal without a silanating step. Post position control is greatly improved as a result of its enhanced radiopacity.

3M™ Scotchbond™ Universal Adhesive provides a reliable bond of the core build-up material to the post and tooth structure – with no need for an additional primer. One application – done!

3M™ Filtek™ One Bulk Fill Restorative is ideally suited as core build-up material, especially for big lesions in combination with a post: The material is sculptable but also easy to adapt to both tooth structure and post. For easier preparation, the cured material cuts like dentin. Core build-ups can be done faster and easier.



Ordering information

Item no. Product information

3M™ RelyX™ Fiber Post 3D Glass Fiber Post Kits

56957	Introductory kit: 20 posts sizes 0, 1, 2, 3 (5 posts of each size), 5 drills sizes 0, 1, 2, 3, Universal (1 drill of each size)
56958	Trial Kit with RelyX™ Unicem 2 Automix Cement: 10 posts sizes 1, 2 (5 posts of each size), 3 drills sizes 1, 2, Universal (1 drill of each size), 1 RelyX™ Unicem 2 Automix Syringe A2 – 8.5g, 10 mixing tips, 10 endo tips
56960	Post and core build-up kit with RelyX™ Unicem 2 Automix Cement: 1 Trial Kit RelyX™ Fiber Post 3D/RelyX™ Unicem 2 Automix Cement, 20 capsules Filtek™ Bulk Fill Posterior Restorative A3 – 0.2g each, 1 vial Scotchbond™ Universal Adhesive – 5 ml

3M™ RelyX™ Fiber Post 3D Glass Fiber Posts and Drills

56948	10 posts, size 0 (1.1 mm dia), white
56949	10 posts, size 1 (1.3 mm dia), yellow
56950	10 posts, size 2 (1.6 mm dia), red
56951	10 posts, size 3 (1.9 mm dia), blue
56873	1 drill, size 0 (1.1 mm dia), white
56864	1 drill, size 1 (1.3 mm dia), yellow
56865	1 drill, size 2 (1.6 mm dia), red
56866	1 drill, size 3 (1.9 mm dia), blue
56867	1 Universal drill

3M™ RelyX™ Unicem 2 Automix Self-Adhesive Resin Cement

	Refills – each containing 1 syringe – 8.5g, 10 mixing tips (regular), 5 mixing tips (wide), 5 endo tips
56846	Shade A2 Universal
56847	Shade A3 Opaque
56848	Shade translucent
	Value packs – each containing 3 syringes – 8.5g, 30 mixing tips (regular), 20 mixing tips (wide), 15 endo tips, 5 intra-oral tips
56849	Shade A2 Universal
56857	Shade A3 Opaque
56858	Shade Translucent

Item no. Product information

Accessories – Cement mixing tips

56851	30 mixing tips (regular)
56919	15 mixing tips (wide)
56918	15 mixing tips (wide), 15 intra-oral tips

3M™ Scotchbond™ Universal Adhesive

41254	Introductory Kit: 1 vial Scotchbond™ Universal Adhesive – 5ml, 50 disposable mixing wells, 50 disposable microbrushes, 1 etchant syringe – 3ml, 25 dispensing tips for etchant syringe
41258	Refill: 1 vial Scotchbond™ Universal Adhesive – 5 ml

3M™ Filtek™ One Bulk Fill Restorative

Refill	Item no.	Item no.
Shade	Syringe incl. 1 syringe – 4g	Capsule incl. 20 capsules – 0.2g each
A1	4866A1	4867A1
A2	4866A2	4867A2
A3	4866A3	4867A3
B1	4866B1	4867B1
C2	4866C2	4867C2

3M Oral Care
3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough
Leicestershire LE11 5RB
www.3M.co.uk

www.3M.com/RelyX3D

GD023.
J382343.
Ref: The Dental Advisor Vol. 34, No. 01 Jan/Feb 2017.
3M, ESPE, Filtek, RelyX and Scotchbond are trademarks of the 3M company. © 3M 2018. All rights reserved.

