
# Science. Applied to Life.™

## **3M™ Lava™ Esthetic**

Fluorescent Full-Contour Zirconia

# Natural esthetics.

Whatever the case.





Now one disc has it all. The first zirconia with **built-in fluorescence** for all shades.\* Are you still struggling with staining? Get ready for the next level of highly esthetic, full-contour zirconia restorations.

#### Shade match.

Built-in shade gradient, designed for excellent shade match with all VITA classical shades including Bleach.

#### High translucency.

Optimized for esthetic full-contour restorations.

#### **Durability.**

Durable full-contour with 800 MPa through and through. Indicated for crowns, 3-unit bridges, inlays/onlays, veneers.\*\*

### Productivity.

Benefit from a streamlined mill-sinterglaze process. No shading liquids. No drying time.



<sup>\*</sup> Fluorescence determined with light sources simulating natural UV light

<sup>\*\* 3-</sup>point bending strength according to ISO 6872:2015; qualified for Type II, class 4; indications: crowns, bridges with a maximum of one pontic that must be supported on each side by a crown (prosthesis not to exceed three units), inlays, onlays and veneers.

#### Ordering information



#### 3M<sup>™</sup> Lava<sup>™</sup> Esthetic Fluorescent Full-Contour Zirconia

1 disc per pack	Disc height 14 mm	Disc height 18 mm	Disc height 22 mm
Shade	Item No.	Item No.	Item No.
Bleach	69319	69327	69335
A1	69320	69328	69336
A2	69321	69329	69337
A3	69322	69330	69338
A3.5	69323	69331	69339
A4	69352	69361	69370
B1	69624	69332	69340
B2	69353	69362	69371
B3	69354	69363	69372
B4	69355	69364	69373
C1	69325	69333	69341
C2	69356	69365	69374
C3	69357	69366	69375
C4	69358	69367	69376
D2	69326	69334	69342
D3	69359	69368	69377
D4	69360	69369	69378

#### Reliable cementation made easy.

3M recommends using 3M™ RelyX™ Unicem 2 Self-Adhesive Resin Cement combining high bond strength with simple handling.



3M™ RelyX™ Unicem 2 Automix Cement		
Shade		
A2 Universal		
A3 Opaque		
Translucent		

Each refill contains 1 syringe – 8.5 g, 15 mixing tips (regular), 5 mixing tips (wide), 5 endo tips

Ask your 3M sales representative for your copy of the

3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia Crown Quick Start Guide

or download from 3M.com/LavaEsthetic



### **3**M

**3M Oral Care** 2510 Conway Avenue St. Paul, MN 55144-1000 USA 1-800-634-2249

### 3M.com/Dental

**3M Canada** Post Office Box 5757 London, Ontario N6A 4T1 Canada 1-888-363-3685

3M, Lava and RelyX are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies.

© 3M 2018. All rights reserved.