

**3M** Science.  
Applied to Life.™



**3M™ Lava™ Esthetic**  
Fluorescent Full-Contour Zirconia

Available from  
your lab.  
**In all VITA  
classical  
shades.**

# Natural esthetics. Whatever the case.



Ask your lab for 3M™ Lava™ Esthetic Zirconia – the high-quality full-contour zirconia from 3M. Featuring a one-of-a-kind shading technology, Lava Esthetic zirconia is the first zirconia with built-in fluorescence for all shades\* that combines the esthetics of glass ceramics with high full-contour strength.

## Shade match.

Built-in shade gradient, designed for excellent shade match with all VITA classical shades including Bleach.

## High translucency.

Optimized for esthetic full-contour restorations.

## Durability.

Durable full-contour with 800 MPa through and through. Indicated for crowns, 3-unit bridges, inlays/onlays, veneers.\*\*

## From a brand you can trust.

With your Lava Esthetic Zirconia crown or bridge you receive a digital high-precision restoration from 3M – made by your lab.



Before



After



Under UV light



Light source Benda NU-6 KL, UV-A setting


\* Fluorescence determined with light sources simulating natural UV light.

\*\* 3-point bending strength according to ISO 6872:2015; qualified for Type II, class 4; indications: crowns, bridges with a maximum of one pontic that must be supported on each side by a crown (prosthesis not to exceed three units), inlays, onlays and veneers.

# Excellent bond strength. Fewer steps.

For final cementation of 3M™ Lava™ Esthetic Zirconia restorations, 3M recommends using 3M™ RelyX™ Unicem 2 Automix Self-Adhesive Resin Cement offering high bond strength and ease of use without compromising reliability and esthetics. And, as there is no need to etch or prime, this easy-to-use cement will save you precious time and work steps.

Please refer to the instructions for use for the proper cementation procedure.



Shear bond strength after 24 hours, 12,000 thermocycles and 90 days storage, respectively. Source: Rosentritt et al., University of Regensburg, Dentistry Journal 2017, 5, 32

## 3M™ RelyX™ Unicem 2 Automix Cement

Refills	Shade
56846	A2 Universal
56847	A3 Opaque
56848	Translucent


Each refill contains 1 syringe – 8.5 g, 15 mixing tips (regular), 5 mixing tips (wide), 5 endo tips



Ask your 3M sales representative for your copy of the

## 3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia Crown Quick Start Guide

or download from [3M.com/LavaEsthetic](http://3M.com/LavaEsthetic)



[3M.com/Dental](http://3M.com/Dental)

**3M Oral Care**  
2510 Conway Avenue  
St. Paul, MN 55144-1000  
USA  
1-800-634-2249

**3M Canada**  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1-888-363-3685

3M, Lava and RelyX are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies.  
© 3M 2018. All rights reserved.