


3M™ Dual Arch Impression Tray

3M™ Dual Arch Impression Tray is a disposable plastic impression tray that can be used for small impression cases (1-2 units). Using 3M impression materials with 3M Dual Arch Impression Tray allows for a robust procedure with easier handling and less waste.

3M Dual Arch Impression Tray can be used with Polyether or VPS impression materials of medium or heavy body tray consistency for:

- Crown/Bridge Arch Impression
- Bite Registration
- Opposing Arch Impression



Ordering information

Item No.	Product Information
71622	3M™ Dual Arch Impression Tray — Posterior Contains: 48 Posterior Trays
71623	3M™ Dual Arch Impression Tray — Anterior Contains: 35 Anterior Trays

3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3M.com/oralcare

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada

Phone 1-888-363-3685

Customer Care Center: 1-800-634-2249
3M.com/greatimpressions

3M is a trademark of 3M. Used under license in Canada. Please recycle.
Printed in U.S.A. © 3M 2018. All rights reserved.