


3M™ Clarity™ Ultra Self-Ligating Brackets

Practise confidently. Smile beautifully.




Create success for patients and your practice.

As you know, the orthodontic business and treatment models are changing. Patients have more choices than ever, and they're playing an active role in decision-making when it comes to the products they purchase. At the same time, emerging technology is driving innovation and new efficiencies in the practice of orthodontics.

In the midst of this change, your goal remains the same—creating beautiful, confident smiles that will last your patients a lifetime.

Showcase your expertise. Deliver precise results.



Doors withstand more than typical treatment use*

Reliable self-ligating doors withstand more than 2 times the opening and closing cycles of typical treatment.

Typical 2-year treatment

2x

3M[™] Clarity[™] Ultra Bracket door reliability

Empowering your efficiency and results.

You and your patients want the same thing: effective treatment options with positive results. To achieve this, you need control and predictability from start to finish. 3M™ Clarity™ Ultra Brackets give you that power—allowing for successful treatment planning to generate precise outcomes.



Easy wire changes without ligation and ligatures




Wide ceramic door that assists rotational control with a broad mesial-distal span



Passive or active on demand



Rotational control without adding accessories



Simple, efficient bonding process.

Powered by technology, the 3M™ APC™ Flash-Free Adhesive Coated Appliance System is the most efficient bonding system in orthodontics. 3M's innovative design creates a unique adhesive bonding substrate on the base of the bracket that saves you time and improves patient satisfaction.



Available precoated with 3M™ APC™ Flash-Free Adhesive for superior bonding confidence.





Saves Time and Cost

- Faster and more efficient bonding than non-precoated brackets
- Eliminates the flash removal step from bonding
- Simple, easy clean-up at debonding

Superior Bonding Efficiency

- Bond time is reduced up to 40% per bracket
- No flash removal reduces accidental bracket movement after positioning, reducing the need to reposition
- Helps protects tooth enamel under adhesive





3M[™] Victory Series[™] Superior Fit Buccal Tubes



3M[™] Clarity[™] Ultra Self-Ligating Brackets



3M[™] Forsus[™] Class II Correctors

One methodology. Optimized practice efficiency.

Powered by 3M science, 3M™ Clarity™ Self-Ligating Brackets and other 3M products – including 3M™ APC™ Flash-Free Adhesive, 3M™ Forsus™ Class II Correctors and the 3M™ MBT™ Versatile+ Appliance System – are designed to work together to give you efficient and custom control to achieve results with a single treatment methodology.



The clear solution that inspires smiles.

Designed to blend in with natural tooth colour, 3M™ Clarity™ Ultra brackets provide a fully esthetic solution. Translucent and stain-resistant, the bracket has no visible metal parts, giving it an attractive appearance.

The absence of elastomeric ligatures that may stain during treatment also enhances the visual appeal. Smooth, rounded, ceramic surfaces and edges are designed to deliver comfort over the course of treatment that both you and your patients will appreciate.



A perfect companion for healthy, beautiful smiles.

Paired with Clarity™ Ultra
Brackets, 3M™ Unitek™ Esthetic
Archwires are a perfect
complement to help your
patients smile confidently.

3M™ Clarity™ Ultra Self-Ligating Brackets

3M™ MBT™ Appliance System Prescription



Maxillary			In/Out		Part Number022 in	
Tooth	Torque	Angulation	in	mm	L	R
Central	+17°	4°	0.026	0.66	003-301	003-302
Lateral	+10°	8°	0.035	0.89	003-303	003-304
Cuspid Hk	0°	8°	0.021	0.53	003-307	003-308
Bicuspid Hk	-7°	0°	0.029	0.74	003-310	003-311

Mandibular			In/Out		Part Number022 in	
Tooth	Torque	Angulation	in	mm	L	R
Anterior	-6°	0°	0.045	1.14	003	-350
Cuspid Hk	0°	3°	0.020	0.51	003-353	003-354
1st Bicuspid Hk	-12°	2°	0.030	0.76	003-357	003-358
2 nd Bicuspid Hk	-17°	2°	0.034	0.86	003-361	003-362
					Add Prefix "4" for 3	M™ APC™ Flash-Free

Adhesive System. Example: 4003-301

Patient Kits

	Part Number022 in
U/L Lt/Rt 5×5 Hk	003-111
U Lt/Rt 5×5 Hk	003-116
	Add Prefix "4" for 3M™ APC™ Flash-Free Adhesive System.

Instruments for the Ultra SL System

	Part Number
3M™ Unitek™ Debonding Instrument	804-175
3M™ Unitek™ Open-Close Instrument	804-155







Practise with clarity.

Our full portfolio of naturally esthetic treatment options are informed by 3M science and fueled by creativity. Innovation drives our leading-edge technology and delivers the choice and control you seek.

Our commitment is to collaborate with orthodontists to provide custom, flexible, esthetic treatment solutions that showcase your expertise and help ensure precise outcomes. Together, we can change lives by creating healthy, beautiful smiles.



3M™ Clarity™ Ultra Self-Ligating Brackets

Visit **3M.ca/Ortho**, or ask your 3M Oral Care representative for more information.





3M Oral Care 3M Canada P.O. Box 5757 London, Ontario N6A 4T1 1 800-364-3577 3M.ca/Ortho

3M Unitek Orthodontic Products 2724 Peck Rd Monrovia, CA 91016 USA 1-800-423-4588

3M, 3M Science. Applied to Life., Clarity, APC, Forsus, MBT, Victory Series and Unitek are trademarks of 3M. Used under license in Canada. © 2018, 3M. All rights reserved. 1807-12514 E BA-18-27129