





















	Order code	Dispensing system	Viscosity low  high	Impression technique	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
Super Quick Setting Materials						
3M™ Impregum™ Penta™ Super Quick (Heavy Body)	69406			1-step	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Medium Body)	69385			Monophase or 1-step	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)	69379			1-step	0:45	2:00
Quick Setting Materials						
3M™ Impregum™ Penta™ Soft Quick (Heavy Body)*	31769			1-step	1:00	3:00
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	31765			Monophase or 1-step	1:00	3:00
3M™ Impregum™ Garant™ Soft Quick (Medium Body)	31788			Monophase or 1-step	1:00	3:00
3M™ Impregum™ Soft Quick (Light Body)	31767			1-step	1:00	3:00
Regular Setting Materials						
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body)	31740			1-step	2:30	3:30
3M™ Impregum™ Penta™ Soft (Medium Body)	31730			Monophase or 1-step	2:45	3:15
3M™ Impregum™ Garant™ Soft (Medium Body)	31785			Monophase or 1-step	1:45	4:15
3M™ Impregum™ Penta™ (Medium Body)	31644			Monophase	2:45	3:15
3M™ Impregum™ Penta™ L DuoSoft™ (Light Body)	31745			1-step	3:15	3:15
3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)	31751			1-step	2:00	3:30

* 3M™ Impregum™ Penta™ Soft Quick (Heavy Body) is not available in New Zealand.

Indications	Impression technique	Single-unit crown	Small implant (≤ 2 units)	Inlay/onlay (≤ 2 units)	Small bridge (≤ 3 units)	Veneers	Bridge (> 3 units)	Implant (> 2 units)	Pick-up impression	Functional impression
Super Quick Setting Materials										
3M™ Impregum™ Penta™ Super Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ Super Quick (Medium or Heavy Body) 3M™ Impregum™ Super Quick (Light Body)	1-step	✓	✓	✓	✓					
Quick Setting Materials										
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	Monophase									
3M™ Impregum™ Penta™ Soft Quick (Medium or Heavy Body) 3M™ Impregum™ Soft Quick (Light Body)	1-step	✓	✓	✓	✓	✓				
Regular Setting Materials										
3M™ Impregum™ Penta™ Soft (Medium Body) 3M™ Impregum™ Penta™ (Medium Body)	Monophase									
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body) 3M™ Impregum™ Penta™ Soft (Medium Body) 3M™ Impregum™ Penta™ or Garant™ L DuoSoft™ (Light Body)	1-step					✓	✓	✓	✓	✓

Always follow the Instructions for Use (IFU) and refer to IFU for full indications, precautions and warnings.

Claims supported by 3M data on file. 3M, Garant, Impregum and Penta are trademarks of 3M. © 3M 2019. All rights reserved.

3M Australia Pty Limited | Building A, 1 Rivett Road, North Ryde, NSW 2113 | 1300 363 454 | 3M.com.au/ImpregumSQ

3M New Zealand Limited | 94 Apollo Drive, Rosedale, Auckland 0632 | 0800 80 81 82 | 3M.co.nz/ImpregumSQ