

A full line of restorative solutions to fit your unique needs.






With a rich legacy of oral care products and services, 3M is one of the world's most trusted companies. Filtek™ restoratives have been used in more than 1 billion restorations worldwide. Our family of Filtek™ dental restoratives lets you choose just the right option to fit any patient situation – from highly aesthetic, to highly simple for an easier day.

Aesthetics	Simplicity	Efficiency
<p>3M™ Filtek™ Supreme XTE Universal Restorative</p> <ul style="list-style-type: none"> • Next level aesthetics • Natural-looking fluorescence • Excellent wear resistance • Exceptional handling • Outstanding strength for anterior and posterior • Wear resistance similar to enamel* • TRUE nanotechnology • 36 shades, 4 opacities • Excellent clinical performance^ 	<p>3M™ Filtek™ Universal Restorative</p> <ul style="list-style-type: none"> • Predictable, repeatable aesthetics • Long-lasting polish • Creamy handling • Low-stress monomer technologies • Excellent wear resistance • TRUE nanotechnology • 8 designer shades plus an Extra White (XW) cover the 19 VITA classical and bleach shades** • Pink Opaquer with excellent masking ability 	<p>3M™ Filtek™ One Bulk Fill Restorative</p> <ul style="list-style-type: none"> • Fast and easy one-step placement up to 5 mm • Low shrinkage stress/stress relief • Excellent adaptation • Excellent handling and sculptability • High radiopacity • Unique tip design for easier access to deep cavities • TRUE nanotechnology • 5 shades 

*In vitro data on file. **3M internal data.

^ 2018 Dental Advisor's 15 Year Clinical Evaluation. Please request your copy from your 3M Representative.

We have you covered for every class of direct restoration.

	Aesthetics	Simplicity	Efficiency
	3M™ Filtek™ Supreme XTE Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ One Bulk Fill Restorative (Bulk up to 5 mm)
Class I 	✓	✓	✓ (up to 4 mm)
Class II 	✓	✓	✓ (up to 5 mm)
Class III 	✓	✓	
Class IV 	✓	✓	
Class V 	✓	✓	
Minimally Invasive	✓	✓	✓
Stress Bearing	✓	✓	✓

✓ = Recommended use

Refer to Instructions for Use for full list of indications.

3M.com.au/Filtek
3M.co.nz/Filtek



3M Oral Care

3M Australia Pty Limited Building A, 1 Rivett Road, North Ryde, NSW 2113 | 1300 363 454 | 3M.com.au
3M New Zealand Limited 94 Apollo Drive, Rosedale, Auckland 0632 | 0800 808 182 | 3M.co.nz



Find us on Facebook

3M and Filtek are trademarks of 3M. Printed in Australia.
© 3M 2019. All rights reserved. AV011477237

