











A full line of restorative solutions to fit your unique needs.

With a rich legacy of oral care products and services, 3M is one of the world's most trusted companies. Filtek™ restoratives have been used in more than 1 billion restorations worldwide. Our family of Filtek™ dental restoratives lets you choose just the right option to fit any patient situation – from highly aesthetic, to highly simple for an easier day.

Aesthetics	Simplicity	Efficiency		
<p>3M™ Filtek™ Supreme XTE Universal Restorative</p> <ul style="list-style-type: none"> • Next level aesthetics • Natural-looking fluorescence • Excellent polish retention • Exceptional handling • Outstanding strength for anterior and posterior • Wear resistance similar to enamel¹ • Can be safely warmed² and injection molded • TRUE nanotechnology • 36 shades, 4 opacities • Excellent clinical performance[^] 	<p>3M™ Filtek™ Supreme Flowable Restorative</p> <ul style="list-style-type: none"> • Virtually no bubbles or material “run-on” • Easy-to-bend cannula for improved access • Ergonomic design is easy to hold and inject • Excellent adaptation, polish retention and wear resistance • Can be safely warmed³ and injection molded • Identical shade match to 3M™ Filtek™ Universal Restorative and 3M™ Filtek™ Supreme XTE Universal Restorative • 12 shades 	<p>3M™ Filtek™ Universal Restorative</p> <ul style="list-style-type: none"> • Predictable, repeatable aesthetics • Long-lasting polish • Creamy handling • Low-stress monomer technologies • Excellent wear resistance • Can be safely warmed² and injection molded • TRUE nanotechnology • 8 designer shades plus an Extra White (XW) cover the 19 VITA classical and bleach shades⁴ • Pink Opaquer with excellent masking ability 	<p>3M™ Filtek™ One Bulk Fill Restorative</p> <ul style="list-style-type: none"> • Fast and easy one-step placement up to 5 mm • Low shrinkage stress/ stress relief • High radiopacity • Excellent adaptation • Excellent handling and sculptability • Unique tip design for easier access to deep cavities • Can be safely warmed¹ and injection molded • TRUE nanotechnology • 5 shades 	<p>3M™ Filtek™ Bulk Fill Flowable Restorative</p> <ul style="list-style-type: none"> • Virtually no bubbles or material “run-on” • Easy-to-bend cannula for improved access • Ergonomic design is easy to hold and inject • Injectable, one-step placement up to 4 mm • Excellent adaptation, polish retention and low shrinkage • Can be safely warmed³ and injection molded • 4 shades 

1. In vitro data on file. 2. Warm capsules to 70°C (158°F) for up to 1 hour. 3. Up to 70°C (158°F) for up to 1 hour, up to 25 times. 4. 3M internal data. ^ 2018 Dental Advisors 15 Year Clinical Evaluation. Please request your copy from your 3M Representative.

We have you covered for every class of direct restoration.

	Aesthetics		Simplicity	Efficiency	
	3M™ Filtek™ Supreme XTE Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Supreme Flowable Restorative (Incremental up to 2 mm)	3M™ Filtek™ Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ One Bulk Fill Restorative (Bulk up to 5 mm)	3M™ Filtek™ Bulk Fill Flowable Restorative (Incremental up to 4 mm)
Class I 	✓	✓ (Liner/Base)	✓	✓ (up to 4 mm)	✓ (Liner/Base)
Class II 	✓	✓ (Liner/Base)	✓	✓ (up to 5 mm)	✓ (Liner/Base)
Class III 	✓	✓	✓		✓
Class IV 	✓		✓		
Class V 	✓	✓	✓		✓
Minimally Invasive	✓	✓	✓	✓	✓
Stress Bearing	✓		✓	✓	

✓ = Recommended use.

Refer to Instructions for Use for full list of indications.

www.3M.com.au/Filtek

3M Australia Pty Limited
Building A, 1 Rivett Road
North Ryde, NSW 2113
1300 363 454
3M.com.au

3M New Zealand Limited
94 Apollo Drive
Rosedale, Auckland 0632
0800 808 182
3M.co.nz

3M and Filtek are trademarks of 3M. © 3M 2020. All rights reserved.
Please recycle. AV011477419