



**Littmann®**  
Stethoscopes

Cardiology

# Built to hear more.

Sturdy and refined, the **3M™ Littmann® Cardiology IV™ Diagnostic Stethoscope** helps isolate subtle sounds\* so you can hear changes in your patients in the most critical and challenging environments. Designed to help you listen for hard-to-hear sounds, including aortic stenosis, S3 gallop murmurs and faint pulmonary anomalies.



**Diagnose and make decisions with confidence.**

- Useful in critical care and challenging environments such as the ED, ICU, Cardiac ICU, step-down unit and other dynamic locations
- Ergonomic and high profile construction; 40% larger chestpiece and 60% deeper bell than 3M™ Littmann® Classic III™ Monitoring Stethoscope
- Designed for use with adult and pediatric patients
- Tunable, dual-sided stainless steel chestpiece with open or closed bell
- Non-chill bell sleeve for greater patient comfort
- Soft-sealing eartips provide optimal comfort and excellent sound occlusion
- Made in the USA without natural rubber latex

## Product Features

<b>Acoustic Performance</b>	9 (1–10 scale)	
<b>Chestpiece Design/Technology</b>	Machined stainless steel, two-sided	
<b>Headset Material</b>	Lightweight, aircraft-grade, aluminum alloy	
<b>Length</b>	22"/56 cm, 27"/69 cm	
<b>Total Weight</b>	175, 185 grams	
<b>Diaphragm Details/Diameter</b>	Two tunable diaphragms 1.7"/4,3 cm, 1.3"/3,3 cm	
<b>Warranty</b>	7 years	

**3M Health Care**  
3M Medical Solutions Division  
2510 Conway Ave., St. Paul, MN 55144 USA  
1-800-228-3957 | Littmann.com

\* Better audibility of high-frequency sounds (adult side, firm pressure) compared to 3M™ Littmann® Cardiology III™ Stethoscope

3M, Littmann, the L Littmann logo, Classic III, Cardiology III and Cardiology IV are trademarks of 3M. Used under license in Canada. Please recycle. Printed in U.S.A. © 3M 2018. All rights reserved. 70-2011-7185-0