



Lab. Order Form

Guide, 2017



3M Incognito
Please submit impressions to:
TOP-Service für Lingualtechnik GmbH
Schleddehauer Straße 81
49152 Bad Essen - Germany
Web: www.incognito.net
Mail: info-incognito@mmm.com
Tel.: +49 (0) 59 71 14 0
Fax: +49 (0) 59 71 14 0

Patient (Surname) _____ Patient (First Name) _____ Date of Bonding _____ To be filled out by TOP-Service

Orthodontic Practice _____ Orthodontist _____ E-Mail _____ Telephone Number _____



Billing Address (Postal Code, City, Country). Please inform us immediately if the shipping address is different.

Please fully complete the treatment plan.

Impression Setup	Digital Impression Digital Setup ¹	Silicone Impression Digital Setup ²	Silicone Impression Manual Setup
Incognito® Clear Precision Tray	☐	☐	☐
Thermoplastic Tray	☐	☐	☐
2-Phase Silicone	☐	☐	☐

Trays

☐ I want to review setup
Please be aware that dispatch date is confirmed after setup approval

☐ Formu® Class II Correctors
2 Lower Arch Attachments

☐ Patient with Ni allergy

☐ Surgery Case

☐ Crossbite Correction

☐ Herbst® Appliance

☐ Wires Recommended by TOP-Service³

Stripping MAX

☐ Essential _____ mm

☐ Yes, if necessary

☐ No stripping

MAX

☐ To be bonded

☐ Setup only
eg. Pasted with other technique or prior on with the Incognito® System

Archwire Sequence MAX

Wires required (please tick)	straight	Indic.	Indic.	Wires required (please tick)	straight	Indic.	Indic.	Wires required (please tick)	Indic.	Indic.
	lat. sec.	lat. sec.	lat. sec.		lat. sec.	lat. sec.	lat. sec.		lat. sec.	lat. sec.
12 Ø				10 × 22				17,5 × 17,5		
14 Ø				10 × 22 ET*				17 × 25		
16 Ø				10 × 24				18,2 × 18,2	E	N
18 × 22	E		N	10 × 24 ET**	E		N	18,2 × 25		
17 × 25				18 × 25						
18 × 25	E		N	18 × 25 (red)						

Archwire Series

ET* Extra torque of 10° on 11/21

ET** Extra torque of 10° on 3-3

red. Laterally reduced

E Recommended wire for no-case

N Recommended wire for non-no-case

Bracket Series / Extras

Ribbowide VH

MAX TIP¹ 3SP TH TL TS

Bracket Series

TIP Tip Be Bracket 3-3

3SP Side plane on 3's

TH Tube with Hook

TL Tube Extra Long

TS Tube Short

T1 Standard

Archwire Sequence MDB

Wires required (please tick)	straight	Indic.	Indic.	Wires required (please tick)	straight	Indic.	Indic.	Wires required (please tick)	Indic.	Indic.
	lat. sec.	lat. sec.	lat. sec.		lat. sec.	lat. sec.	lat. sec.		lat. sec.	lat. sec.
12 Ø				10 × 22				17,5 × 17,5		
14 Ø	E		N	10 × 24	E		N	17 × 25		
16 Ø				18 × 25				18,2 × 18,2	E	N
18 × 22	E		N	18 × 25 (red)				18,2 × 25		
17 × 25										
18 × 25	E		N							

Archwire Series

E Recommended wire for no-case

N Recommended wire for non-no-case

red. Laterally reduced

Bracket Series / Extras

Ribbowide VH

DB TIP¹ TH TL TS

Bracket Series

3-3SL Standard

TIP Tip Be Bracket 3-3

TH Tube with Hook

TL Tube Extra Long

TS Tube Short

T1 Standard

Stripping MDB

☐ Essential _____ mm

☐ Yes, if necessary

☐ No stripping

MDB

☐ To be bonded

☐ Setup only
eg. Pasted with other technique or prior on with the Incognito® System

Remarks

Stamp, Date and Signature

I agree with the General Terms of Sale covered on the back of this document.

¹ This is a recommendation from TOP-Service. It is not a recommendation from the manufacturer of the appliance or the appliance itself.

² See notice on reverse side.

³ In some cases when it is not clinically possible to design an Incognito® Bracket with Tip Box, 3M reserves the right to design a conventional Incognito Bracket.

© 3M 2017. All Rights Reserved. FB-751-60e-Laborschutz Incognito EN_s_01/02/2017

- Fill in all the data
- Allow SIX WEEKS for back delivery
- Sign!!



 Science. Applied to Life.™		Please submit <u>impressions</u> to: TOP-Service für Lingualtechnik GmbH		Schiedehauser Straße 81 49162 Bad Essen · Germany		Web: www.incognito.net Mail: info-incognito@mmm.com		Booking: Mo - Fr 8 am - 5 pm Tel: +49 (0)5472 - 94 91 - 0			
Patient (Surname)			Patient (First Name)			Date of Bonding			To be filled out by TOP-Service		
Orthodontic Practice			Orthodontist			E-Mail			Telephone Number		
Billing Address (Postal Code, City, Country). Please inform us immediately if the shipping address is different.											

Stamp, Date and Signature

I accept the General Terms of Sale covered on the back of this document.

- † This is a recommendation from TOP-Service. The decision on the treatment plan is up to the treating orthodontist.
 ‡ See notice on reverse side.
 § In some cases when it is not clinically possible to design an Incognito™ Bracket with Tip Bar, 3M reserves the right to design a conventional Incognito Bracket.


■ Tray choice:

- CLEAR PRECISION TRAY CPT: The OPTION
- Available from PVS impression or Scanner
- Leads only to DIGITAL SET UP



	Impression Setup	Digital Impression Digital Setup ²	Silicone Impression Digital Setup ²	Silicone Impression Manual Setup
Trays	Incognito™ Clear Precision Tray	☒	☒	
	Thermoplastic Tray	☐	☐	☐
	2-Phase Silicone	☐	☐	☐

The Incognito™ Clear Precision Tray



Inner, transparent silicone tray with digitally defined bracket spaces

Incognito brackets with sandblasted surface in digitally defined position

Hard, transparent outer shell for stable positioning

■ Miscellanea green box

- I ALWAYS want to review the SET UP
- Choose class II correctors, allergies
- Crossbite correction



I want to review setup

Please be aware that dispatch date is confirmed after setup approval



Forsus™ Class II Correctors

2 Lower Arch Attachments



Patient with Ni allergy



Surgery Case



**Wires Recommended
by TOP-Service¹**



Herbst® Appliance



Crossbite Correction

■ Upper dental arch treatment

Stripping MAX

☐ Essential _____ mm


☐ Yes, if necessary

☐ No stripping

MAX

☐ To be bonded

☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System



Please fill in:

B = Bracket | Circle Tooth = Casted Ring | EX = To Be Extracted | F = Fixed
 L = Lower Arch Attachment | Shade In = Occlusal Pad | T = Tube | X = Missing

For an ex-case please indicate space closure or not.

■ Stripping:

- Essential: your choice. (Bolton?)
- Yes, if necessary: **default choice**
- No stripping

Stripping MAX

☐ Essential _____ mm

☒ Yes, if necessary

☐ No stripping

■ Upper dental arch treatment

Stripping MAX

☐ Essential _____ mm


☐ Yes, if necessary

☐ No stripping

MAX

☐ To be bonded

☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System



Please fill in:

B = Bracket | Circle Tooth = Casted Ring | EX = To Be Extracted | F = Fixed
 L = Lower Arch Attachment | Shade In = Occlusal Pad | T = Tube | X = Missing

For an ex-case please indicate space closure or not.

■ Incognito Treatment:

- **To be bonded:** Incognito will be manufactured
- **Set up only:** any other appliance will be used, and incognito lower


MAX

☒ To be bonded

☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System

■ Upper dental arch treatment

■ Appliance description



Bracket Series / Extras

Ribbonwise VH

MAX	TIP ³	3BP	TH	TL	TS
	☐	☐	☐	☐	☐

Bracket Series

TIP	Tip Bar Bracket 3-3
3BP	Bite plane on 3's
TH	Tube with Hook
TL	Tube Extra Long
TS	Tube Short
TI	Standard

Please fill in:

B = Bracket | **Circle Tooth** = Casted Ring | **EX** = To Be Extracted | **F** = Fixed
L = Lower Arch Attachment | **Shade In** = Occlusal Pad | **T** = Tube | **X** = Missing

For an ex-case please indicate space closure or not.

- Appliance design:
 - Brackets, tubes, half occlusal pads.
 - Missing teeth, extractions, IMPLANTS



- Appliance design: Half occlusal pads




Half occlusal pads

- Increase bonding surface in short clinical crowns
- Indications:
 - Second molars: all cases
 - Premolars: adolescents
 - Upper and lower



- Bite planes:
 - Upper canines: classII



Bracket Series / Extras

Ribbonwise VH

MAX	TIP [®]	3BP	TH	TL	TS
	☐	☒	☐	☐	☐

Bracket Series

TIP	Tip Bar Bracket 3-3
3BP	Bite plane on 3's
TH	Tube with Hook
TL	Tube Extra Long
TS	Tube Short
TI	Standard

Incognito Anterior Brackets TIP BAR

- Available for upper 3 to 3 and lower canines




Bracket Series / Extras

Ribbonwise VH

MAX	TIP [®]	3BP	TH	TL	TS
	X				

Bracket Series

TIP	Tip Bar Bracket 3-3
3BP	Bite plane on 3's
TH	Tube with Hook
TL	Tube Extra Long
TS	Tube Short
TI	Standard



Incognito Anterior Brackets TIP BAR

- INDICATIONS:
 - Extraction mechanics
 - Sagittal intermaxillary elastics
 - Upper anterior diastemas
 - Midline shifting
 - DEFAULT OPTION??
 - Crowding interferes with bracket design
 - Also available for LITE



■ Tubes:

- With hook?
- Easy insertion: Standard



Bracket Series / Extras

Ribbonwise VH


MAX	TIP ³	3BP	TH	TL	TS
	☐	☐	☒	☐	☐

Bracket Series

TIP	Tip Bar Bracket 3-3
3BP	Bite plane on 3's
TH	Tube with Hook
TL	Tube Extra Long
TS	Tube Short
TI	Standard



■ Easy Insertion Tubes TI



■ Tubes: combination possibilities



Bracket Series / Extras

Ribbonwise VH

MAX	TIP [®]	3BP	TH	TL	TS
	☐	☐	☒	☒	☐

Bracket Series

TIP	Tip Bar Bracket 3-3
3BP	Bite plane on 3's
TH	Tube with Hook
TL	Tube Extra Long
TS	Tube Short
TI	Standard




■ Lower dental arch treatment



For an ex-case please indicate space closure or not.

Finishing Class	
RIGHT	LEFT
☐	Canine
☐	Molar



Stripping MDB

☐ Essential _____ mm

☐ Yes, if necessary

☐ No stripping

MDB

☐ To be bonded

☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System

■ Finishing class:

- Molar will be defined by canine

Finishing Class	
RIGHT	LEFT
☐	Canine
☐	Molar

■ Lower dental arch treatment


■ Stripping:

- Essential: your choice. (Bolton?)
- Yes, if necessary: **default choice**
- No stripping



For an ex-case please indicate space closure or not.

Finishing Class	
RIGHT	LEFT
☐	Canine ☐
☐	Molar ☐



Stripping MDB

☐ Essential _____ mm

☐ Yes, if necessary

☐ No stripping

MDB

☐ To be bonded

☐ Setup only
e.g. if treated with other
technique or later on with
the Incognito™ System

Stripping MDB

- ☐ Essential _____ mm
- ☒ Yes, if necessary
- ☐ No stripping

■ Lower dental arch treatment

■ Incognito Treatment:

- **To be bonded:** Incognito will be manufactured
- **Set up only:** any other appliance will be used, and incognito Upper



For an ex-case please indicate space closure or not.

Finishing Class	
RIGHT	LEFT
☐	Canine
☐	Molar

Stripping MDB

☐ Essential _____ mm

☐ Yes, if necessary

☐ No stripping

MDB


☐ To be bonded

☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System

MDB

- ☒ To be bonded
- ☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System

- Both arches Incognito
 - Upper & lower set up and Incognito appliance



MAX

- ☒ To be bonded
- ☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System

MDB

- ☒ To be bonded
- ☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System

Upper Incognito, lower labial:

- Upper and lower treatment planning and set up
- Indicate arch form preference: Orthoform III...



MAX

- ☒ To be bonded
- ☐ Setup only
e.g. if treated with other technique or later on with the Incognito™ System

MDB

- ☐ To be bonded
- ☒ Setup only
e.g. if treated with other technique or later on with the Incognito™ System



Upper Incognito, no lower treatment :

- Set up & tray only of the upper arch
- Send antagonist plaster model or impression



MAX



To be bonded



Setup only

e.g. if treated with other technique or later on with the Incognito™ System

MDB



To be bonded




Setup only

e.g. if treated with other technique or later on with the Incognito™ System

■ Lower dental arch treatment

■ Appliance description



Bracket Series / Extras					
Ribbonwise VH					
MAX	TIP ³	3BP	TH	TL	TS
	☐	☐	☐	☐	☐

Bracket Series
 TIP Tip Bar Bracket 3-3
 3BP Bite plane on 3's
 TH Tube with Hook
 TL Tube Extra Long
 TS Tube Short
 TI Standard

Please fill in:

B = Bracket | **Circle Tooth** = Casted Ring | **EX** = To Be Extracted | **F** = Fixed
L = Lower Arch Attachment | **Shade In** = Occlusal Pad | **T** = Tube | **X** = Missing

For an ex-case please indicate space closure or not.

No Extractions

- Individual lateral segments
- SE Ni-Ti for alignment
- SSteel for sliding or elastics??
- Beta III for finishing

5 wires per arch are included in the price

ET* = Extra torque of 15° on 11/21

ET** = Extra torque of 13° on 3-3

red. = laterally reduced

E = recommended wire for Ex-cases

N = recommended wire for Non-Ex-cases

■ = not recommended

Set of archwires MAX			
Wires required (please tick)		straight lat. sec.	indiv. lat. sec.
SE Ni-Ti	12 Ø		
	14 Ø		X ^N
	16 Ø	E	
	16 x 22	E	X ^N
	17 x 25		
	18 x 25	E	X ^N
Steel	16 x 22		
	16 x 22 ET*		
	16 x 24		
	16 x 24 ET**	E	X [?]
	18 x 25		
	18 x 25 (red.)		
Beta III Titanium	17,5 x 17,5		
	17 x 25		
	18,2 x 18,2		E X ^N
	18,2 x 25		

No Extraction: individual lateral sections

Archwire Series

ET* Extra torque of 15° on 11/21

ET** Extra torque of 13° on 3-3

red. Laterally reduced

F Recommended wire for ex-cases

N Recommended wire for non-ex-cases

- Archwire Sequence Design

No Extraction: individual lateral sections



■ Archwire Sequence Design

Extraction: Straight lateral sections



Archwire Sequence MAX								
SE NiTi			STEEL			Beta III Titanium		
Wires required (please tick)	straight lat. sec.	indiv. lat. sec.	Wires required (please tick)	straight lat. sec.	indiv. lat. sec.	Wires required (please tick)	indiv. lat. sec.	
12 Ø	☐	☐	16 × 22	☐	☐	17,5 × 17,5	☐	
14 Ø	☐	N	16 × 22 ET*	☐	☐	17 × 25	☐	
16 Ø	E X	☐	16 × 24	☐	☐	18,2 × 18,2	E X	
16 × 22	E X	N	16 × 24 ET**	E	X N	18,2 × 25	☐	
17 × 25	☐	☐	18 × 25	☐	☐			
18 × 25	E X	N	18 × 25 (red.)	☐	☐			

Archwire Series

ET* Extra torque of 15° on 11/21

ET** Extra torque of 13° on 3-3

red. Laterally reduced

E Recommended wire for ex-cases

N Recommended wire
for non-ex-cases

Posterior first order bend recovery on last BIII Wire

- Archwire Sequence Design

Extraction: Straight lateral sections



Straight lat. seg.

Customized **Lower** Archwire
0.016" x 0.024" Stainless Steel



Customized **Lower** Archwire
0.0182" x 0.0182" Beta III Titanium



- **Appliance design and treatment planning**
 - Sending the models is not enough
 - Treatment planning is our job: EXPLAIN

Remarks

14/24 extractions

Objective: Canine class I, molar class II

Alignment

■ Appliance design and treatment planning

- Compulsory
- Defined by us, not a Top Service job!
- Explain treatment planning

Remarks

Upper incognito treatment

Lower labial treatment:

-MBT prescription

-Orthoform III arch form



Lab. Order Form

Guide, 2017