

	Dispenserings- system	Viskositet lav  høy	Herdehastighet	Maksimal arbeidstid (23 °C) min:sek 	Maksimal intraoral sprøytetid (37 °C) min:sek 	Total intraoral herdetid (37 °C) min:sek 
Skjematerialer						
3M™ Imprint™ 4 Penta™ Putty		Putty	Regular	1:30	–	2:30
3M™ Imprint™ 4 Penta™ Heavy		Heavy Body	Regular	2:00	–	2:00
3M™ Imprint™ 4 Penta™ Super Quick Heavy		Heavy Body	Super Quick	1:15	–	1:15
Sprøytematerialer						
3M™ Imprint™ 4 Super Quick Ultra-Light			Super Quick	1:15	eller 0:35	1:15
3M™ Imprint™ 4 Light			Regular	2:00	eller 1:00	2:00
3M™ Imprint™ 4 Super Quick Light			Super Quick	1:15	eller 0:35	1:15
3M™ Imprint™ 4 Regular			Regular	2:00	eller 1:00	2:00
3M™ Imprint™ 4 Super Quick Regular			Super Quick	1:15	eller 0:35	1:15