

Cement selection made easy.

	RelyX™ Ultimate Adhesive Resin Cement	3M™ RelyX™ U200 Self-Adhesive Resin Cement	RelyX™ Luting 2 Resin Modified Glass Ionomer Cement	RelyX™ Veneer Cement
Metal/Metal Based				
Inlays/Onlays	+	++	+	-
Crowns/Bridges	+	++	++	-
Endodontic Posts	+	++	++	-
Maryland Bridges	++	+	-	-
On Implant Abutments	+	++	++	-
Glass Ceramics (incl. Li. Disilicate) (e.g. e.max®, VITA Mark II, IPS Empress 2)				
Inlays/Onlays/Table Tops	++	+	-	-
Crowns/Bridges	++	++	-	-
Veneers	+	-	-	++
On Implant Abutments	+	+	-	-
Oxide Ceramics (e.g. 3M™ Lava™ Plus High Translucency Zirconia, Brux Zir®, Procera®)				
Inlays/Onlays	++	++	+	-
Crowns/Bridges	+	++	+	-
Endodontic Posts	+	++	+	-
Maryland Bridges	++	+	-	-
On Implant Abutments	+	++	++	-
Resin Nano Ceramics (e.g. 3M™ Lava™ Ultimate CAD/CAM Restorative)				
Inlays/Onlays	++	-	-	-
Veneers	+	-	-	++
Resin Composites				
Inlays/Onlays	++	+	-	-
Crowns	++	++	-	-
Endodontic Posts	+	++	-	-
Veneers	+	-	-	++
On Implant Abutments	++	++	-	-

++ Highly recommended* + Recommended - Not indicated

* Either better performance for this indication or easier handling with equal performance.

Note: The cement selection recommendations in this guide list the key indications for a particular cement, and are intended as a general guide only. 3M RelyX cements are highly versatile products that can be used for multiple indications. Your selection is dependent on individual patient circumstances and your preferred procedural techniques.



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3M.com/dentalcements

3M, Lava and RelyX are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by their respective owners. Please recycle. © 3M 2019. All rights reserved. 70-2013-7175-7