


A complete solution.

3M™ Sample Collection solutions aligned with
3M™ Petrifilm™ Plates

Sample Collection Products



3M™ Swab-Samplers

- Self-contained swab
- Ideal for sampling hard-to-reach areas with long, pre-hydrated 9cm swab
- Offered in multiple volumes for qualitative or quantitative use
- Offered in multiple buffers



3M™ Quick Swabs

- Ready and easy-to-use
- Quickly delivers 1mL sample onto a 3M™ Petrifilm™ Plate
- Labor savings in dilution mixing and pipetting



3M™ Environmental Scrub Sampler with 10mL Wide Spectrum Neutralizer

- Multilayered non-woven composite material maintains integrity when sampling rough areas and scrub dot technology disrupts biofilm and enhances sample collection
- Stick is metal detectable and features an easy-to-use 1-step twist detachment mechanism
- 3M Wide Spectrum Neutralizer is formulated to effectively neutralize common sanitizers and is designed to be compatible with 3M™ Petrifilm™ Plates
- AOAC® Performance Tested Method™ Certificate #022104



3M™ Hydrated-Sponges and 3M™ Sponge-Sticks

- Ready-to-use, pre-hydrated sponges for easy and convenient sample collection
- Biocide-free cellulose sponge maintains viability for a wide range of organisms
- Durable sponge withstands scrubbing
- Offered in multiple buffers
- Allows sampling without direct handling of the sponge



3M™ Sample Bags

- Durable, quality construction
- Easy-to-open bag
- Wire ties for strong seals and reduced leakage

Sample Preparation and Processing Products




3M™ Flip-Top Dilution Bottles

- Leak-resistance, strong, lightweight clear PET plastic
- Pre-filled, eliminate making and sterilizing media
- Hinged lid stays open for easy sample addition
- Offered in multiple buffers



3M™ Mini Flip-Top Dilution Vials

- Hinged lid stays open for easy sample addition
- Requires less sample material than regular Flip-Top dilution bottles
- Offered in multiple buffers



3M™ Homogenizer Bags

- Quality construction to survive your paddle blender



3M™ Dilution Pouches

- Flat bottom stand-up pouch ensures stability
- Wire tie holds top open for easy sample addition
- Multi-laminated pouch is puncture resistant



3M™ Dehydrated Media

- Convenient
- Consistent
- Cost-effective

Create a seamless sample collection program with 3M's integrated product offerings.

Product	Description	Catalog #
3M™ Quick Swab	1mL Lethen Broth	6432
	1mL Lethen Broth	6433
3M™ Swab-Sampler	1mL Lethen Broth	RS9601LET
	4mL Lethen Broth	RS9604LET
	5mL Lethen Broth	RS9605LET
	10mL Lethen Broth	RS96010LET
	10mL Buffered Peptone Water	RS96010BPW
3M™ Environmental Scrub Sampler Stick	10mL 3M™ Wide Spectrum Neutralizer	ESS10WSN
3M™ Environmental Scrub Sampler and Gloves	10mL 3M™ Wide Spectrum Neutralizer	HES10WSN2G
3M™ Hydrated-Sponge	10mL Buffered Peptone Water	HS10BPW
	10mL Buffered Peptone Water, with 2 gloves	HS10BPW2G
	10mL Lethen Broth	HS10LET
3M™ Sponge-Stick	10mL Lethen Broth, with 2 gloves	HS10LET2G
	Without Media	SSL100
	10mL Buffered Peptone Water	SSL10BPW
	10mL Lethen Broth	SSL10LET
3M™ Dry-Sponge	10mL Lethen Broth, with 2 gloves	SSL10LETG
	Dry sponge, with 18 oz. sample bag	BP133ES
3M™ Letheen Broth	Dry sponge with 18 oz. sample bag, with 2 gloves	BP237SPG
	500g	BP0245500
3M™ Mini Flip-Top Vials with Lethen	9mL Lethen Broth	BPPFV9LB
	10mL Lethen Broth	BPPFV10LB
3M™ Sample Bags, Plain	3" x 7", 4 oz., 2.5 mil	BP111S
	4.5" x 7.5", 12 oz., 2.5 mil	BP134S
	4.5" x 9", 15 oz., 2.5 mil	BP112S
	5.5" x 9", 22 oz., 3 mil	BP113S
	4.5" x 12", 27 oz., 2.5 mil	BP118S
	4.5" x 15", 36 oz., 2.5 mil	BP119S
	5.5" x 15", 50 oz., 3 mil	BP121S
	12" x 12", 120 oz., 4 mil	BP41212
	10" x 15", 135 oz., 4 mil	BP1015LP
12" x 18", 228 oz., 4 mil	BP41218	
3M™ Sample Bags, Printable	3" x 5", 2 oz., 2.5 mil	BP123S
	3" x 7", 4 oz., 2.5 mil	BP114S
	4.5" x 9", 15 oz., 2.5 mil	BP115S
	5.5" x 9", 22 oz., 3 mil	BP127S
	7" x 12", 55 oz., 3 mil	BP701
	10" x 15", 135 oz., 4 mil	BP1015

Product	Description	Catalog #
3M™ Mini Flip-Top Vials	9mL Butterfields Buffer	BPPFV9BFD
	9mL Buffered Peptone Water	BPPFV9BPW
	10mL Buffered Peptone Water	BPPFV10BPW
	9mL MRD, Max Rec Diluent	BPPFV9MRD
	9mL MRS Broth	BPPFV9MRS
3M™ Flip-Top Dilution Bottles	90mL Buffered Peptone Water (ISO)	FTBPW9060
	90mL Buffered Peptone Water	FTBPW9060
	90mL Butterfields Buffer	FTBFD9060
	99mL Butterfields Buffer	FTBFD9960
	90mL Peptone Water	FTPW9060
	99mL Peptone Water	FTPW9960
	90mL Phosphate Buffer	FTPHB9060
3M™ Sample Bags, Blender	99mL Phosphate Buffer	FTPHB9960
	No wire tie, 7" x 12", 55 oz., 3 mil	BP912S
3M™ Sample Bags, Blender with Filter	No wire tie, 7" x 12", 55 oz., 4 mil	BP913S
	No wire tie, 7" x 12", 55 oz., 3 mil	BP914S
	No wire tie, 15" x 20", 405 oz., 3 mil	BP1520F
3M™ Quick Dilution Pouch	No wire tie, 6.75" x 12", 55 oz., 3 mil	6469
	99mL Butterfields Buffer	QDBFD99
	225mL Butterfields Buffer	QDBFD225
	225mL Butterfields Buffer, with filter	QDBFD225F
3M™ Bird Rinse	450mL Butterfields Buffer	QDBFD450
	400mL Butterfields Buffer	BRBFD40030
	400mL Buffered Peptone Water	BRBPW40030
3M™ Buffered Peptone Water Broth	500g	BP0115500
	2.5kg	BP0115005
3M™ Buffered Peptone Water, ISO	500g	BPW500
	2.5kg	BPW025
	10kg	BPW010
3M™ Salmonella Enrichment Base	500g	SEB500
	2.5kg	SEB025
3M™ Salmonella Enrichment Supplement	1g/vial	SESUP001



3M Food Safety
 3M Center, Building 275-5W-05
 St. Paul, MN 55144-1000 USA

1-800-328-6553
 3M.com/foodsafety

3M and Petrifilm are trademarks of 3M. All other trademarks are property of their respective owners. Please recycle. Printed in USA.
 © 3M 2021. All rights reserved. 70-2011-5188-6 (Rev-07-2021)