

Upgrade to Dual Performance

Introducing the New 3M[™] Attest[™] dual Auto-reader BI System with Any-Well Technology for Steam and VH2O2





The Power of 2-in-1

Improve the efficiency of how you process BIs with 3M's latest innovation, so you can do more with less. Introducing our 3M[™] Attest[™] Auto-reader technology that incubates and delivers BI results for both steam and vaporised hydrogen peroxide sterilisation in one reader.

Any-Well Technology

The new 24-minute Dual System works with Any-Well technology, which means Attest[™] 1492V (for steam) or 1295 (for vaporised hydrogen peroxide sterilisation — VH2O2) BIs can be incubated in any well of the reader at the same time.

An Upgrade that Gives You More

Get the peace of mind and confidence with a dual Auto-reader system that's:

Simple

Current customers can have their existing Attest™ Auto-reader model 490 (steam) or 490H (VH202) upgraded to start incubating both steam and VH2O2 biological indicators in one unit.



Consistent

3M offers 24-minute biological indicator results for both steam and vaporised hydrogen peroxide (VH2O2) providing the ease of consistency.

Familiar & Cost-Sensitive

Use the same trusted Attest[™] systems with no change to 3M BIs or test packs. The upgrade is free, with no new hardware to buy.

Connected

Use the same 3M[™] Web Application you use now for real-time BI status and results and for digital record keeping.

Your reader will still connect to leading Instrument Tracking Software Systems.

24 minute Bl results*



For more details on how to upgrade your current 3M[™] Attest[™] 490 or 490H Auto-Reader, contact your 3M Sales Representative.

Current Attest[™] Auto-reader 390 users or new 3M customers can contact a 3M sales rep to purchase a new model 490 Auto-reader pre-configured with dual functionality.

Steam Sterilisation Monitoring, Now Optimised

3M[™] Attest[™] Super Rapid Readout Biological Indicators use the same technology you've trusted for years, now optimised to deliver an even faster result.

Standard Compliance

Attest[™] Super Rapid Readout BIs meet the performance requirements of ISO 11138-1:2017 and ISO 11138-3:2017.

Trusted Steam Sterilisation Monitoring Portfolio

- Get the performance of the same trusted technology you are familiar with by using 3M[™] Attest[™] Super Rapid Readout process challenge devices (PCDs).
- 3M[™] Attest[™] system is the #1 utilised steam sterilisation test pack technology in Australia and New Zealand, and is now available with the speed of 24-minute results.
- 3M[™] Comply[™] Chemical indicators, Bowie-Dick test packs and indicator tape provide the core four steam sterilisation monitoring products.





*When used with 3M[™] Attest[™] Auto-reader 490 or the Attest[™] Auto-reader 490H with software version 4.2.7 or greater.

Know. Then Go.

Now Every Load Monitoring is not just possible, it's practical. Take control with 24-minute results and get the peace of mind that comes from knowing that every patient will receive a high standard of care.

\checkmark	

	`	
	` 🗸	
\bullet		

Results You Can Count on for Vaporised Hydrogen Peroxide Sterilisation

The 3M[™] Attest[™] system for VH2O2 provides confidence, in just 24 minutes. The 1295 Biological Indicator system is cleared for use in the following machines and cycles:

STERRAD[®] Vaporised Hydrogen Peroxide Systems

For use with STERRAD® 100S, STERRAD® 100NX (Standard, Flex, Express and Duo cycles), STERRAD® NX (Standard and Advanced cycles), STERRAD® 100NX with ALLClear™ Technology (Standard, Flex, Express and Duo cycles), STERRAD® NX with ALLClear™ Technology (Standard and Advanced cycles).

Amsco[®] V-PRO[®] Low Temperature Sterilisation System

For use with V-PRO® maX (Lumen, Non Lumen, and Flexible cycles), V-PRO® 1 plus (Lumen and Non Lumen cycles), V-PRO® 1 (Lumen cycles), V-PRO® 60 (Lumen, Non Lumen, and Flexible cycles).

Successfully validated and assessed for compatibility with the sterilisers outlined below: HMTS-80E, HMTS-142, HMTS-142D and HMTS -30E.

3M Provides a Complete VH2O2 System for Monitoring

- Use the key three monitoring products from 3M— BIs, CIs, and indicator tape—to leverage the full system of sterilisation assurance products.
- A consistent blue-towards-pink color change for all VH2O2 indicator products including the CI on the biological indicator, internal process CI and external indicator tape, makes it easy to read and interpret results for all of your VH2O2 monitoring needs.



3M Attest™ EEE 1295 Date Ster.


The 3M[™] Attest[™] Brand Advantage



The 3M[™] Attest[™] brand advantage gives you a standardised, simplified and supported system for sterilisation across steam, VH2O2 and EO modalities.

Fight SSIs from Every Angle

Our 3M portfolio of sterilisation assurance products and sterilisation monitoring expertise is your first line of defense in your fight to help protect patients from surgical site infections.





Pack Control Exposure Control

Product	Description	Packaging	Shipping
490/490H	3M [™] Attest [™] Auto-reader 490/490H	1 unit/box	1 box/case
1295	3M [™] Attest [™] Rapid Readout Biological Indicator	30 each/bag	4 bags/case
1248	3M [™] Comply [™] Chemical Indicator Strip	250 each/bag	4 bags/case
1228	3M [™] Comply [™] Indicator Tape	1,9cm roll	24 rolls/case

For more details on how to upgrade your current 3M[™] Attest[™] 490 or 490H Auto-Reader, contact your 3M Sales Representative.

For more information or support visit <u>go.3M.com/3MAdvantage</u> or call the 3M Health Care Help Line at Australia: 1300 363 878 or New Zealand: 0800 80 81 82.



3M Medical Solutions Division 3M Australia Pty Limited Bldg A, 1 Rivett Road North Ryde NSW 2113 1300 363 878 www.3M.com.au/healthcare

3M New Zealand Limited 94 Apollo Drive Rosedale Auckland 0632 0800 80 81 82 www.3M.com/healthcare

3M, Attest, Comply and SteriGage are trademarks of 3M. All other trademarks are property of their respective owners. Please recycle. Printed in Australia. © 3M 2018. All rights reserved. PB8794.06.18