





# The Power of 2-in-1

Improve the efficiency of how you process BIs with 3M's latest innovation, so you can do more with less. Introducing our 3M™ Attest™ Auto-reader technology that incubates and delivers BI results for both steam and vaporised hydrogen peroxide sterilisation in one reader.

## Any-Well Technology

The new 24-minute Dual System works with Any-Well technology, which means Attest™ 1492V (for steam) or 1295 (for vaporised hydrogen peroxide sterilisation — VH<sub>2</sub>O<sub>2</sub>) BIs can be incubated in any well of the reader at the same time.



# Steam Sterilisation Monitoring, Now Optimised

3M™ Attest™ Super Rapid Readout Biological Indicators use the same technology you've trusted for years, now optimised to deliver an even faster result.

## Standard Compliance

Attest™ Super Rapid Readout BIs meet the performance requirements of ISO 11138-1:2017 and ISO 11138-3:2017.

## Trusted Steam Sterilisation Monitoring Portfolio


- Get the performance of the same trusted technology you are familiar with by using 3M™ Attest™ Super Rapid Readout process challenge devices (PCDs).
- 3M™ Attest™ system is the #1 utilised steam sterilisation test pack technology in the Australia and New Zealand, and is now available with the speed of 24-minute results.
- 3M™ Comply™ Chemical indicators, Bowie-Dick test packs and indicator tape provide the core four steam sterilisation monitoring products.



\*When used with 3M™ Attest™ Auto-reader 490 or the Attest™ Auto-reader 490H with software version 4.2.7 or greater.

## Know. Then Go.

Now Every Load Monitoring is not just possible, it's practical. Take control with 24-minute results and get the peace of mind that comes from knowing that every patient will receive a high standard of care.



# Results You Can Count on for Vaporised Hydrogen Peroxide Sterilisation

The 3M™ Attest™ system for VH2O2 provides confidence, in just 24 minutes. The 1295 Biological Indicator system is cleared for use in the following machines and cycles:

## STERRAD® Vaporised Hydrogen Peroxide Systems

For use with STERRAD® 100S, STERRAD® 100NX (Standard, Flex, Express and Duo cycles), STERRAD® NX (Standard and Advanced cycles), STERRAD® 100NX with ALLClear™ Technology (Standard, Flex, Express and Duo cycles), STERRAD® NX with ALLClear™ Technology (Standard and Advanced cycles).

## Amsco® V-PRO® Low Temperature Sterilisation System

For use with V-PRO® maX (Lumen, Non Lumen, and Flexible cycles), V-PRO® 1 plus (Lumen and Non Lumen cycles), V-PRO® 1 (Lumen cycles), V-PRO® 60 (Lumen, Non Lumen, and Flexible cycles).

## 3M Provides a Complete VH2O2 System for Monitoring

- Use the key three monitoring products from 3M— BIs, CIs, and indicator tape—to leverage the full system of sterilisation assurance products.
- A consistent blue-towards-pink color change for all VH2O2 indicator products including the CI on the biological indicator, internal process CI and external indicator tape, makes it easy to read and interpret results for all of your VH2O2 monitoring needs.




## The 3M™ Attest™ Brand Advantage



The 3M™ Attest™ brand advantage gives you a standardised, simplified and supported system for sterilisation across steam, VH2O2 and EO modalities.

# Fight SSIs from Every Angle


Our 3M portfolio of sterilisation assurance products and sterilisation monitoring expertise is your first line of defense in your fight to help protect patients from surgical site infections.



**STEAM**

Load Control  
Pack Control  
Exposure Control  
Equipment Control

Product	Description	Packaging	Shipping
490/490H	<b>3M™ Attest™ Auto-reader 490/490H</b>	1 unit/box	1 box/case
1492V	<b>3M™ Attest™ Super Rapid Readout Biological Indicator</b>	50 each/box	4 box/case
41482V	<b>3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack</b>	24 each/box	1 box/case
1243	<b>3M™ Comply™ (SteriGage™) Chemical Integrators</b>	500 each/bag	2 bags/case
		100 each/bag	10 bags/case
1355	<b>3M™ Comply™ Indicator Tapes</b>	18mm roll	28 rolls/case
		24mm roll	20 rolls/case
4208/4109	<b>3M™ Electronic Test System</b>	1 each/case	1 each/case
1300	<b>3M™ Comply™ Bowie-Dick Test Pack 1300</b>	20 packs/box	1 boxes/case
1301	<b>3M™ Comply™ Bowie-Dick Test Pack 1301</b>	20 packs/box	1 boxes/case



**VH2O2**

Load Control  
Pack Control  
Exposure Control

Product	Description	Packaging	Shipping
490/490H	<b>3M™ Attest™ Auto-reader 490/490H</b>	1 unit/box	1 box/case
1295	<b>3M™ Attest™ Rapid Readout Biological Indicator</b>	30 each/bag	4 bags/case
1248	<b>3M™ Comply™ Chemical Indicator Strip</b>	250 each/bag	4 bags/case
1228	<b>3M™ Comply™ Indicator Tape</b>	1,9cm roll	24 rolls/case

**For more details on how to upgrade your current 3M™ Attest™ 490 or 490H Auto-Reader, contact your 3M Sales Representative.**

For more information or support visit [go.3M.com/3MAdvantage](http://go.3M.com/3MAdvantage) or call the 3M Health Care Help Line at Australia: 1300 363 878 or New Zealand: 0800 80 81 82.



**3M Infection Prevention Division**  
3M Australia Pty Limited  
Bldg A, 1 Rivett Road  
North Ryde NSW 2113  
1300 363 878  
[www.3M.com.au/healthcare](http://www.3M.com.au/healthcare)

3M New Zealand Limited  
94 Apollo Drive  
Rosedale Auckland 0632  
0800 80 81 82  
[www.3M.com/healthcare](http://www.3M.com/healthcare)

3M, Attest, Comply and SteriGage are trademarks of 3M. All other trademarks are property of their respective owners.

Please recycle. Printed in the Australia.  
© 3M 2018. All rights reserved.  
PB8794.06.18