



**Choose the right stethoscope to help you perform at your best.**

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the reliable acoustic performance, high quality and comfort you need with a Littmann stethoscope.



	Lightweight II S.E. Stethoscope	Classic II Pediatric and Infant Stethoscopes	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	CORE Digital Stethoscope
	Ideal for taking blood pressure readings and making limited patient assessments.	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.		Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more with amplification, noise cancellation and mobile connectivity.
<b>Approximate Length</b>	28"/71cm	28"/71cm	27"/69cm	22"/56cm, 27"/69cm	27"/68cm	27"/69cm
<b>Approximate Weight</b>	118g	105g, 95g	150g	167g, 177g	185g	232g
<b>Warranty</b>	2 years	3 years	5 years	7 years	7 years	2 years
<b>Diaphragm</b>	Tunable	Non-Tunable (Floating)	Tunable (2)	Tunable (2)	Tunable	Tunable (2)
<b>Approximate Chestpiece Diameter</b>	Double-sided 2.0 x 1.8"/5.2 x 4.5cm 1.4 x 1.3"/3.5 x 3.2cm	Double-sided Pediatric 1.4"/3.4cm 1.1"/2.7cm Infant 1.1"/2.8cm 0.7"/1.7cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Single-sided 1.9"/4.8cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm
<b>Features</b>	Teardrop-shaped chestpiece Fits under blood pressure cuff	Specially designed for the smallest patients	Convertible open bell	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Up to 40x amplification* Switch from analog listening to digital  Connects to Eko software on a smart device** to capture, store and share data



**The power to learn, in the palm of your hand.**

3M™ Littmann® Learning App – Download free version in the App Store® or Google Play™ Store

**Need more help? Visit [Littmann.com.au/FindMyStethoscope](http://Littmann.com.au/FindMyStethoscope)**  
**Need more help? Visit [Littmann.co.nz/FindMyStethoscope](http://Littmann.co.nz/FindMyStethoscope)**

1300 363 878 | [www.littmann.com.au](http://www.littmann.com.au)  
0 800 80 8182 | [www.littmann.co.nz](http://www.littmann.co.nz)

\*At peak frequency, vs. analog mode  
\*\*Not included



# 3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide

**Components**  
(Chestpiece/Eartube/Stem)

Lightweight II SE Stethoscope	Classic II Infant Stethoscope	Classic II Pediatric Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	CORE Digital Stethoscope
-------------------------------	-------------------------------	----------------------------------	-------------------------------------	---------------------------------------	--------------------------------	--------------------------

**Tube Colour**

All Standard Finish	2450 Black	2114 Black	2113 Black	5620 Black	5627 Burgundy	6156 Plum	2160 Black	
	2451 Burgundy	2114R Red	2113R Red	5621 Gray	5633 Pearl Pink	6152 Black	2163 Burgundy	
	2452 Caribbean	2124 Caribbean	2119 Caribbean	5622 Navy	5831 Plum	6154 Navy	2164 Navy	
	2454 Ceil Blue		2122 Raspberry	5623 Caribbean	5832 Lavender	6155 Hunter Green	2167 Plum	
	2456 Pearl Pink			5630 Ceil Blue	5835 Turquoise			
				5839 Lemon-Lime				
All Black			5803 Black	5868 Burgundy	6163 Black	6168 Navy	2161 Black	8480 Black
All Brass Finish								
All Copper Finish			5809 Chocolate					
All Smoke			5811 Black					
Rainbow/Brass/Brass		2157 Raspberry		5806 Raspberry	5807 Caribbean			
Brass/Smoke/Smoke								
Rainbow/Black/Black			5870 Black					
Mirror/Black/Black								8890 Black
Mirror/Smoke/Smoke				5863 Navy				
			5959 Ceil Blue					
Mirror/Standard/Mirror					6177 Black	6170 Burgundy		
Champagne/Smoke/Smoke			5861 Black	5864 Burgundy	6179 Black			
Champagne/Black/Champagne					6176 Burgundy	6190 Caribbean		
High-Polish Copper/Black/Black								8870 Black
High-Polish Rainbow/Black/Black						6242 Navy		8570 Black
High-Polish Rainbow/Smoke/Smoke						6240 Black	6241 Raspberry	
<b>COLOURED STEM SERIES</b>	<b>Components:</b> Chestpiece/Eartube/Stem			<b>COLOURED STEM SERIES: (square inset represents stem colour)</b>				
Black/Black/Colour	Standard Finish		Rainbow			6200 Black	6203 Black	
	Black		Mirror			6201 Black		
High-Polish Smoke/Black/Colour	Brass Finish		Champagne			6202 Navy	6204 Black	
	Copper Finish		High-Polish Smoke			6232 Black	6234 Caribbean	
High-Polish Smoke/Smoke/Colour	Smoke		High-Polish Rainbow			6238 Gray		
	High-Polish Rainbow		High-Polish Copper					
Rainbow/Black/Colour						6205 Plum		
Smoke/Smoke/Colour				5872 Turquoise	5875 Lime			
				5873 Gray				
Mirror/Smoke/Colour				5960 Plum	5962 Pearl Pink			

Colours shown are approximate.

1300 363 878 | [www.littmann.com.au](http://www.littmann.com.au)  
0800 80 81 82 | [www.littmann.co.nz](http://www.littmann.co.nz)

© 3M 2021. All rights reserved. 3M, Littmann, the L Littmann logo, Classic III, Cardiology IV, and Master Cardiology are trademarks of 3M. Eko trademarks, including CORE, are marks and/or registered marks of Eko Devices Inc. All other trademarks are property of their respective owners. 70-2011-8157-8, Rev. 28/01/22