

Soft Putty

Vinyl Polysiloxane Impression Material



Soft but Strong

Features and Benefits:

- Easy Mix with soft handling
- **Multiple Pours** - Allows for multiple pours compared to alginate material
- **Patient Friendly** -Blueberry mint flavor with Xylitol is more pleasing to patients
- Accurate & Dimensionally stable
- Cost effective

A greater level of confidence with a new generation of VPS impression material !

After more than 40 years of developing innovative impressing solutions, 3M ESPE has a product portfolio that meets virtually all dentists' needs and preferences.

A world leader in Impressing materials, 3M ESPE now brings you Soft Putty, a value based VPS impressing material . Its good balance of tear strength, dimensional stability and elasticity makes it a great choice for your daily impressing needs.

Indications :

- Conventional Crown and Bridge impressions
- Preliminary impressions for complete denture cases
- Provisional crown and bridge impressions
- Study model impressions
- Opposing dentition impressions
- Impressions for orthodontic models
- Impressions for bleaching tray and mouth guards



Important clinical advantages over c-silicones

	C-Silicone (e.g., Speedex)	VPS Wash Material
Dimensional stability	—	✓
Very high tear strength	—	✓
Multiple pouring	✓	✓
Long-term storage	—	✓
Hydrophilicity	✓	✓
Automatic mixing	—	✓

Condensation silicone (cc-silicone) materials display inferior dimensional stability and continue to shrink after setting.

Great Compatibility with light body impression materials

Soft Putty has great compatibility with existing Light body impression materials. Separation doesn't occur and the color is clearly distinguishable. That enables you to make precise impressions.



Fresh Flavour

Friendly flavour from blueberry mint and Xylitol removes chemical odor of rubber impressions and make patients cooperative to dentists.






Excellent stability

Soft Putty passed a stability test whose standard has been currently enhanced by Korean FDA. This guarantees excellent stability within shelf life. (Safety test completed in Jun. 2012 at Kyung-Hee University in South Korea.)

Soft, but Strong

Soft Putty is soft during mixing, but strong after mixing. High strength and elasticity prevent impressions from damage and distortion.



Graph 1. Prevents impressions from being torn out when removing from teeth

Graph 2. Minimizes distortion when removing impressions from teeth

Graph 3. Minimizes distortion after setting.

*Each test followed ISO 4823:2000 (Dentistry-Elastomeric impression materials)
*Source: 3M Korea Laboratory test

Ordering Information

Item Number - 70326

Description - 3M ESPE Soft Putty™

Components - 1 Base Paste (300ml), 1 Catalyst Paste (300ml), 2 Scoops (One Orange Base Scoop and One Purple Catalyst Scoop) and 1 Instructions For Use

3M

3M India Limited

Concorde Block, UB City, 24,

Vittal Mallya Road,

Bangalore - 560001

Tel : +91-80- 22231414/66595759

Fax : +91-80-22231450

e-mail: 3mespe.in@mmm.com

www.3mespeindia.com

Toll Free: 1-800-425-3030



3M ESPE India

Please recycle. Printed in India
© 3M 2013. All rights reserved.

