

3M Science.
Applied to Life.™

Flexible



3M™ Medipore™ H Soft Cloth Surgical Tape

Made to move with your patients.



The choice for stretchable, conformable securement.

Medipore H Tape features a soft, conformable backing with bi-directional stretch to accommodate swelling and movement, while minimizing the risk for skin injury. It's now available in individually packaged single-use rolls to help reduce the risk of cross-contamination.



Conformability

Bi-directional stretch accommodates swelling, distention and movement.



Breathability

Highly breathable, soft cloth tape.



Easy to use

Easy tear at perforations with no scissors necessary.



Helps reduce risk

Individually packaged single-use rolls can help protect patients by reducing the risk of cross-contamination.*

Helping you focus on
patient care and safety



74%
of tape rolls

In one study, 74% of partially used tape rolls from various sites within the hospital had some bacterial growth.¹

*Individually packaged single-use rolls help prevent tape from being exposed to environmental contaminants, minimize contact with surfaces and equipment, and minimize exposure to healthcare worker hands.

3M™ Medipore™ H Soft Cloth Surgical Tape

Clinical Application: **Flexible Securement**

- Dressings and added pressure
- Secondary securement for:
 - ▶ IV lines and tubing
 - ▶ Central venous catheters
 - ▶ Chest tubes
 - ▶ Surgical drain tubes
- When swelling or movement is anticipated
- When repeated taping is needed

Ideal for patients who need to maintain movement or are experiencing swelling or distension.

Tape material: Perforated nonwoven soft cloth



Ordering information

Catalog Number	Size	Individually packaged single-use rolls	Rolls	Case	HCPCS Code
2860S-1	1 in. x 2 yd. (2,5 cm x 1,8 m)	✓	72 per bag	1 bag	A4452
2860S-2	2 in. x 2 yd. (5 cm x 1,8 m)	✓	48 per bag	1 bag	A4452
2860S-4	4 in. x 2 yd. (10,1 cm x 1,8 m)	✓	24 per bag	1 bag	A4452
2860S-6	6 in. x 2 yd. (15,2 cm x 1,8 m)	✓	16 per bag	1 bag	A4452
2861	1 in. x 10 yd. (2,5 cm x 9,1 m)	—	24 per box	1 box	A4452
2862	2 in. x 10 yd. (5 cm x 9,1 m)	—	12 per box	1 box	A4452
2863	3 in. x 10 yd. (7,6 cm x 9,1 m)	—	12 per box	1 box	A4452
2864	4 in. x 10 yd. (10,1 cm x 9,1 m)	—	12 per box	1 box	A4452
2866	6 in. x 10 yd. (15,2 cm x 9,1 m)	—	12 per box	1 box	A4452
2868	8 in. x 10 yd. (20,3 cm x 9,1 m)	—	6 per box	1 box	A4452



For more information about 3M™ Medipore™ H Soft Cloth Surgical Tape and other 3M Medical Securement Solutions, visit [3M.com/MedicalSecurement](https://www.3m.com/MedicalSecurement).

Medical Solutions Division
3M Health Care
2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-228-3957
Web 3M.com/medical

¹Lavelle BE. Reducing the risk of skin trauma related to medical adhesives, *Managing Infection Control*, June 2004.

3M and Medipore are registered trademarks of 3M. Please recycle.
Printed in USA. © 3M 2020. All rights reserved. 70-2011-7899-6