

3M

Tegaderm™

Brand Wound Care Solutions



When monitoring an incision, you need the whole picture.

3M™ Tegaderm™ Absorbent Clear Acrylic Dressings

Maintaining an ideal healing environment and monitoring surgical incisions are critical steps in reducing post-operative surgical site infections. For small to large incisions, Tegaderm™ Absorbent Clear Acrylic Dressings offer a waterproof¹, transparent² barrier – allowing active monitoring and early intervention without disturbing the incision³. Plus, extended wear time⁴ allows patients to shower and return more quickly to daily activities¹.



Available in five sizes for complete coverage of incisions up to 28 cm (11 in).

3M™ Tegaderm™ Absorbent Clear Acrylic Dressings

At 3M, we partner with providers to develop surgical solutions that have the power to transform patient outcomes through science, including the post-operative protection and care of surgical incisions. Choose from our full line of transparent Tegaderm™ Absorbent Clear Acrylic Dressings for a variety of applications or procedures.



Enhanced monitoring

Designed with a clear film and acrylic absorbent pad composed of a hydrophilic polymer⁵, Tegaderm absorbent clear acrylic dressings give providers and patients the ability to monitor incisions for early signs of infection.²



Extended protection




Wear time of up to 28 days⁴ gives patients extended protection beyond their initial hospital stay, and supports early showering¹ without removing the dressing.



Optimal healing environment

Helps protect from external contaminants, bodily fluids, bacteria, viruses and water ingress¹, plus is designed to absorb moderate exudate levels² to promote an optimal healing environment.⁶

Dressing size selection: Choose the right dressing size for your application.

	Catalogue Number	Dressing Pad Size	Overall Dressing Size	Incision Size	Dressings/Box	Boxes/Case
	90800 Oval	3.8 x 5.7 cm 1.5 in x 2.25 in	7.6 x 9.5 cm 3 in x 3.75 in	≤ 2.5 cm ≤ 1 in	5	6
	90801 Oval	6.0 x 7.6 cm 2.375 in x 3 in	11.1 x 12.7 cm 4.375 in x 5 in	≤ 4.0 cm ≤ 1.5 in	5	6
	90803 Oval	8.5 x 10.7 cm 3.375 in x 4.25 in	14.2 x 15.8 cm 5.625 in x 6.25 in	≤ 10.0 cm ≤ 4 in	5	6
	90815 Rectangle	4.5 x 20.0 cm 1.75 in x 8 in	9.0 x 25.0 cm 3.5 in x 10 in	≤ 18.0 cm ≤ 7 in	30	4
	90817 Large Rectangle	4.5 x 30.0 cm 1.75 in x 11.75 in	9.0 x 35.0 cm 3.5 in x 13.75 in	≤ 28.0 cm ≤ 11 in	30	4

Application tip: Dressing pad should extend past ends of incision by 1 cm (.375 in).

For more information, contact your 3M Health Care sales representative, call the 3M Health Care Customer Helpline at 1-800-364-3577 or visit 3M.ca/Medical.



Available in Canada from

3M Medical Solutions Division
3M Canada
 300 Tartan Drive
 London, ON N5V 4M9
 Canada
 1-800-364-3577
www.3M.ca/Medical

3M Health Care
 2510 Conway Avenue
 St. Paul, MN 55144
 USA
 1-800-228-3957

¹ 3M data on file. EM-05-014684.

² 3M data on file. EM-05-000138.

³ 3M data on file. EM-05-000003.

⁴ 3M data on file. EM-05-014692, EM-05-014725.

⁵ 3M data on file. EM-05-137984.

⁶ 3M data on file. EM-05-001314.