

3M Science.
Applied to Life.™

Non-adhesive



3M™ Coban™ NL Non-Latex Self-Adherent Wrap

**Designed to adhere
to itself, not skin.**



Comfortable and versatile hold.

Give patients the protection they deserve with Coban NL Wraps. Unlike tape and some other wraps, Coban NL Wraps adhere to themselves without adhesive or clips. They're also available in individually packaged single-patient-use rolls to help reduce the risk of cross-contamination.*



Reliable hold

Maintains compression and is an excellent choice when you don't want to stick to skin.



Skin performance

Delivers the securement power you need while minimizing damage to skin.



Easy to use

Can be torn by hand and is easily adjusted during application.



Helps reduce risk


Individually packaged single-patient-use rolls can help protect patients by reducing the risk of cross-contamination.*

**Gentle hold.
Minimal skin irritation.**

**27% less
irritating**

Coban NL Wrap is 27% less irritating to skin than CoFlex NL Wrap.¹

24-hour self-assessment of irritation



Scale:

0 = less irritation to 100 = more irritation

3M™ Coban™ NL Non-Latex Self-Adherent Wrap

Clinical Application: Non-adhesive Securement

- Blood draws
- Dressings
- Immobilization
- Secondary securement for difficult to dress areas (head, fingers, toes)
- Support and mild compression for soft tissue injuries (i.e. strains, sprains)**

Suitable for all patients, including those with fragile and at-risk skin.

Wrap material: Nonwoven, elastic

Ordering information

Catalog Number	Size		Individually packaged single-patient-use rolls	Sterile	Rolls	HCPCS Code
2082-1X	2 in. x .75 yd. (50 mm x 0,6 m) 2 in. x 2 yd. (50 mm x 1,5 m)	↔ ↔	✓	—	36 per bag	A6453
2083-1X	3 in. x .75 yd. (75 mm x 0,6 m) 3 in. x 2 yd. (75 mm x 1,5 m)	↔ ↔	✓	—	24 per bag	A6454
2084-1X	4 in. x .75 yd. (100 mm x 0,6 m) 4 in. x 2 yd. (100 mm x 1,5 m)	↔ ↔	✓	—	18 per bag	A6454
2086-1X	6 in. x .75 yd. (150 mm x 0,6 m) 6 in. x 2 yd. (150 mm x 1,5 m)	↔ ↔	✓	—	12 per bag	A6455
2082S	2 in. x 5 yd. (50 mm x 4,5 m)		—	✓	36 per bag	A6453
2083S	3 in. x 5 yd. (75 mm x 4,5 m)		—	✓	24 per bag	A6454
2084S	4 in. x 5 yd. (100 mm x 4,5 m)		—	✓	18 per bag	A6454
2086S	6 in. x 5 yd. (150 mm x 4,5 m)		—	✓	12 per bag	A6455
2081	1 in. x 5 yd. (25 mm x 4,5 m)		—	—	30 per box	A6453
20815	1.5 in. x 5 yd. (37,5 mm x 4,5 m)		—	—	48 per box	A6453
2082	2 in. x 5 yd. (50 mm x 4,5 m)		—	—	36 per box	A6453
2083	3 in. x 5 yd. (75 mm x 4,5 m)		—	—	24 per box	A6454
2084	4 in. x 5 yd. (100 mm x 4,5 m)		—	—	18 per box	A6454
2086	6 in. x 5 yd. (150 mm x 4,5 m)		—	—	24 per box	A6455
2081C (5 rolls of each)	1 in. x 5 yd. (25 mm x 4,5 m) ●●●●●		—	—	30 per box	A6453
20815C (8 rolls of each)	1.5 in. x 5 yd. (37,5 mm x 4,5 m) ●●●●●●●●		—	—	48 per box	A6453
2082C (6 rolls of each)	2 in. x 5 yd. (50 mm x 4,5 m) ●●●●●●		—	—	36 per box	A6453
2083C (4 rolls of each)	3 in. x 5 yd. (75 mm x 4,5 m) ●●●●●●		—	—	24 per box	A6454
2084C (3 rolls of each)	4 in. x 5 yd. (100 mm x 4,5 m) ●●●●●●		—	—	18 per box	A6454
2084L	4 in. x 6.5 yd. (100 mm x 6 m)		—	—	12 per box	A6454



Skin performance:
The feather designates products that deliver the securement power you need while minimizing damage to skin.

For more information about 3M™ Coban™ NL Non-Latex Self-Adherent Wrap and other 3M Medical Securement Solutions, visit [3M.com/MedicalSecurement](https://www.3m.com/MedicalSecurement).



Medical Solutions Division
3M Health Care
2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-228-3957
Web 3M.com/medical

3M medical securement products are intended for use on intact skin.

¹Data on file with 3M. EM-05-012882.

^{*}Individually packaged, single-use rolls help prevent tape from being exposed to environmental contaminants, minimizes contact with hospital surfaces and equipment, and exposure to healthcare worker hands.

^{**}May use as a component in a compression wrap system support only under supervision of a wound care specialist.

3M and Coban are registered trademarks of 3M. Please recycle. Printed in USA. © 3M 2020. All rights reserved. 70-2011-7898-8