




A full line of restorative solutions to fit your unique needs.

With a rich legacy of oral care products and services, 3M is one of the world's most trusted companies. Filtek™ restoratives have been used in more than 1 billion restorations worldwide. Our family of Filtek™ dental restoratives lets you choose just the right option to fit any patient situation—from highly esthetic, to highly simple for an easier day.

Versatility & Esthetics		Esthetics & Simplicity		Efficiency & Simplicity	
<p>3M™ Filtek™ Supreme Ultra Universal Restorative</p> <ul style="list-style-type: none"> • Exceptional esthetics • Natural-looking fluorescence • Excellent polish retention • Exceptional handling • Outstanding strength for anterior and posterior • Wear resistance similar to enamel¹ • Can can be safely warmed² and injection molded • TRUE nanotechnology • 36 shades, 4 opacities 	<p>3M™ Filtek™ Supreme Flowable Restorative</p> <ul style="list-style-type: none"> • Virtually no bubbles or material “run-on” • Easy-to-bend cannula for improved access • Ergonomic design is easy to hold and inject • Excellent adaptation, polish retention and wear resistance • Can can be safely warmed³ and injection molded • 12 shades 	<p>3M™ Filtek™ Universal Restorative</p> <ul style="list-style-type: none"> • Predictable, repeatable esthetics • Long-lasting polish • Excellent, creamy handling • Low shrinkage stress/stress relief • Excellent wear resistance • Can can be safely warmed² and injection molded • TRUE nanotechnology • 8 designer shades plus an Extra White (XW) cover the 19 VITA classical and bleach shades⁴ • Pink Opaquer with excellent masking ability 	<p>3M™ Filtek™ One Bulk Fill Restorative</p> <ul style="list-style-type: none"> • Increased opacity for improved esthetics • One-step placement up to 5 mm • Excellent handling • Excellent adaptation • Stress relief • Excellent wear resistance • Can can be safely warmed² and injection molded • TRUE nanotechnology • 5 shades 	<p>3M™ Filtek™ Bulk Fill Flowable Restorative</p> <ul style="list-style-type: none"> • Virtually no bubbles or material “run-on” • Easy-to-bend cannula for improved access • Ergonomic design is easy to hold and inject • Injectable, one-step placement up to 4mm • Excellent adaptation, polish retention and low shrinkage • Can can be safely warmed³ and injection molded • 4 shades 	
					






¹ In vitro data on file.

² Warm capsules to 70°C (158°F) for up to 1 hour.

³ Up to 70°C (158°F) for up to 1 hour, up to 25 times.

⁴ 3M internal data.

We have you covered for every class of direct restoration.

	Versatility & Esthetics		Esthetics & Simplicity	Efficiency & Simplicity	
	3M™ Filtek™ Supreme Ultra Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Supreme Flowable Restorative (Incremental up to 2 mm)	3M™ Filtek™ Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ One Bulk Fill Restorative (Bulk up to 5 mm)	3M™ Filtek™ Bulk Fill Flowable Restorative (Incremental up to 4 mm)
Class I 	✓	✓ (Liner/Base)	✓	✓ (up to 4 mm)	✓ (Liner/Base)
Class II 	✓	✓ (Liner/Base)	✓	✓ (up to 5 mm)	✓ (Liner/Base)
Class III 	✓	✓	✓		✓
Class IV 	✓		✓		
Class V 	✓	✓	✓		✓
Minimally Invasive	✓	✓	✓	✓	✓
Stress Bearing	✓		✓	✓	

✓ = recommended use

Refer to Instructions for Use for full list of indications.



www.3M.com/Filtek

3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
Phone 1-888-363-3685

3M and Filtek are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2019. All rights reserved. (Rev. 10/19) 70-2013-0807-2