



Filtek™

Bulk Fill
Flowable Restorative

The bulk fill flowable that's easier to inject.

3M™ Filtek™ Bulk Fill Flowable Restorative is available in a new ergonomic syringe designed for improved comfort and easier injecting. This gives you more control and visibility for precise placement in deep restorations.

Reduce bubbles and waste

Virtually no bubbles or material “run-on” during dispensing gives you better control.

No more guessing

The green syringe barrel indicates the remaining material volume. When your syringe is empty, only the white plunger shows.



Bendable cannula

Easy-to-bend without kinking for better access.

Easy to hold and inject

Triangle finger plate and plunger make it quick and easy to personalise your grip.



Key benefits

- Virtually no bubbles or material “run-on”
- Easy-to-bend cannula for improved access
- Ergonomic design is easy to hold and inject
- Injectable, one-step placement up to 4 mm
- Excellent adaptation and low shrinkage stress

Can be warmed for injection molding[^]

Warm up to 70°C (158°F) for up to one hour, up to 25 times.



Posterior restoration made simple



Prepare the tooth and apply an adhesive e.g. 3M™ Scotchbond™ Universal Adhesive.



Dispense up to 4 mm of 3M™ Filtek™ Bulk Fill Flowable Restorative into the prepared cavity and cure, e.g. with 3M™ Elipar™ DeepCure-S LED Curing Light.



Place a 2 mm (minimum) capping layer of an aesthetic composite, such as a 3M™ Filtek™ Universal Restorative, and cure.

Clinical photos courtesy of Dr. Reich, Seefeld, Germany.

Visit 3M.com.au/FiltekFlowables



3M Oral Care

3M Australia Pty Limited
Building A, 1 Rivett Road
North Ryde, NSW 2113
1300 363 454
3M.com.au

3M New Zealand Limited
94 Apollo Drive
Rosedale, Auckland 0632
0800 808 182
3M.co.nz

Indications

(Incremental up to 4 mm)



Base
(Class I – II)



Liner
(Class I – II)



Minimally Invasive
Dentistry (MID)¹



Class III



Class V



Repair of Indirect
(small defects)



Repair²



Undercut
Block-out



Pit & Fissure
Sealant



Core
Build-up³

Refer to Instructions
for Use for full list of
indications.

Ordering information



3M™ Filtek™ Bulk Fill Flowable Restorative Refills

Syringe refills include: 2 – 2 g syringes;
20 dispensing tips.

SHADE	ITEM NO.
A1	4862A1
A2	4862A2
A3	4862A3
Universal	4862U

ITEM NO.	PRODUCT INFO.
3740T	3M™ Filtek™ Bulk Fill Flowable Restorative Dispensing Tips, 19-gauge, Green (20)



[^] Data on file.

1. Including small, non stress-bearing occlusal restorations.
2. Repair of acrylic temporary and composite resin.
3. Where at least half the coronal tooth structure is remaining to provide structural support for the crown.

Please always refer to the 3M Product Instructions For Use. 3M, Elipar, Filtek and Scotchbond are trademarks of 3M. Please recycle. Printed in Australia © 3M 2020. All rights reserved. AV011477393.