

Solutions scaled to fit your needs.



Not shown true to size and scale.



Now available! The 3M™ Attest™ Mini Dual Auto-reader.



3M™ Attest™ Dual Auto-reader System 490 and Mini Dual Auto-reader 490M

Get results for both steam and vaporized hydrogen peroxide (VH2O2) in one reader in just 24 minutes.* 3M now offers its dual auto-readers in two sizes to meet your facility's needs and budget.



Increase ease of use

Place any Attest 1491 or 1492V (for steam) or 1295 (for VH2O2) BI in any well of the reader at the same time.



Raise the standard of care

Help protect patients from the risk of infection with every load monitoring. With fast 24-minute BI results — it's possible.



Be connected

Connect with 3M product apps to get real time BI results and quickly integrate data into your facility's tracking system.

Innovation you can trust



3M is a leader in rapid self-contained BI technology for steam, vaporized hydrogen peroxide, and ethylene oxide.

850+
million
loads monitored



*When used with 3M™ Attest™ Auto-reader 490 or 490H with software version 4.2.7 or greater.

Maximizing space and affordability

With 3M™ Attest™ Dual Auto-reader size options you can standardize practice across multiple facilities for efficiency and training ease.

Our dual auto readers are compatible with the following 3M™ Attest™ biological indicators:

- 1491 Steam*
- 1492V Steam*
- 1295 VH2O2*



| | Catalog No. | Product | Description | Packaging | Shipping |
|-------|-------------|------------------------------------------------------|-------------------------------------------------------------------------------------------------------|-------------|--------------|
| STEAM | 490 | 3M™ Attest™ Dual Auto-reader | 22.2(W) x 8.25(H) x 22.2(D) cm | 1 unit/box | 1 box/case |
| | 490M | 3M™ Attest™ Mini Dual Auto-reader | 15.5(W) x 8.25(H) x 5.25(D) cm | 1 unit/box | 1 box/case |
| STEAM | 41482V | 3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack | Dynamic-air-removal (Pre-vacuum and SFPP): 132°C (270°F) 4 min. 135°C (275°F) 3 min. | 24 each/box | 1 box/case |
| | 1491 | 3M™ Attest™ Super Rapid Readout Biological Indicator | Gravity: 132°C (270°F) 3 min. and 10 min. 135°C (275°F) 3 min. and 10 min. | 50 each/box | 4 boxes/case |
| | 1492V | 3M™ Attest™ Super Rapid Readout Biological Indicator | Dynamic-air-removal (Pre-vacuum and SFPP): 132°C (270°F) 3 min. and 4 min. 135°C (275°F) 3 min. | 50 each/box | 4 boxes/case |
| VH2O2 | 1295 | 3M™ Attest™ Rapid Readout Biological Indicator | Vaporized Hydrogen Peroxide (VH2O2) | 30 each/bag | 4 bags/case |

Together we can raise the standard of patient care. For more information, contact your 3M representative or visit go.3M.com/sterilization.



Medical Solutions Division
3M Health Care
2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-228-3957
Web 3M.com/medical

*Refer to the indications for use provided in the biological indicator Instructions for Use. The indications for use should align with the sterilization cycle to be monitored.

© 2020 3M. All rights reserved. 3M and Attest are marks and/or registered marks of 3M. Unauthorized use prohibited. Please recycle. Printed in USA. EN_GL_70-2011-7940-8