

**3M** Science.  
Applied to Life.™



# Know now.

**Protect your patients  
and your practice with  
biological indicator  
results in 24 minutes\*.**



**3M™ Attest™ Mini Auto-reader 490M**  
for sterilization monitoring

**Now available! The 3M™ Attest™ Mini Auto-reader 490M.**

In-office monitoring (spore testing) can help reduce the risk of patient exposure to contaminated dental instruments. Sized and priced for office use, the Attest™ Mini Auto-reader 490M lets you bring rapid sterilization monitoring to your practice – because your patients deserve the highest standard of care.

- In-house sterilization monitoring results in just 24 minutes\*
- Know results before you release instruments for use – offering your patients a high standard of care
- Designed to be sized and priced for your practice
- Same expert 3M sterilization assurance technology as used in hospitals

\*When used with Attest™ Auto-reader 490 or 490H with software version 4.2.7 or greater.

# Your partner in patient safety.

- Easy to read display
- Small footprint for space-constrained areas
- Connect to a PC for digital documentation of biological indicator results



Catalogue No.	Product	Other Notes	Packaging
490M	<b>3M™ Attest™ Mini Auto-reader</b>	Dimensions (W x H x D): 15.5 cm x 8.25 cm x 5.25 cm (6.1 in x 3.2 in x 2.1 in)	1 unit/box
1491†	<b>3M™ Attest™ Super Rapid Readout Biological Indicator (blue cap) for Steam</b>	Gravity: 132-135°C (270-275°F) 3 min. and 10 min.	50 each/box (4 boxes/case)
1492V†	<b>3M™ Attest™ Super Rapid Readout Biological Indicator (brown cap) for Steam</b>	Dynamic-air-removal (Pre-vacuum and SFPP): 132-135°C (270-275°F) 3 min. and 4 min. 135°C (275°F) 3 min.	50 each/box (4 boxes/case)
1243A	<b>3M™ Attest™ Steam Chemical Integrator (CSA Type 5<sup>1</sup>)</b>	Dimensions: 5.1 cm x 1.9 cm (2 in x 0.75 in)	500 each/bag (2 bags/case)
1243B			100 each/bag (10 bags/case)
1322-18mm	<b>3M™ Comply™ Lead Free Steam Indicator Tapes for Disposable Wraps (CSA Type 1<sup>1</sup>)</b>	Dimensions: 18 mm x 55 m (0.70 in x 60 yds)	28 rolls/case
1322-24mm			20 rolls/case

To get started, contact your local 3M Canada Oral Care Sales Representative, authorized distributor, or call the 3M Oral Care Customer Care Centre at **1-888-363-3685**.

†Refer to the indications for use provided in the biological indicator Instructions for Use. The indications for use should align with the sterilization cycle to be monitored.

<sup>1</sup> As defined by CAN/CSA Z314-18.



Available in Canada from

**3M Canada Company**  
3M Medical Solutions Division  
300 Tartan Drive  
London, Ontario N5V 4M9  
Canada  
1-800-364-3577  
3M.ca

**3M Company**  
2510 Conway Avenue  
St. Paul, MN 55144  
USA  
1-800-228-3957

3M, 3M Science. Applied to Life., Attest and Comply are trademarks of 3M. Used under license in Canada. © 2020, 3M. All rights reserved. 2005-17352 E