

# PREVENA™

INCISION MANAGEMENT SYSTEM

## Protect Surgical Incisions, Minimise Postoperative Risk

PREVENA™ Therapy manages surgical incisions and surrounding tissue in patients at risk for developing post-operative complications, including surgical site infection.



### PREVENA™ 125 Therapy Unit

**Included with:**  
PREVENA™ 13cm, PREVENA™  
20cm and PREVENA DUO™  
System Kits

#### PEEL & PLACE™ Dressing – 13cm

- For linear incisions up to 13cm in length



#### PEEL & PLACE™ Dressing – 20cm

- For linear incisions up to 20cm in length



### PREVENA PLUS™ 125 Therapy Unit

**Included with:**  
PREVENA™ 35cm and  
PREVENA CUSTOMIZABLE™  
System Kits

#### PEEL & PLACE™ Dressing – 35cm

- For linear incisions up to 35cm in length



#### CUSTOMIZABLE™ Dressing

- For linear, non-linear and intersecting  
incisions up to 90cm in length



For further information, contact your local representative,  
or call Customer Services on **0800-980-8880**

# PREVENA™ Therapy System Kits

## PREVENA™ Therapy System Kit



## PREVENA CUSTOMIZABLE™ System Kit



- Ideal for mobile patients and those available for discharge
- Portable, single-use therapy for up to 7 days
- Patient carry case included
- Shower friendly dressings
- Audible and visual alarms

SIZE	3M+KCI CODE	NHSSC CODE	CONTENTS
13cm	PRE1101	ELZ850	1 x PREVENA™ 125 Therapy Unit, 1 x 13cm PREVENA PEEL & PLACE™ Dressing, Patch Strips, V.A.C.® Connector
20cm	PRE1001	ELZ283	1 x PREVENA™ 125 Therapy Unit, 1 x 20cm PREVENA PEEL & PLACE™ Dressing, Patch Strips, V.A.C.® Connector
35cm	PRE3201	ELZ984	1 x PREVENA™ PLUS Therapy Unit, 1 x 35cm PREVENA PEEL & PLACE™ Dressing, Patch Strips, V.A.C.® Connector
90cm	PRE4001	ELZ821	1 x PREVENA™ PLUS Therapy Unit, 1 x 90cm PREVENA CUSTOMIZABLE™ Dressing with SENSAT.R.A.C.™
DUO™ 13cm / 13cm	PRE1121	ELZ983	1 x PREVENA™ PLUS Therapy Unit, 2 x 13cm PREVENA PEEL & PLACE™ Dressings, 1 x V.A.C.® Y-Connector

# PREVENA™ Therapy Dressing Kits

## V.A.C.ULTA™ Therapy and ACTI.V.A.C.™ Therapy Sold Separately



- Compatible with V.A.CULTA™ Therapy Units in the acute setting
- Compatible with ACTI.V.A.C.™ Units for ambulatory patients

SIZE	3M+KCI CODE	NHSSC CODE	CONTENTS
13cm	PRE1155	ELZ833	5 x 13cm PREVENA PEEL & PLACE™ Dressings
20cm	PRE1055	ELZ647	5 x 20cm PREVENA PEEL & PLACE™ Dressings
35cm	PRE3255	ELZ985	5 x 35cm PREVENA PEEL & PLACE™ Dressings
90cm	PRE4055	ELZ822	5 x 90cm PREVENA CUSTOMIZABLE™ Dressings with SENSAT.R.A.C.™

# PREVENA™ Therapy Accessories

From left to right:  
45ml canister  
150ml canister  
V.A.C.® Connector



- Order additional canisters for high exudate wounds
- Replacement V.A.C.® Connectors

SIZE	3M+KCI CODE	NHSSC CODE	CONTENTS
45ml Canister	PRE1095	ELZ832	5 x 45ml PREVENA™ Canister
150ml Canister	PRE4095	ELZ823	5 x 150ml PREVENA PLUS™ Canister
V.A.C.® Connector	PRE9090	ELZ824	10 x PREVENA™ Therapy V.A.C.® Connector

**NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a physician and product instructions for use prior to application. This material is intended for healthcare professionals.**