

# SNAP™

THERAPY SYSTEM

## Attention users of Negative Pressure Wound Therapy (NPWT)

During these unprecedented times we are still available to offer you wound care solutions.

**Transition appropriate\* wound patients to SNAP™ Therapy System - It may help facilitate an easier and earlier<sup>1</sup> discharge to a home setting!**



**vac therapy**

**SNAP™**  
THERAPY SYSTEM

SNAP™ Therapy System's **off-the-shelf** availability, **simple application** process and **ultra-portability** are advantages in an outpatient care setting as compared to electrically powered negative pressure wound therapy.<sup>2</sup>


\* The SNAP™ Therapy System is indicated for patients who would benefit from wound management via the application of negative pressure, particularly as the device may promote wound healing through the removal of excess exudate, infectious material and tissue debris. For a detailed list of indications, please see Instructions for Use.

#### References


1. Harding K, et al. The role of mechanically powered disposable negative pressure wound therapy (dNPWT) in practice. Wounds International 2017. Available to download from [www.woundsinternational.com](http://www.woundsinternational.com).
2. Lerman B, Oldenbrook L, Ryu J, Fong KD, Schubart PJ. The SNAP™ wound care system: A case series using a novel ultraportable negative pressure wound therapy device for the treatment of diabetic lower extremity wounds. Journal of Diabetes Science and Technology. 2010 Jul 1;4:825-830.



- SNAP™ Therapy System combines the **simplicity** of advanced wound dressings with the **proven benefits of negative pressure wound therapy** in a **discreet design**.<sup>3</sup>
- It is **small, silent** and **lightweight** to preserve **patient quality of life**.<sup>3</sup>
- The **SNAP™ Therapy System** employs a familiar NPWT mechanism of action with reticulated open cell foam and -125mmHg pressure.<sup>4</sup>
- Designed for chronic wounds with exudate less than 180cc/week and less than 13cm x 13cm in size.



## Considerations for transitioning from V.A.C. VERAFLOR™ Therapy or V.A.C.® Therapy in acute setting to NPWT in the community.



For more information, please contact your local representative

#### References

3. Armstrong DG, Marston WA, Reyzelman AM, Kirsner RS. Comparative effectiveness of mechanically and electrically powered negative pressure wound therapy devices: a multicenter randomized controlled trial. *Wound Rep Reg.* 2012; 20(3):332-341
4. Fong KD, Hu D, Eichstadt S et al. The SNaP system: biomechanical and animal model testing of a novel ultraportable negative-pressure wound therapy system. *Plastic and Reconstructive Surgery.* 2010 May;125(5):1362-71.

**NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a physician and product instructions for use prior to application. This material is intended for healthcare professionals only.**

©Copyright 2020 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. PRA-PM-SG-00029 (07/20)

