

3M Science.
Applied to Life.™

3M™ Scotchbond™ Universal Plus
Adhesive


**More than just
a strong bond.**





3M™
Scotchbond™
Universal Plus
Adhesive takes
universal adhesives to
the next level.
It keeps all the benefits of the
original 3M™ Scotchbond™
Universal Adhesive, the first truly
universal adhesive, while adding
even greater control and
predictability – because a
good adhesive should
provide more than just a
strong bond.


From excellent...

Built on the clinically proven and world's most researched universal adhesive.

 **All direct and indirect indications**
light-curing, self-curing or dual-curing
methacrylate-based composites / cements

 **All etching techniques**
total-etch, self-etch, and selective
enamel etch

 **Bonds to all restorative materials**
without the need for separate primers

 **Available in vial and unit dose**
for efficient hygiene management


 **Virtually no post-operative sensitivity**
in field evaluation:


99% of dentists reported
that 3M™ Scotchbond™ Universal Plus
Adhesive shows virtually no post-
operative sensitivity.*

*Field evaluation with 309 dentists and over 20,000 applications

...to exceptional


Formulation carefully optimized to achieve unique features.

 **First radiopaque universal adhesive**
reducing the risk of X-ray misdiagnosis
and invasive overtreatment


 **Bonds and seals caries-affected dentin**
maximizing preservation of natural tooth structure,
supporting minimally invasive preparation guidelines

 **Gold standard adhesion to all dental substrates**
including glass ceramics

 **Full dual- and self-cure compatibility**
no need for an additional dual-cure activator bottle

 **BPA free formulation**
free of BPA derivatives like
BisGMA




 **Fully aligned system with 3M™ RelyX™
Universal Resin Cement**
excellent bond strength for virtually all resin
cement indications. Adhesive cured by cement.



Added radiopacity


On average, dentists see at least one hard to interpret X-ray per day due to radiolucency under an existing restoration. This could lead to misdiagnosis and invasive overtreatment. 3M™ Scotchbond™ Universal Plus Adhesive is the first radiopaque universal adhesive with dentin-like radiopacity, helping to overcome this issue.



Non-Radiopaque Adhesive

3M™ Scotchbond™ Universal Plus Adhesive

Intentional adhesive pooling under Class I restorations – X-rays of extracted human molars, filling material 3M™ Filtek™ One Bulk Fill Restorative. 3M internal data, H. Loll et al., J Dent Res 99 (Spec Iss A): No. 757, 2020



“Radiolucent adhesive pooling underneath a restoration doubles the risk of restoration replacement.”

**Dr. John O. Burgess,
University of Alabama, USA**

Ordering information

Item no. Product information

3M™ Scotchbond™ Universal Plus Adhesive – Vial

41293	Intro Kit Vial 1 Vial 5ml, 50 Disposable Brushes, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41294	Refill Vial – 1 × 5ml
41295	Refill Vial – 3 × 5ml
41263	3M™ Scotchbond™ Universal Etchant Etching Gel 2 Syringes x 3ml, 50 Dispensing Tips



Vial
with smooth and even surface –
easy to wipe off excess

3M™ Scotchbond™ Universal Plus Adhesive – Unit Dose

41297	Intro Kit Unit Dose 50 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41298	Unit Dose 100 100 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41299	Unit Dose 200 200 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips
41304	Unit Dose 400 400 Unit Doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips



Unit Dose
for efficient hygiene management –
reduces the risk of cross
contamination



3M.com

3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000
USA
Phone 1-800-634-2249
3M.com/dental

3M, Filtek, RelyX and Scotchbond are trademarks of 3M Company or 3M Deutschland GmbH.
© 3M 2020. All rights reserved.