


Preventing infections in intensive care units.

While vascular catheters provide the advantage of prolonged venous access, they present a risk of infectious complications. In fact, 60% of hospital-acquired bloodstream infections originate from some form of vascular access.¹ These infections can be acquired at the time of the initial insertion or anytime throughout the duration of the venous access.

Reducing infection risk at all access points with antimicrobial protection.



3M
Tegaderm™
Chlorhexidine Gluconate (CHG)
I.V. Securement Dressing

3M
Curoso™
Disinfecting Caps

Extraluminal protection

Intraluminal protection

3M™ Tegaderm™ CHG Chlorhexidine Gluconate I.V. Securement Dressing

Cleared and proven to reduce Catheter-related Bloodstream Infection (CRBSI) through immediate and consistent antimicrobial activity while providing continuous site visibility, consistent application and catheter securement.

3M™ Curoso™ Disinfecting Caps

Consistent use of Curoso disinfecting caps on IV needleless connectors is associated with decreased Central Line Associated Bloodstream Infection (CLABSI). Disinfect and protect needleless connectors, open female luers and male luer devices to help reduce the risk of contaminants entering the catheter post-insertion.

3M™ Tegaderm™ Chlorhexidine Gluconate (CHG) I.V. Securement Dressing

The gel pad provides 2% CHG to the skin surface immediately, without requiring moisture to activate. The integrated design ensures consistent application, aligning with evidence-based guidelines and practice standards.



Antimicrobial I.V. securement dressing



3M™ Curoso™ Disinfecting Cap for Needleless Connectors

Provides quick and verifiable disinfection of needleless connectors. Fits commonly used needleless connectors.

3M™ PICC/CVC Securement Device + Tegaderm™ CHG I.V. Securement Dressing

An engineered stabilisation device (ESD) plus antimicrobial (CHG) dressing cleared and proven to reduce CRBSI.



Sutureless securement device




3M™ Curoso™ Tips Disinfecting Cap for Male Luers

Disinfects and protects the distal end of I.V. tubing and other male luer devices.

3M™ Tegaderm™ CHG Chlorhexidine Gluconate I.V. Port Dressing

Antimicrobial (CHG) gel pad plus I.V. port dressing specifically designed to protect single or double implanted venous ports and non-coring 'Huber' needles from pathogens most commonly found in CRBSIs.



I.V. securement dressing

CHG gel pad



3M™ Curoso™ Stopper Disinfecting Cap for Open Female Luers

Designed to fit on a wide range of stopcocks and catheter hubs. Their unique design maintains pressure and disinfects the critical areas with 70% (w/v) isopropyl alcohol.



3M™ Curoso™ Disinfecting Cap for Tego® Hemodialysis Connectors

This specially designed Curoso disinfecting cap is compatible with Tego® Needlefree Hemodialysis Connector.



3M is here to support you and your ICU teams.

Ordering information

Product code	NHS code	Product description	Size in cm	Dressings per box
1657R	ELW294	3M™ Tegaderm™ CHG I.V. Securement Dressing	8.5 × 11.5	25
1658R	ELW625	3M™ Tegaderm™ CHG I.V. Securement Dressing	10 × 12	25
1659R	ELW295	3M™ Tegaderm™ CHG I.V. Securement Dressing	10 × 15.5	25
1660R	ELW366	3M™ Tegaderm™ CHG I.V. Securement Dressing	7 × 8.5	25
1877R-2100	ELW944	3M™ Tegaderm™ PICC/CVC Securement Device + CHG	8.5 × 11.5	20
1879R-2100	ELW945	3M™ Tegaderm™ PICC/CVC Securement Device + CHG	10 × 15.5	20
1665R		3M™ Tegaderm™ I.V. Port Dressing + 6.2 × 4.9 gel pad device	12 × 12	25

Product code	NHS code	Product description	Dispenser	Units per box
CFF1-270R	FSW582	3M™ Curoso™ Disinfecting Caps for Needleless Connectors	Individual caps	270
CFF10-250R	FSW583	3M™ Curoso™ Disinfecting Caps for Needleless Connectors	Strips (10 caps)	25 strips
CTG1-270R	FSW1355	3M™ Curoso™ Disinfecting Caps for Tego® Hemodialysis Connectors	Individual caps	270
CSV1-270R	FSW1376	3M™ Curoso™ Stopper Disinfecting Caps for Open Female Luers	Individual caps	270
CSV5-250R	FSW1377	3M™ Curoso™ Stopper Disinfecting Caps for Open Female Luers	Strips (5 caps)	50 strips
CM5-200R		3M™ Curoso™ Tips Disinfecting Cap for Male Luers	Strips (5 tips)	40 strips

To find out more:

 [3M.co.uk/medical](https://www.3m.co.uk/medical)

 **0330 053 8938**

 C3SD@mmm.com

References

- 1 Crnich CJ, Maki DG. The Role of Intravascular Devices in Sepsis. *Curr Infect Dis Rep.* 2001 Dec;3(6):496–506.
- 2 Timsit JF, Mimoz O, Mourvillier B, et al. Randomized Controlled Trial of Chlorhexidine Dressing and Highly Adhesive Dressing for Preventing Catheter-related Infections in Critically Ill Adults. *American Journal of Respiratory and Critical Care Medicine.* 2012; 186: 1272–1278.

3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough
LE11 5RB
01509 611611

3M Ireland
The Iveagh Building
Carrickmines Park
Carrickmines
Dublin 18
00 353 (01) 280 3555

Tego® is a trademark of ICU Medical, Inc. 3M, Curoso and Tegaderm are trademarks of the 3M company.
© 3M 2020. All rights reserved. OMG103039.

