

V.A.C. FREEDOM™

T H E R A P Y S Y S T E M

A portable system for managing wounds



Easy to Operate

- LCD with menu interface allows for efficient selection of negative pressure wound therapy
- SENSAT.R.A.C.™ Technology monitors and maintains wound site target pressure
- Smart Alarms help ensure patient safety and assist staff in providing care
- 300mL canister minimizes canister changes and is easily removed and replaced

Patient Comfort

- Long battery life (~12 hours) and lightweight design enable patients to be mobile and active
- Carrying case allows discreet delivery of therapy on the go
- Adjustable rate of dressing draw down intensity for increased comfort

Therapy Options

Negative Pressure:

- 50-200 mmHg

Mode of Operation:

- Continuous and Intermittent

Intensity:

- Variable

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a physician and product instructions for use prior to application. Rx only.

©Copyright 2020 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.
LIT# 29-A-438 • PRA-PM-US-02378 (03/20)

Specifications

- V.A.C. Freedom™ Therapy Unit: P/N 320000

Dimension:

- 6.5"(16.51cm) W x 3.125"(7.9cm) H x 7.5"(19.05cm) D

Weight:

- 3.20 lbs. (1.45 kg)

Canister volume:

- 300 mL

Battery Life:

- Approximately 12 hours
Electrical Charger/External Supply Input:
100-240 VAC, 50/60 Hz, 1 Amp

Charger/External Supply Output:

- 18 VDC, 2.2 Amp

Maximum Electrical Leakage:

- <100 Microamps

IEC Classification:

- Medical Equipment
- Equipment not suitable for use in Presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide

Storage Conditions

Temperature Range..... -4°F (-20°C) to 140°F (60°C)
Relative Humidity Range.. 10 to 95% non-condensing

Operating Conditions

Temperature Range..... 50°F (10°C) to 85°F (30°C)
Relative Humidity Range.. 30% to 75%